3M ESPE

Protemp™ 4

Temporization Material
Временный материал
Материал за временна употреба
Materijal za izradu privremenih restauracija
Ideiglenes
Materiał tymczasowy
Materiál pentru protezări temporare
Materiál pre dočasné náhrady
Materiál za izdelavo začasnih restavracij
Materiál pro dočasné náhrady
Geçici Restorasyon Materyali
Ajutište restauratsioonide komposiitmaterjal
Materiāls pagaidu protezēšanai
Laikinoji medžiaga

Instructions for Use
Инструкция по применению
Указания за употреба
Upute za uporabu
Használati utasítás
Instrukcja użycia
Návod na použitie
Navodila za uporabo
Kasutusjuhend
Lietošanas pamācība
Vartojimo instrukcija
Návod na použitie
Návod k použití
Kullanma Kılavuzu
Lietošanas pamācība
Vartojimo instrukcija

3M ESPE AG
Dental Products
D-82229 Seefeld – Germany

3M ESPE
Dental Products
St. Paul, MN 55144-1000

3M, ESPE, Filtek, Garant, Penta, Protemp, RelyX and Soft-Sec are trademarks of 3M or 3M ESPE AG.
© 3M 2008. All rights reserved.

64002073/0708

0123

ENGLISH

Product Description
Protemp™ 4, manufactured by 3M ESPE, is a composite material for the fabrication of temporary restoration at the chairside. Based on multi-functional methacrylate esters, this two-component system has filling material-like properties and is available in the following colors: A1, A2, A3, B3, Bleach. With the Garant Dispenser, manufactured for 3M ESPE, the material can be applied directly from the cartridge without causing any bubbles to form. Protemp™ 4 temporary restorations can be modified and individualized using low-viscosity composites (e. g., from the Filtek™ Flowable Composite product line, manufactured by 3M ESPE).

☞ For details on all mentioned products please refer to the corresponding Instructions for Use. Instructions for Use should not be discarded for the duration of product use.

Application Areas
- Fabrication of temporary crowns, bridges, inlays, onlays and veneers
- Fabrication of long-term provisional restorations
- Lining material for prefabricated temporary crowns made of composite (e. g., Protemp™ Crown, manufactured by 3M ESPE) and metal (e. g., Iso-Form crowns, manufactured by 3M ESPE)

Precautionary Measures
MSDSs for the mentioned 3M ESPE products can be obtained at www.mmm.com or from your local subsidiary.

Preparatory Work
- Take an alginate, silicone or polyether impression from the intact row of teeth (e. g., using a product manufactured by 3M ESPE).
 – Instead of an impression, a laboratory-made vacuum-formed template or, for single crowns, a preformed crown can be used for shaping.
- To enhance the strength of the temporary restoration, remove proximal grooves in the impression. If required due to tight spacing conditions, enlarge the interior of the mold.
- For easy replacement in the mouth, relieve undercuts.
- With missing teeth or gaps in the molar area, cut a groove in the impression to produce a stable, bar-shaped connection for improved temporary restoration.
- If required, prior to taking the impression close any gaps between front teeth with denture teeth acting as spacers; interlock multiple denture teeth with wax.

Times
00:00–00:40 min:sec	Place Protemp™ 4 in the impression and reposition in the mouth
00:40–01:40 min:sec	Setting in the mouth
01:40–02:50 min:sec	Best point in time for removal of the impression
Up to 05:00 min:sec	Complete setting
From 05:00 min:sec	Removal of the temporary restoration from the impression, finishing

The processing times are applicable to products which have been stored and processed at a temperature of 23°C/74°F and at 50% relative humidity. Higher temperatures shorten these times, while lower temperatures extend them.

Dosing and Mixing
The pastes are dosed and statically mixed in the Garant dispenser.
Use the enclosed Garant™ Mixing tips blue only! Garant Mixing tips are manufactured for 3M ESPE.

During the initial Dosing
- **Remove and discard the cartridge lid.** Under no circumstances should the cartridge lid be used to reseal the cartridge as this could result in the inclusion of air bubbles.
- Check whether both cartridge openings are free from obstruction and remove obstructive material, if any, with a suitable instrument.
- Before a mixing tip is applied for the first time, squeeze paste out onto a mixing block until the base and catalyst pastes are the same. Discard the paste which has been squeezed out.

Application
It takes the material a certain amount of time to emerge from the newly attached mixing tip. The flow of material is interrupted as soon as the pressure on the handle is reduced. Material that has set inside the mixing tip should not be extruded by force, as this may damage both the cartridge and the mixing tip.
Core build-ups which are fabricated using composites, compomers, or resin-modified glassionomer cements must be sealed off before the application of Protemp™ 4, e. g., using Vaseline.
- First, install a new mixing tip.
- Prior to each application, extrude a small (pea-sized) amount of paste from the mixing tip onto a mixing pad and discard this material.
- Subsequently, load the dried impression or vacuum-formed template in the relevant spaces from bottom up.
- Reposition the impression or vacuum-formed template in the mouth.
- The material attains a hard-elastic consistency within 40 sec after the onset of mixing. The material and the impression or vacuum-formed template must be removed from the mouth within 2 min 50 sec after the onset of mixing.
- Check the setting process by means of the excess material in the mouth.
- Keep the filled mixing tip on the cartridge to serve as a closure until the next use.
- Remove the smear layer and finish the residues, see "Finishing".

Cementation
- If the definitive restoration has been sealed off (e. g., using Vaseline), remove the rest of the sealant thoroughly before the temporary cementation.
- Insert the temporary restoration with a commercial temporary cement, e. g., RelyX™ Temp E or RelyX™ Temp NE, both products manufactured by 3M ESPE.
 – Eugenol-containing cements may impair the setting of composite luting cements to be used at a later time in the final luting step! If a composite cement is to be used for the final luting step, seat the temporary restoration with a eugenol-free cement.

After Cementation
- Thoroughly examine and explore the sulcus and surrounding dentition.
- Remove any residuals of the temporary material or cement from the mouth.

Repair
- Provide for mechanical retention by roughening the break surfaces and adjacent areas.
- Clean with alcohol and dry with air.
- Join together using Protemp™ 4 or a Filtek Flowable Composite.
 – Cure Protemp™ 4 for 90 sec. Filtek Flowable Composite as described in the relevant Instructions for Use.
 – Once the material is fully set, finish as usual.

Cleaning of the Garant Dispenser
- Remove the Garant™ cartridge.
- Remove paste that has not set using an alcohol-soaked cloth.
- The Dispenser handle and plunger can be autoclaved up to a temperature of 135°C/275°F; remove the plunger beforehand.
- Glutaraldehyde-based solutions can be used for disinfection.

Storage and Stability
Store the product at 15–25°C/59–77°F.
Do not store in the refrigerator.
Do not use after the expiration date.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of this product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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