3M | Impregum™ Garant™ L DuoSoft™ ESPE

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3M Deutschland GMBH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss – Germany

3M, ESPE, DuoSoft, Garant, Impregum and Penta are trademarks of 3M or 3M Deutschland GmbH. © 2020, 3M. All rights reserved.

MD CE 0123

44000118057/02

en ENGLISH

Product Description
Impregum™ Garant™ L DuoSoft™ is a light-bodied consistency polyether impression material supplied in the Garant cartridge. The mixing ratio is 2 volumes base paste:1 volume catalyst paste.
Due to its lower shore hardness, the impression material is suitable for use with impressions with slight undercuts even without additional blocking out.

☞ These Instructions for Use should be kept for the duration of product use. The product may only be used when the product documentation is available. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

Intended Purpose
Intended purpose: dental impression material
Intended users: educated dental professionals who have theoretical and practical knowledge on usage of dental products
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use
Clinical benefit: create a precise reproduction of the clinical situation within dental workflows

Indications
- Impression of inlay, onlay, crown, bridge and veneer preparations
- Functional impressions
- Fixation and Implant impressions

Contraindications
Polyether impression materials are contraindicated for use with patients who have a known allergy to polyether impression materials or contained allergens, as listed in section Precautionary Measures.

Precautionary Measures

For Patients
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue use of the product in this patient and seek medical attention as needed.
Known allergens are 1-Dodecylimidazole, mint flavors and Sulfonium salt.

For Dental Personnel
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In case of prolonged, avoid exposure to uncured material. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.

Potential Undesirable Side Effects and Complications
Potential effects are local irritations, or local or systemic allergic reactions. In rare cases impression material can be swallowed or aspired. Swallowing of impression material can cause harm to the gastrointestinal tract. As the impression material cures in the mouth it gets hard. It is possible that loose teeth or restorations are removed with the impression material.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (ies) of local regulatory authority.

3M Safety Data Sheets/ Safety Information Sheets can be obtained from www.3M.com or contact your local subsidiary.

Preparation

Garant dispenser/cartridge:
▶ Upon use of a new cartridge and with each new mixing process, start without attaching the mixing tip by extruding a small quantity of paste until both components emerge uniformly.
▶ Remove clogged paste, if necessary.
▶ Subsequently, install the Garant Mixing Tip White and Garant Intraoral Tip White.

Impression trays:
All impression trays generally used for precision impressions are suitable.
▶ For sufficient adhesion, apply a thin layer of Polyether Adhesive/Polyether Contact Tray Adhesive to the tray and allow to dry completely according to the Instructions for Use.

Retraction
Solutions based on aluminum hydroxide chloride or aluminum sulphate, are suitable retraction agents. Retraction with threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulphate or iron (III) sulphate may impair the setting of polyether impression materials.
▶ Areas from which impressions are to be taken should be kept dry.
▶ In subgingival preparations, hemostatic threads or rings may be used.
▶ **Prior to taking the impression, completely remove any residual retraction agent by rinsing and drying.**

Times

	Processing time from start of mixing* min:sec	Setting time from start of mixing* min:sec	Residence time in the mouth** min:sec
Impregum Garant L DuoSoft	2:00	5:30	3:30
Impregum Penta H DuoSoft	2:30	6:00	3:30

The working times shown are in conformity with ISO 4823. Higher temperatures will shorten the working time, while lower temperatures will prolong it. We recommend a working time of 1:45 min:sec and compliance with the indicated setting time from the start of mixing to assure an ideal impression result at all times under clinical conditions.
* Start of mixing = entry of paste into the mixing tip.
** For product combinations, longer times are valid.

Impression Taking

One-step technique
▶ Load the selected material onto the impression tray according to the corresponding manufacturer instructions.
▶ At the same or a later time, apply Impregum Garant L DuoSoft around the dry preparation from bottom up, keeping the tip of the intraoral tip immersed in the material at all times, and in contact with the surface of the tooth.
– If required, fill a separate Elastomer syringe directly from the mixing tip **without** attaching the intraoral tip.
▶ Without delay, place the loaded tray in the mouth and hold it in place without applying pressure until setting is complete.
▶ To remedy initial adhesion ("setting the valve"), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior position. If this proves difficult it may be necessary to carefully blow some air or water between the impression and the gingiva.

Functional impression
▶ After taking a functional impression with a heavy-bodied or medium-bodied consistency polyether material, use Impregum Garant L DuoSoft for precise impression taking (correction).

After Impression Taking
▶ Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene
▶ Place the impression in a standard disinfectant solution commonly used for impressions, for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
▶ After disinfection, rinse the impression under running water for approx. 15 sec.

Model Preparation
▶ Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
▶ To avoid introducing bubbles into the model, briefly pre-rinse with water, and dry with air before pouring.
▶ **Do not use surfactants as these impair the quality of polyether impressions and are not necessary!**
▶ Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

Cleaning and Disinfection of the Garant Cartridge
Caution
Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:
- Oxidation agents (e.g., hydrogen peroxides)
- Oils
- Glutaral/aldehydes

Use only validated methods for cleaning and disinfection. Use solely cleaning and disinfectant agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

Limitations on Reprocessing
The Garant cartridge has been tested 18 times without suffering visible damage. The Garant cartridge may be reused, provided that it is undamaged and clean.

General
This cartridge is non-sterile when delivered and must be cleaned and disinfected before initial use and prior to every further use.

Preparation for Decontamination
The cartridge should be reprocessed as soon as possible after every use.

Reprocessing
The reprocessing procedure must be conducted before the initial use and after every further use. The cartridge must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

Manual Cleaning by Wiping
▶ Clean the cartridge before its first use and after every further use with ready-to-use disinfectant wipes (CaviWipes®) or equivalent cleaning agents. The wipes contain the disinfectant CaviCide® (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use.
▶ Continue wiping the cartridge with the wipes until no contamination is visible.
▶ Special care must be exercised in cleaning the grooves.

Manual Disinfection by Wiping
▶ Disinfect the cartridge with ready-to-use disinfectant wipes (CaviWipes) or equivalent disinfection agents. The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.
▶ Disinfect the cartridge for three minutes using at least one new disinfectant wipe. Keep all surfaces of the entire cartridge moist during the three-minute period.
▶ Special care must be exercised in disinfecting the grooves.

Inspection, Maintenance, Testing
▶ Examine the cartridge for damage, discoloration, and contamination before every use.
▶ Do not use damaged cartridges under any circumstances.
▶ If the cartridge is not visibly clean, repeat the reprocessing procedure.

Cleaning the Garant Dispenser and Elastomer Syringe
Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Garant Dispenser and the Elastomer Syringe.

Notes
▶ Garant Mixing Tips White, Garant Intra-oral Tips White, Elastomer Syringe Tips are single-use devices and would be damaged when removing set impression material. When reattaching Garant Intra-oral Tip White, the proper fit of the tip on the mixing tip cannot be ensured.
▶ Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.
▶ **Turning or reattaching the mixing tip without subsequent impression taking may lead to carry-over of paste and clogging.**
▶ At storage temperatures below 8 °C/46 °F, the paste may contract such way that precise dosing cannot be guaranteed even at room temperature.
▶ Direct exposure to sunlight and damp storage conditions may damage the impression material.
▶ Polyether impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modelling.
▶ 3M Polyether materials may only be combined with 3M polyether materials. Doing so the working and setting times of the individual products must necessarily be kept.

Disposal
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

Storage and Stability
Store the product at 18–25° C/64–77° F. **Do not refrigerate!**
Do not use after the expiration date.
Store impressions dry and below 30° C/86° F in the dark.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark	CE 0123	Indicates conformity to European Union Medical Device Regulation or Directive with notified body involvement.
Reference Number and Symbol Title	Symbol	Description of Symbol
CE Mark	CE	Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device	MD	Indicates the item is a medical device.
Recycle	21 PAP	Indicates packaging is recyclable with mixed paper.)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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