For better dentistry

DENTSPLY

prime&bond® one

ETCH & RINSE

Universal Etch&Rinse Adhesive

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ETCH & RINSE

Universal Etch&Rinse Adhesive

CAUTION: For dental use only.
USA: Rx only.

Content Page

1 Product description **2**
2 Safety notes **4**
3 Step-by-step instructions **6**
4 Hygiene **9**
5 Lot number and expiration date **10**

1 Product description

prime&bond® one etch&rinse adhesive is a universal etch-and-rinse one-bottle dental adhesive designed to bond resin-based materials to enamel and dentin as well as to metals and ceramic.

prime&bond® one etch&rinse adhesive offers excellent performance due to high bond strength to enamel and dentin, easy and comfortable application and a high degree of technique robustness.

prime&bond® one etch&rinse adhesive, when used on its own, is designed to be used with light-cured resin-based restorative materials and cements.

prime&bond® one etch&rinse adhesive, when mixed with Self Cure Activator (available separately), is designed to be used with dual-cure/self-cure resin materials. For the use in conjunction with Self Cure Activator please refer to the Instructions for Use “Self Cure Activator”.

1.1 Delivery forms

- Easy-squeeze bottle

1.2 Composition

- Carboxylic acid modified dimethacrylate (TCB resin)
- Phosphoric acid modified acrylate resin (PENTA)
- Urethane dimethacrylate (UDMA)
- Triethyleneglycol dimethacrylate (TEGDMA)
- 2-Hydroxyethylmethacrylate (HEMA)
- Butylated benzenediol (stabilizer)
- Ethyl-4-(dimethylamino)benzoate
- Camphorquinone
- Functionalised amorphous silica
- Tertiary butanol

1.3 Indications

Self-priming adhesive for:

- Direct, light-cured composite and compomer restoratives.
- Indirect restorations; light-cured, resin cemented inlays, onlays, crown and bridge retainers.
- Composite/compomer, ceramic and amalgam repairs.
- Cavity varnish for use with fresh amalgam.

For the use of prime&bond® one etch&rinse for the following indications please refer to the Instructions for Use "Self Cure Activator":

- Direct, dual-cure or self-cure composite restorations and core build-ups.
- Indirect restorations; dual-cured and self-cured resin cemented inlays, onlays, crown and bridge retainers.
- Dual-cured and self-cured resin cemented endodontic post cementation.

For the use of prime&bond® one etch&rinse for adhesive bonding of direct amalgam restorations please refer to the Instructions for Use "Amalgam Bonding Acessory Kit" [1].

1.4 Contraindications

- Use with patients who have a history of severe allergic reactions to acrylate or methacrylate resins or any of the other components.
- Direct application to dental pulp (direct pulp capping).

1.5 Compatible tooth conditioners

prime&bond® one etch&rinse adhesive is used following phosphoric acid conditioning with conventional dentin/enamel conditioners including DENTSPLY 34% and 36% conditioners (see Instructions for Use of selected conditioner).

1.6 Compatible restorative materials

prime&bond® one etch&rinse adhesive is used in conjunction with suitable restorative or cement material and is chemically compatible with conventional methacrylate-based composite restorative materials including DENTSPLY visible light-cured composite restorative or cement materials (see Instructions for Use of selected restorative or cement material).

[1] Not available in all countries.

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

The material contains polymerizable acrylate and methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Tightly close bottle immediately after use.
- Use in a well-ventilated area. Avoid inhaling vapor.
- Flammable: Keep away from sources of ignition. Take precautionary measures against static discharges.
- Avoid the product saturating the gingival retraction cord. If it soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.
- Applicator tips are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.

- Interactions:
 - Do not use eugenol- and hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening of the product.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
 - prime&bond® one etch&rinse adhesive is a light-cured material. Protect from ambient light.
 - Insufficient data exist to support the use of desensitizing agents and/or cavity cleansing agents with prime&bond® one etch&rinse adhesive. Recommended pre-treatment is pumice or prophylaxis paste with a rubber cup. (See Step-by-Step Instructions)
 - If H_2O_2 has been used to clean the cavity, proper rinsing is essential. Higher concentration H_2O_2 may interfere with the setting of polymerizable material and should not be used prior to the application of prime&bond® one etch&rinse adhesive.

2.3 Adverse reactions
- Eye contact: Irritation and possible corneal damage.
- Skin contact: Sensitization. Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.
- Keep the product out of direct sunlight and store in a well ventilated place.
- Store the product at temperatures between 10 °C and 28 °C (50 °F and 82 °F). Use the product at room temperature.
- Protect from moisture.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 Direct restoration (light-cured composite resin and compomers)

3.1.1 Cleaning

1. Clean uninstrumented enamel and dentin with a rubber cup and pumice or a non-fluoride cleaning paste such as Nupro® prophylaxis paste.
2. Wash thoroughly with water spray and air dry.
3. Clean freshly instrumented enamel and dentin with water spray, then air dry.

3.1.2 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Instructions for Use) leaving the rest of the cavity surface free for bonding.

3.1.3 Tooth conditioning/dentin pre-treatment

When used as a bonding agent for direct restorative materials, it is recommended to follow the total etch technique described in this chapter. Prior to use, please see the Instructions for Use of the respective conditioner.

Application of 34%-36% phosphoric acid conditioner

1. After application of rubber dam or other suitable isolation technique, apply the respective conditioner.
2. Apply the respective conditioner to the cavity surfaces starting at the enamel margins.
3. For best results, condition enamel for at least 15 seconds and dentin for 15 seconds.

Rinsing and blot drying

1. Remove conditioner with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.
2. Remove rinsing water completely by blowing gently with an air syringe or by blot drying with a cotton pellet. Do not desiccate dentin.
3. Proceed immediately to application of prime&bond® one etch&rinse adhesive.

Contamination.

Restoration failure.

1. Once the surfaces have been properly treated, they must be kept uncontaminated.
2. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure for 5 seconds only. Rinse and dry as described above.

3.1.4 Dosage

1. To prevent the bottle from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves.
2. Dispense prime&bond® one etch&rinse adhesive directly onto a clean applicator tip, making sure that the bottle does not come in direct contact with the applicator tip, or place one to two drops of prime&bond® one etch&rinse adhesive into a clean CliXdish™ mixing well or standard dappen dish.
3. Replace cap promptly.

Technique Tip: Material in a closed CliXdish™ mixing well will remain usable for up to 15 minutes. Material dispensed to a standard mixing well must be used immediately.

3.1.5 Application and curing

1. Using the applicator tip, apply the product wetting all cavity surfaces uniformly. Avoid pooling. This surface should remain undisturbed for 20 seconds.
2. Evaporate solvent by thoroughly drying with clean, dry air from a dental syringe for at least 5 seconds. Use a vacuum aspirator to prevent spatter of the adhesive onto the mucosa. Surface should have a uniform glossy appearance. If not, repeat application and air dry.
3. Cure prime&bond® one etch&rinse adhesive for 10 seconds [2] using a curing light.

Insufficient curing.

Inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output.

4. Place restorative material over the cured prime&bond® one etch&rinse adhesive as per restorative material manufacturer's Instructions for Use.

3.2 Light-cured composite, ceramic and amalgam repairs

3.2.1 Preparation

1. Roughen and clean fractured restoration surface with a fine diamond bur and create mechanical retention, if possible.

3.2.2 Tooth conditioning/dentin pre-treatment

1. Etch fractured tooth surface with the respective conditioner as outlined in 3.1.3 Tooth conditioning/dentin pre-treatment.

3.2.3 Treatment of restoration

1. Etch fractured restoration surface with the respective conditioner for at least 15 seconds as outlined in 3.1.3 Tooth conditioning/dentin pre-treatment.

3.2.4 Application and curing

1. Apply adhesive to the fractured restoration surfaces and light-cure [2] the adhesive [3.1.5 Application and curing].

3.2.5 Completion

1. Complete repair with placement and curing of desired shade(s) of light-cured composite restorative following manufacturer's Instructions for Use.

3.3 Cementation of light transmitting indirect restoration with a light-cured resin cement

3.3.1 Treatment of restoration

1. Treat internal bonding surface of restoration according to manufacturer's or dental laboratory's Instructions for Use, i.e., etching, mechanical roughening and/or Silane Coupling Agent application.

[2] Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 500 and 800 mW/cm².

2. Apply a single coat of prime&bond® one etch&rinse adhesive to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.

Light-curing of applied prime&bond® one etch&rinse adhesive is not necessary.

3.3.2 Treatment of tooth preparation

1. Clean and condition tooth preparation, apply and light-cure the material as described in 3.1 Direct restoration (light-cured composite resin and compomers).

3.3.3 Cementation

1. Prepare and apply light-cured resin cement according to manufacturer's Instructions for Use.

3.4 Cavity varnish for use with fresh amalgam

When used as a cavity varnish, prime&bond® one etch&rinse adhesive is not an amalgam adhesive.

3.4.1 Preparation

1. Finish preparation by removing existing restorations and/or caries.

3.4.2 Cleaning

Clean preparation and place a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Instructions for Use) if needed, according to:

1. Cleaning [3.1.1 Cleaning].
2. Pulp protection [3.1.2 Pulp protection].

3.4.3 Tooth conditioning/dentin pre-treatment

1. Rinse and carefully air dry cavity preparation, but do not desiccate exposed dentin.
2. Condition enamel and dentin following section 3.1.3 Tooth conditioning/dentin pretreatment.

3.4.4 Application

1. Apply and light cure prime&bond® one etch&rinse adhesive as described for direct restorations, see section 3.1.3 Tooth conditioning/dentin pretreatment.

3.4.5 Completion

1. Place and condense amalgam (e.g. Dispersalloy® Dispersed Phase Alloy, see Instructions for Use) as per manufacturer's directions.

4 Hygiene

Cross-contamination.

Infection.

1. Do not reuse single use products. Dispose in accordance with local regulations.
2. Reprocess reusable products as described below.

4.1 CliXdish™ mixing well

4.1.1 Manual cleaning

1. Discard used gloves according to local regulations and disinfect hands with an appropriate bactericidal, virucidal and fungicidal hand disinfectant solution according to local regulations and use according to disinfectant solution manufacturer's Instructions for Use.
2. Use a new clean pair of examination gloves.
3. Disassemble device.
4. Remove residual adhesive with tissue.
5. Clean device with a soft brush, warm water, and an appropriate cleaning solution until free of visible contamination. Pay special attention to device seams and insertions.
6. Remove cleaning solution residue with a cloth soaked with water.
7. Dry device with a lint-free single-use cloth.

4.1.2 Manual disinfection

1. After cleaning [4.1.1 Manual cleaning] thoroughly wipe all device surfaces with a single-use cloth in combination with a water-based, bactericidal, virucidal and fungicidal disinfection solution approved according to local regulations and use according to disinfectant solution manufacturer's Instructions for Use. Make sure disinfectant solution is compatible with cleaning solution. Pay special attention to device seams and insertions.
2. Remove disinfectant solution residue with a cloth soaked with water.
3. Dry device with a lint-free single-use cloth.

Device can be manually disinfected for up to 200 times.

4.1.3 Automatic Washer-Disinfector

As an alternative to manual cleaning and disinfection, the device can also be processed in a washer-disinfector [3]. Use only properly maintained, calibrated, and approved washer-disinfector according to ISO 15883-1.

1. Disassemble device.
2. Remove residual adhesive with tissue.
3. Place device in the washer-disinfector allowing water and detergent to enter into and drain off from device openings during washer-disinfection program and follow washer-disinfector's Instructions for Use.
4. Run washer-disinfection program with A0 value ≥ 3000 (e.g. 5 min at ≥ 90 °C) using appropriate detergents [3].

[3] For example Miele Washer-Disinfector G7835 CD using program vario TD with neodisher® MediClean (alkaline detergent) and neodisher® Z (acid neutralization and cleaning detergent); both detergents from Dr. Weigert, Hamburg, Germany (no products or registered trademarks of DENTSPLY International, Inc.).

Device can be washer-disinfected for up to 500 times.

4.1.4 Sterilization (optional)

The device can also be autoclaved.

1. Steam autoclave device unwrapped at 134 °C/ 2 bar = 273 °F/ 29 psi for 3 minutes 30 seconds.

Device can be disinfected for up to 200 times and can be sterilized for up to 20 times.

5 Lot number (LOT) and expiration date ()

1. Do not use after expiration date. ISO standard is used: “YYYY-MM”
2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number
 - Expiration date

© DENTSPLY DeTrey 2014-02-12
[These Instructions for Use are based on Master Version 01]

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DENTSPLY DeTrey GmbH
DENTSPLY Limited