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Simboli

ENGLISH

RAPID PGLA
Synthetic absorbable undyed braided multifilament surgical suture polyglactin 910 coated with a co-polymer polyglactin 370 + Calcium Stearate

Description
RAPID PGLA suture is a braided synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide. The characteristic rapid loss of strength occurs due to the lower molecular weight than the normal PGLA. RAPID PGLA suture is available undyed in a wide variety of length(45cm-150cm) – diameter(USP 6/0-USP 2) combinations, with or without needle of various sizes and types, manufactured from high quality medical grade stainless steel.
RAPID PGLA meets all requirements established by the European Pharmacopeia for sterile Synthetic Absorbable Braided Sutures and with the essential requirements of the directive 93/42/EEC.

Indications
RAPID PGLA Synthetic Absorbable Surgical Suture is indicated only for use in soft tissue approximation of the skin and mucosa, plastic surgery, where only short-term wound support is required. It is not indicated for use in ligation, ophthalmic, cardiovascular, or neurological procedures.

Actions
RAPID PGLA Synthetic Absorbable Surgical Suture elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of RAPID PGLA Synthetic Absorbable Sutures occurs by means of hydrolysis, where the polymer degrades to glycolic and lactic acids which are subsequently absorbed and metabolized in the body. RAPID PGLA retains approximately 50% of the original tensile strength at 5 days post implantation. All of the original tensile strength is lost by approximately 10 to 14 days. Absorption is essentially complete by 42 days.

Contraindications
This suture, being absorbable, should not be used where extended approximation of tissue is required.

Warnings / Precautions / Interactions
RAPID PGLA surgical sutures should be used only from professionals who are familiar with surgical procedures and techniques involving absorbable sutures and wound closure techniques, as the risk of wound dehiscence may vary with the site of application and the suture material used. When selecting a suture its in vivo performance should be taken into consideration. The use of RAPID PGLA sutures may be inappropriate in elderly, malnourished and debilitated patients or patients suffering from other conditions that may delay wound healing.
As with any other absorbable surgical suture RAPID PGLA may have the temporary act of a foreign body. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or inflected wounds.
Topical skin sutures which must remain in place longer than 7 days may cause localized irritation and the external portion should be cut and snipped off or removed as indicated.
Adequate knot security is achieved with the standard surgical technique of flat and square ties, with additional throws depending on the surgical circumstances and the experience of the surgeon. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
Special care should be taken in handling surgical needles. The needle holder should grasp the needle in an area which should not exceed the 1/3 to 1/2 of the total length of the needle from the suture attachment point. Grasping the needle from its penetration point could impair seriously its penetration performance or even cause a needle fracture. An attempt to reshape a bended needle could cause loss of its strength and its resistance to breaking. Used needles should be discarded in special waste containers.

Adverse Reactions
The use of RAPID PGLA suture in some patients may cause temporary local irritation or temporary inflammatory response to foreign body or even erythema and induration in cases of subcuticular applications. Like every other foreign body RAPID PGLA suture may potentiate an existing infection.

Sterilization
RAPID PGLA surgical sutures are sterilized with Ethylenoxide gas. The sterilization method is mentioned on each single package. It is intended to be used only once and it should be discarded if its package is damaged or opened. Unused open sutures must be discarded.
Sutures should not be re-sterilized.

Storage
Store in cool and dry place away from direct heat and moisture. Never use after expiration date.

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