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
Scatola contenente 12 sutura

Scatola contenente 24 sutura

Rx Only
Federal law restricts this device to sale by or on the order of a surgeon

Fabbricante
CAMARK SA
Industrial Park Axioupoli - 61400 Axioupoli - Greece

ENGLISH

PGLA
Synthetic absorbable violet braided multifilament surgical suture polyglactin 910 coated with Poly(glycolide-co-lactide) (30/70) + Calcium Stearate

Description
PGLA suture is a braided synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide. PGLA suture is provided coated with Poly(glycolide-co-lactide) (30/70) + calcium stearate.
PGLA suture is available undyed or dyed (violet) in a wide variety of length(45cm-150cm) – diameter(USP 7/0-USP 2) combinations, with or without needle of various sizes and types, manufactured from high quality medical grade stainless steel.
PGLA meets all requirements established by the European Pharmacopeia for sterile Synthetic Absorbable Braided Sutures and with the essential requirements of the directive 93/42/EEC.

Indications
PGLA absorbable surgical sutures are indicated for use in all types of soft tissue approximation and/or ligation, including dental and general surgeries. These devices are not indicated for use in cardiovascular surgery, microsurgery and for the application in neural tissue or ophthalmic.

Actions
PGLA Synthetic Absorbable Surgical Suture elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGLA Synthetic Absorbable Surgical Suture occurs by means of hydrolysis, where the polymer degrades to glycolic and lactic acids which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength without appreciable loss of mass. PGLA retains approximately 50% of the average E.P. tensile strength requirement at the end of the 3rd post implantation week. The absorption of the suture is essentially complete between 56-70 days.

Contraindications
This suture, being absorbable, should not be used where extended approximation of tissue is required.

Warnings / Precautions / Interactions
PGLA surgical sutures should be used only from professionals who are familiar with surgical procedures and techniques involving absorbable sutures and wound closure techniques, as the risk of wound dehiscence or wound rupture may vary with the site of application and the suture material used. Adequate knot security is achieved with the standard surgical technique of flat and square ties, with additional throws depending on the surgical circumstances and the experience of the surgeon.
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
Special care should be taken in handling surgical needles. The needle holder should grasp the needle in an area which should not exceed the 1/3 to 1/2 of the total length of the needle from the suture attachment point. Grasping the needle from its penetration point could impair seriously its penetration performance or even cause a needle fracture. An attempt to reshape a bended needle could cause loss of its strength and its resistance to breaking. Used needles should be discarded in special waste containers.
When selecting a suture, its in vivo performance should be taken into consideration. The use of PGLA sutures may be inappropriate in elderly, malnourished and debilitated patients or patients suffering from other conditions that may delay wound healing. Furthermore, the degradation rate may vary following the tissue type (i.e. oral cavity) or concomitant therapy (i.e. radiotherapy).
As with any other absorbable surgical suture PGLA may have the temporary act of a foreign body. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or inflected wounds.
Topical skin sutures which must remain in place longer than 7 days may cause localized irritation and the external portion should be cut and snipped off or removed as indicated.

Adverse Reactions
The use of PGLA suture in some patients may cause temporary local irritation or temporary inflammatory response to foreign body or even erythema and induration in cases of subcuticular applications. Like every other foreign body PGLA suture may potentiate an existing infection.


Sterilization
PGLA surgical sutures are sterilized with Ethylenoxide gas. The sterilization method is mentioned on each single package. It is intended to be used only once and it should be discarded if its package is damaged or opened. Unused open sutures must be discarded. Sutures should not be re-sterilized.

Storage
Store in cool and dry place away from direct heat and moisture. Never use after expiration date.

Symbols


Production date


Expiry date


Lot number

Catalogue number


Ethylene oxide sterilization

Do not re-sterilize

Do not reuse


Caution


Non-pyrogenic


Manufacturer

Consult instructions for use

Rx Only
Federal law restricts this device to sale by or on the order of a surgeon

2803
CE mark and identification number of notified body

Do not use in case of damaged package

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