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ENGLISH

POLYPROPYLENE
monofilament non absorbable
sterile surgical suture

Description:
POLYPROPYLENE is a sterile non absorbable monofilament surgical suture made from Polypropylene (C3H6)n. It is dyed with Phthalocyanine Blue.
The suture is available in a wide variety of lengths and diameters (USP 7/0 – USP 2/0); the needle made from high quality medical-grade stainless steel is available in different dimensions and types. POLYPROPYLENE Surgical Sutures are manufactured according to the requirements of the US and European Pharmacopoeia (USP and EP) for sterile non absorbable sutures and to the essential requirements of the 93/42/EEC directive.

Indications:
POLYPROPYLENE Surgical Suture is indicated for use in all types of soft tissue approximation and/or ligation, including dental and general surgeries. The device is not indicated for use in cardiovascular surgery, ophthalmic surgery, microsurgery and for the application in neural tissue.
As it is relatively biologically inert its use is recommended in cases where minimal tissue reaction is desired. It is especially recommended for use in infected wounds or in wounds with high risk of later infection. POLYPROPYLENE suture is not adherent to the tissues and it can be removed easily and painlessly. Therefore it can be used as a pull out suture.
For Canadian market attention: POLYPROPYLENE Surgical Suture is intended for prolonged contact duration (≥ 24 hours to 29 days).

Contraindications:
The safety and the effectiveness of POLYPROPYLENE Surgical Suture in ophthalmic surgery, microsurgery and peripheral neural tissues have not been established.
Their use is also contraindicated in patients with known sensitivity or allergy to POLYPROYLENE fibers.

Warnings:
POLYPROPYLENE Surgical Suture must only be used by professional physicians who are familiar with surgical procedures and techniques. Users must have a good knowledge of suture handling and surgical knotting techniques. Knots must not be too tight and surgical instruments, such as needle holders and forceps, must be used correctly to prevent damage to the suture or make it ineffective. Avoid crushing or crimping the suture when handling it with forceps or needle holders as this may reduce its tensile strength.
The device is for single use only. Do not re-sterilize.

Precautions:
Misuse of this suture, like any other suture, can result in severe injury or death to the patient. As with any suture, care should be taken to prevent damage when handling. Avoid crushing or crimping the suture with surgical instruments or exposing the suture to sharp edges. In order to minimize needle damage, do not grasp it close to the assembly channel.
Knot security requires standard surgical techniques of flat and square knots with additional throws as indicated by surgical circumstances and the experience of the surgeon. When POLYPROPYLENE suture is used, the knots should be tightened by pulling each top of the suture in the opposite direction with equal force.

Caution:
Do not pull on the needle itself, but grasp the suture with the fingers or surgical instruments.
Care should be taken to avoid using a jerking motion, which could break the suture or cause separation of the suture from the needle. If a well-formed square knot is incorrectly knotted it will result in an unsecure knot. When the POLYPROPYLENE Surgical Suture is properly tensioned and formed, standard surgical knotting techniques will produce a secure knot. Special care should be taken in handling surgical needles. The needle should be grasped with the forceps in an area equal to 1/3 - 1/2 of its total length, from its point in which the thread is held as far as the tip. Grasping the needle from the opposite end area can damage the tip or even cause a needle fracture. Used needles should be safely discarded in special containers, without being handled, to avoid accidental contamination.
POLYPROPYLENE sterile surgical suture selection depends on the condition of the patient, the size of the tissue and wound, the surgical technique and the experience of the surgeon. The user must give the patient all the necessary instructions to ensure a proper post-operative recovery in relation to the specific application.

Sterility:
The device is supplied sterile (EtO). Do not use the device after the expiry date. In any case, make sure that the sterile packaging is intact and handled with care. Dot not use if the package is open or damaged. Unused open sutures must be discarded. The device is for single use only. Do not re-sterilize. Device re-use could cause infections, serious injuries or death of the patient and failure to retain the technical characteristics.

Adverse Reactions:
POLYPROPYLENE sutures are not absorbed, nor is any significant change in tensile strength known to occur in vivo. The use of the suture in some patients can cause initial skin irritation followed by a minimal and temporary inflammatory reaction. As every foreign body, it can enhance an existing inflammation.
Adverse effects associated with the use of suture include wound dehiscence, failure of adequate wound support in closure sites where approximation, stretching or distension occur, minimal acute inflammatory tissue reaction.
Adverse effects associated with the use of suture, in general, include increased bacterial infectivity, calculi formation in urinary and biliary tracts when prolonged contact with saline solutions such as urine and bile occurs, and pain, edema and erythema.

Dosages and administration:
Use as directed by surgical procedure.

Storage:
Store in a cool dry place, away from direct heat and moisture. Do not use the device after the expiry date.

Packaging:
POLYPROPYLENE Surgical Suture is supplied sterile in a box of 12 or 24 individually packaged sutures. Each single suture is packaged in a specific paper pack.

Symbols

Production date
Do not reuse
CE 0051 CE mark and identification number of notified body
Expiry date
Caution
Do not use in case of damaged package
Lot number
Non-pyrogenic
Catalogue number
Manufacturer
Box containing 12 sutures
Ethylene oxide sterilization
Consult instructions for use
Box containing 24 sutures
Do not re-sterilize
Rx Only Federal law restricts this device to sale by or on the order of a surgeon

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