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Fabbricante
CAMARK SA
Industrial Park Axioupoli - 61400 Axioupoli - Greece

ENGLISH

SILK
Braided - wax treated
non absorbable sterile surgical suture

Description:
Silk is a sterile non absorbable surgical suture made from braided fibers of natural silk derived from Bombyx Mori L. The fibers are dyed black with Sulphul Black 1. The number of fibers used in the braiding process depends on the required suture diameter. Silk sutures are coated with a thin layer of wax. The suture is available in a wide variety of length – diameter combinations, with or without needle of various sizes and types, manufactured from high quality medical grade stainless steel.
Silk sterile surgical sutures are manufactured according to the requirements of the European Pharmacopoeia for sterile braided non absorbable sutures and with the essential requirements of the directive 93/42/EEC.

Indications
Braided non absorbable Silk suture is indicated for use in all types of soft tissue approximation and/or ligation, including dental and general surgeries. The device is not indicated for use in cardiovascular surgery, ophthalmic surgery, microsurgery and for the application in neural tissue.

Contraindications:
The safety and the effectiveness of the silk sutures in ophthalmic surgery, microsurgery and peripheral neural tissue have not been verified.
Silk sutures should not be used where permanent tensile strength is required. Their use is also contraindicated in patients with known sensitivity or allergy to silk fibers.

Warnings:
Silk sutures should be used only from members of experienced surgical teams. Users must be familiar with suture handling and knotting techniques. Knot security requires standard surgical technique of flat and square knots according to the experience of the surgeon and the nature of the operation. Unnecessary knot tension and suture handling with surgical instruments such as needle holders or forceps can damage the suture and weaken the suture and therefore must be avoided.
Avoid crushing or crimping the suture when handling it with forceps or needle holders as this may reduce its tensile strength. Tissue invasion: Surgical Suture can result in attachment of the suture to the tissue it penetrates in long-term use. Such attachment may make removal of the suture difficult.

Precautions
Misuse of this suture, like any other suture, can result in severe injury or death to the patient. As with any suture, care should be taken to avoid damage when handling. Avoid crushing or crimping the suture with surgical instruments or exposing the suture to sharp edges. In order to minimize needle damage, do not grasp it close to the assembly channel.
Knot security requires standard surgical techniques of flat and square knots with additional throws as indicated by surgical circumstances and the experience of the surgeon.
When silk suture is used, the knots should be tightened by pulling each top of the suture in the opposite direction with equal force.

Caution
Do not pull on the needle itself, but grasp the suture with the fingers or surgical instruments.
Care should be taken to avoid using a jerking motion, which could break the suture or cause separation of the suture from the needle. If a well-formed square knot is incorrectly knotted it will result in an unsecure knot. When the Silk Surgical Suture is properly tensioned and formed, standard surgical knotting techniques will produce a secure knot.
Special care should be taken in handling of surgical needles. The needle should be grasped with the needle holder in an area 1/3 to ½ of the length from the attachment end to the point. Grasping the needle from the opposite end area can damage the point or even cause a needle fracture. Deformed needles should not be forced to their original shape as this can cause loss of strength or even needle fracture.
Used needles should be safely discarded in special containers, without being handled, to avoid accidental contamination.
Silk sterile surgical suture selection depends on the condition of the patient, the size of the tissue and wound and the experience of the surgeon.
The user must give the patient all the necessary instructions to ensure a proper post-operative recovery in relation to the specific application.
Polyester suture is intended to be used in adult and pediatric populations, regardless of gender.

Sterility
The device is supplied sterile (EtO). Do not use the device after the expiry date. In any case, make sure that the sterile packaging is intact and handled with care. Dot not use if the package is open or damaged. Unused open sutures must be discarded. The device is for single use only. Do not re-sterilize. Device re-use could cause infections, serious injuries or death of the patient and failure to retain the technical characteristics.

Adverse reaction
The use of the Silk suture in tissues results to an initial inflammatory reaction, followed by gradual encapsulation of the suture by connective tissues. Silk sutures are not absorbed, but the protein fibers are progressively degraded with gradual loss of their tensile strength.
The use of the suture in some patients can cause allergic response or temporary local irritation which is followed by temporary inflammatory reaction.
As every foreign body it can enhance an existing infection.
Dehiscence of the suture upon removal may occur in case of improper use.
In case of infected wounds acceptable surgical practices should be followed.

Dosages and administration:
Use as required by surgical procedure.

Storage:
Store in a cool dry place, away from direct heat and moisture.
Do not use the device after the expiry date.

Packaging:
Silk Surgical Suture is supplied sterile in a box of 12 or 24 individually packaged sutures. Each single suture is packaged in a specific paper pack.

Symbols

Production date
Do not reuse
CE 0051 CE mark and identification number of notified body
Expiry date
Caution
Lot number
Non-pyrogenic
Do not use in case of damaged package
Catalogue number
Manufacturer
Box containing 12 sutures
STERILE EO Ethylene oxide sterilization
Consult instructions for use
Box containing 24 sutures
Do not re-sterilize
Rx Only Federal law restricts this device to sale by or on the order of a surgeon

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