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EN Dental etchant

BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY
TO BE USED BY DENTIST ONLY

COMPOSITION: o-phosphoric acid 36%, gel base, dyestuff.

INTENDED USE
BLUE ETCH product is used for etching dentine and enamel prior to composites, fissures and cavity seals.

PREPARATION OF THE PRODUCT BEFORE USE (concerning: variant 10 ml and 50 ml)
1. Unscrew the white cap from the syringe (50 ml or 10 ml) containing the etchant.
2. Screw the connector in place of the white cap.
3. Take a small syringe (dispenser) in your hand, pull back the plunger in it and push it all the way back in.
4. Unscrew the blue plug from the connector and then screw the prepared dispenser in its place.
5. Press the plunger of the syringe with the product and pump the required amount of gel, but no more than 2 ml, into the dispenser.
6. If it is difficult to pump the product into the dispenser, move the plunger of the syringe with the product several times. The preparation will dilute due to its thixotropic properties.
7. Unscrew the filled dispenser and screw the metal applicator onto its tip. The product is ready for use.
8. Leave the connector attached to the multi-unit packaging and secure it with the blue cap.

INSTRUCTIONS FOR USE
Before the treatment, immerse the supplied applicator in ethyl alcohol, allow it to dry an then install the applicator on the syringe. Apply a thin layer of the etchant to a dried tooth surface. Etch no longer than 60 seconds. Rinse with strong air-flow and water (for about 15-20 seconds). Do not dry too much.
In case of enamel impurities repeat the process.
Product intended for multiple use except the applicator which is for single use only. Multiple use of applicator may cause the risk of reinfection.

CONTRAINDICATIONS
Do not use the product for patients with hypersensitivity to the product's components.

PRECAUTIONS
The use of dental etchants may involve with post-operative hypersensitivity, resulting in temporary pain. The BLUE ETCH product has cytotoxic potential.
Causes severe skin burns and eye damage. In case of contact rinse immediately with copious amount of water and seek medical advice.
In case of consumption do not extract emesis. Drink copious amount of water. Seek medical advice.
Use a cofferdam, disposable gloves, safety goggles/visors and protective clothing when handling the product. Store out of reach of children.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Store in original packaging at temperature below 25°C. Expiry date is placed on the main container. Once container is open, the expiry date doesn't change if it is securely closed after each use.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed or returned to the manufacturer.

PACKAGING
Syringe of 2 ml of preparation, disposable applicator.
Syringe of 10 ml or 50 ml with preparation, disposable applicators, connector, dispenser.

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