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EN Gel for root canals preparation

PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE
TO BE USED BY DENTIST ONLY

COMPOSITION: disodium edetate (EDTA 17%), gel base.

INTENDED USE
ENDO-PREP GEL is intended for chemo-mechanical root canal preparation. By flushing out calcium and magnesium ions, disodium edetate (EDTA) in the formulation softens the superficial dentin layer making it easier to remove, facilitating decongestion and shaping the canals. Lubricating substances in the product facilitate the introduction of the tool into the duct and minimize the possibility of breaking the tool in the canal.

INSTRUCTIONS FOR USE
Before the treatment, immerse the supplied applicator in ethyl alcohol, allow it to dry an then install the applicator on the syringe.
Open the tooth cavity and rinse it with sodium hypochlorite. Insert the preparation into tooth cavity using the attached applicator or place it directly on the root canal tool. Prepare the canal mechanically for approx. 1 min, while simultaneously rinsing it with sodium hypochlorite solution. Rinse the canals and cavity thoroughly with purified water.
A reusable product with a single-use applicator. Re-using the applicator may cause the risk of reinfection.

CONTRAINDICATIONS
Do not use for patients who are hypersensitive to the product components. Do not leave product in the canal or the tooth cavity between visits. This can lead to excessive softening of dentin.

PRECAUTIONS
Caution is advised when preparing the canal close to the apex area, as the tool may cause perforation of the canal. The ENDO-PREP GEL product has cytotoxic potential. EDTA can adversely affect the bonding strength of dentin to tooth restoration materials. Endodontic treatment may involve temporary postoperative pain. Causes severe skin burns and eye damage. In case of contact, rinse immediately with plenty of water and seek medical advice. If swallowed, do not induce vomiting. Drink plenty of water. Seek medical advice.
Use a cofferdam, disposable gloves, safety goggles/visors and protective clothing when handling the product. Keep out of reach of children.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Keep in original packaging below 25°C.
Expiry date is placed on the main container. Once container is open, the expiry date does not change if the container is securely closed after each use.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
Syringe containing 2 ml, 5 ml or 10 ml of the preparation, disposable applicators.

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