PL Olejek do zarabiania past

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EN Oil for making pastes

BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY
TO BE USED BY DENTIST ONLY

COMPOSITION: eugenol oil 100%

INTENDED USE
EUGENOL is used as a liquid for making pastes for root canal filling and temporary filling of cavities. It is used with zinc oxide.

INSTRUCTIONS FOR USE
EUGENOL and zinc oxide should be combined in proportions allowing obtaining a paste of the desired consistency. Stir vigorously for 30-50 seconds. It is recommended to mix one spoon of zinc oxide (0,1 g) with one drop of EUGENOL (0.03 g) on a glass plate using a spatula.
Measure the length of the canal, when filling root canals, introduce the paste into the canal, remove excess paste, close cavity tightly.
Setting time: approx. 7 min.
After application, an X-ray should be done to check that the paste adequately filled the canal and did not go beyond the apex.
Leave the filling in the cavity or in the root canal for a period of 1 week but no longer than 30 days.

CONTRAINDICATIONS
Do not use the product for patients with hypersensitivity to the product's components.

WARNINGS AND PRECAUTIONS
It can cause allergic skin reaction.
It causes eye irritation.
Endodontic treatment may cause temporary postoperative pain. The use of temporary materials may be associated with bacterial microleakage.
ZOE paste binds in a moist environment, do not over-dry the root canal / cavity before applying the paste. Do not use ZOE paste under composites, as eugenol can interfere with the polymerisation process. ZOE paste has cytotoxic potential. Care should be taken and the paste should be applied in controlled movements to avoid over-filling and accidental pushing beyond the root apex.
Use cofferdam, disposable gloves, safety goggles/visors and protective clothing when handling the product. In case of contact, rinse immediately with plenty of water. If swallowed, rinse mouth with water. Do not induce vomiting. Avoid creating and inhaling dust. In case of worrying symptoms consult a doctor. Keep out of reach of children.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Store in original packaging at a temperature below 25°C.
Protect from sunlight.
Expiry date is placed on the main container. Once the container is open, the expiry date doesn't change if it the container is securely closed after each use.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
Bottle of 10 ml or 20 ml of preparation with a dropper.

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