Adhesive AC

Gebrauchsinformation DE

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Instruction for use EN

Definition
Coltène Adhesive AC is an autoclavable sprayable adhesive for the permanent bond between addition-polymerised silicones and plastic impression trays.

Material type
Polydimethylsiloxanecopolymer

Indication
Coltène disposable impression tray:
- President Tray AC (autoclavable)
- President Tray (not autoclavable)

Also suitable:
Autoclavable or not autoclavable plastic impression trays made of:
- Acrylates, polystyrene, polycarbonates

Important
Always wear gloves.
Because of the irreversible bond between addition-polymerized silicones and plastic impression trays, the plastic impression trays sprayed with Coltène Adhesive AC after the impression is taken can no longer be used.

Contra-indications
The adhesive is not compatible for use with metal trays or special plastic impression trays, in case of doubt test the acceptability.

Side effects / Interactions
- Autoclaving or not autoclaving has no negative influence on the adhesion values.
- Cloudiness or white discoloration of the tray surface are possible, but will not affect the end result.
- Contains ethyl acetate
- Use only in a well-ventilated areas!
- Keep away from open flames!

Contains ethyl acetate

GHS02 GHS07

Danger
- H225 Highly flammable liquid and vapour.
- H319 Causes serious eye irritation.
- H336 May cause drowsiness or dizziness.
- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313 If eye irritation persists: Get medical advice / attention.

Application
Spray (2 – 3 times) all plastic impression trays from a short distance with a thin, even layer of Coltène Adhesive AC, and let it dry for at least 1 min, or dry it using compressed air.

Cleaning of trays
Disposable plastic impression trays are not intended to be cleaned.

Shelf life and storage
Coltène Adhesive AC fulfils the intended purpose at least until the expiry date, which is displayed on the containers, in a well sealed bottle at a temperature of 15 – 23 °C / 59 – 73 °F and 50 % relative air humidity.

Marking
The expiration date and LOT number are shown on the containers.

Note
Only supplied to dentists or upon their instructions. Keep out of the reach of children.

Caution
Federal law restricts this device to sale by or on the order of a dentist.

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Coltène/Whaledent AG