BRILLIANT EverGlow® DE

Gebrauchsinformation

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Instructions for use — EN

Definition
BRILLIANT EverGlow Flow is a flowable radiopaque submicron hybrid composite and available in several shades.

Composition
Methacrylates
Dental glass
Amorphous silica
Zinc oxide

Technical data
Complies with ISO 4049

Broad dimensions of inorganic filler particles:	0,02–1,5 µm
Average filler content by volume:	37 %
Inorganic filler content by weight:	60 %

Indications
BRILLIANT EverGlow Flow Universal and Translucent shades are suitable for:
- Direct class V fillings (cervical caries, root erosion, wedge-shaped defects)
- Fillings in the anterior tooth region (Class III and IV)
- Small fillings of all cavity classes (not occlusion-bearing)
- Blocking out of undercuts
- Adhesive luting of indirect composite and ceramic restorations in as far as accuracy of fit and light permeability are given
- Repairing of direct and indirect composite restorations
- Extensive resin restorations
- Cavity lining

BRILLIANT EverGlow Flow opaque shade is also suitable for:
- Aesthetic corrections (e.g. deviations in chroma)
- Masking of dark areas

Note: The masking ability is controlled by the coating thickness. To obtain optimal blending, it is recommended to keep the overall coating thickness of opaque layers low and to cover these with Universal shades.

Contra-indications
The use of hypersensitivities to ingredients of BRILLIANT EverGlow Flow. If it is not possible to keep the working area dry during application and curing. In patients with poor oral hygiene.

Side effects
Systemic side effects are not known. Isolated cases of contact allergies with products of similar composition have been described.

Interactions
Materials containing eugenol oil may interfere with the curing of composites. The use of zinc oxide eugenol cements should therefore be avoided in conjunction with BRILLIANT EverGlow Flow. Discolouration may occur with contact to cationic mouthwashes as well as plaque revealers or chlorhexidine.

Application
Cleaning and tooth isolation
Clean the tooth to be restored and the adjacent teeth. Select shade and create a dry working area.

Preparing the cavity
One should employ techniques gentle to the hard tooth substance along the principle of adhesive restoration technique. Taper the enamel margins in the anterior region, in the posterior region, only slightly taper or round the sharp enamel margins. Prepare Class V cavities with a flat taper.

Matrix/Interdental wedges
Apply a thin matrix band when working in the approximal region. Fixate the matrix appropriately with interdental wedges.

Pulp protection
Apply appropriate pulp protection in regions close to the pulp chamber according to manufacturer's instructions for use.

Application of bonding agent
Apply the bonding agent (e.g. ONE COAT 7 UNIVERSAL). Procedure according to manufacturer's instructions for use.

Application of BRILLIANT EverGlow Flow
Apply or insert directly with the application needle. BRILLIANT EverGlow Flow is available with two application needle sizes for more individuality in usage.

Polymerisation
Light cure BRILLIANT EverGlow Flow with the blue light of a polymerisation lamp (430-490 nm).

Maximum thickness of the layer and curing times:

Shade	Coating thickness	800 mW/cm² time/cm²	1600 mW/cm² time/cm²
Universal shades (Bleach, A1/B1, A2/B2, A3/B3, A3.5/B3, A4/C4	1.5 mm	20 s	10 s
Translucent shade Translucent	1.5 mm	20 s	10 s
Opaque shade Opaque A2	1 mm	20 s	10 s

Note: BRILLIANT EverGlow Flow must be cured layer by layer. The inhibition layer may not be removed for bonding. Insufficient exposure bears the risk of incomplete polymerisation.

Contra-indications
The use of BRILLIANT EverGlow Flow is light-sensitive and should therefore be protected from daylight and bright light, in particular surgical lighting and/or sunlight, for too long prior to polymerisation.

Finishing
Finishing diamonds from coarse to fine (40 µm and 15 µm) in a constant wiping motion and ample water jet.
COLTENE composite polishers are recommended for best polishing results.

Emergency measures
In case of direct contact with the oral mucosa, simply rinse with water. In the case with contact with the eyes, rinse thoroughly with water (10 min) and consult a ophthalmologist.

Notes
Only supplied to dentists and dental technicians or upon their instruction. Keep out of the reach of children! Do not use the product beyond its expiry date.
When using direct restorations, application needles are only intended for single use on patients for hygienic reasons. Close syringes after use. To avoid continued flow, pull back the plunger.

Shelf life and labelling
The expiry date and LOT number are found on the containers and external packaging.

Storage
Storage: 4–23 °C / 39-73 °F.
Do not expose to direct sunlight or sources of heat.

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