Fill-Up!™

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Kontra-indications
Known hypersensitivities to ingredients contained in Fill-Up! or in the mentioned complementary products. Inadequate oral opening, When the visibility of the working field is not possible.

Side effects
Ingredients of Fill-Up! or of the mentioned complementary products can cause sensitivity in predisposed persons.

Interactions with other agents
Use of other substances (e.g. zinc oxide eugenol) which inhibit polymerization must not be allowed to come into contact with the etched and cleaned enamel.
Contact with saliva into contact with the etched enamel and dentin or other liquids will impair the adhesive system's bonding effect.

Clinically measured times
Polymerization begins upon first contact between base and catalyst.

	Intraoral 37 °C / 99 °F
Working time	30 sec.
Setting time (incl. working time)	1 1 40 min.

Working time is dependent on the temperature. This decreases at higher temperatures and is prolonged at lower temperatures. Fill-Up! is light curing and should not be exposed to intense light for longer than 30 s prior to curing, particularly to sunlight.

Application
Cavity Bond Coat
Clean the tooth to be restored and the mesial and distal adjacent teeth with a brush and fluoride-free cleaning paste.

Tooth Isolation
A dry working area is the prerequisite for outstanding results. The use of a rubber dam is recommended (e.g. HYGENIC® or RONDO Dental Dam).

Preparing the cavity
When preparing the cavity one should employ techniques gentle to the tooth. Use tooth and stone wheels with good cutting efficiency (principle of selective restoration technique).
Use all common preparations for working on enamel and dentine and finish with 25 µm diamonds. Bevelling of the enamel margin is recommended to increase the bonding area between the tooth and the filling material, to prevent fractures of the enamel margin and thus to improve the marginal fit.

Application with ParaBond and Etchant Gel S *(see Fig. 1)*

Note: Avoid bringing Etchant Gel S into contact with gingiva or mucosa. Etchant Gel S can cause irritation. Rinse gingiva (at least 10 times) and mouth guards are important for safety reasons! Use fluoride-free enamel which has been etched and is not covered by saliva. As an alternative to etching with Etchant Gel S, the self-etching primer Fill-Up! Bond can be applied (Etchant Gel S: Leave on for 15 seconds, rinse for 30 seconds).
Important: The time between applying the adhesive and filling the cavity must not be longer than 5 min. Otherwise repeat the application of Parabond Adhesive.

Application with ONE COAT 7 UNIVERSAL & One Coat 7.0 Activator *(see Fig. 2)*

Note: *In case of a high enamel ratio, it is recommended to etch the enamel with phosphoric acid prior to applying the bonding agent.*

Important: Avoid for safety reasons bringing Etchant Gel S into contact with gingiva or adjacent teeth. Gloves, safety goggles (for both doctor and patient) and mouth guards are important! Use fluoride-free enamel which has been etched and is not covered by composite.

Note: *It is important for adequate healthy tooth substance to be present, to achieve a suitable ferrule effect (1.5 mm apical around the cavity).*

1. To accelerate the curing process and/ or to reduce the inhibition layer, heat the unit to room temp.

Emergency measures
If Fill-Up! or the mentioned complementary products come in contact with the eyes, irrigate immediately with water. If these product gets in the eyes, wash them out thoroughly with water and consult an ophthalmologist immediately.

Notes
Fill-Up! may be provided to dentists and dental technicians or on their behalf. Keep out of the reach of children. Do not use the product beyond the expiry date.

Close the bottle of the bonding agents with the respective lid immediately after application. Store in a cool, dry place.

Shelf life and labelling
Important: to avoid pockets, do not shake during the processing period. The surface of the cavity must be completely covered with the material.

Storage
Store Fill-Up! in the refrigerator (4 – 8 °C / 39 – 46 °F). After the first use, Fill-Up! can be stored at room temperature (maximum 23 °C / 73 °F).
Store Parabond Adhesive, ONE COAT 7 UNIVERSAL and One Coat 7.0 Activator in the refrigerator (4 – 8 °C / 39 – 46 °F). Etchant Gel S can be stored at room temperature (15 – 25 °C). Do not expose materials to direct sunlight or other heat source.

Caution
Federal law restricts this device to sale by or on the order of a dentist.

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