ONE COAT 7 UNIVERSAL

Gebrauchsinformation DE

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Hersteller/Lieferant
03-2018

Instructions for use EN

Definition
ONE COAT 7 UNIVERSAL is a universal, light-curing, one-component bonding agent for adhesive restoration technique, applicable in self-etch, selective-etch or total-etch technique.
Optionally, ONE COAT 7 UNIVERSAL can be used with ONE COAT 7.0 Activator for chemically curing materials. Please follow the instructions for use of One Coat 7.0 Activator.

Composition
Methacrylates
Photoinitiators
Ethanol
Water

Indications
ONE COAT 7 UNIVERSAL is indicated for:
- Direct filling technique with materials or compomers
- Bonding of indirect restorations (inlays, onlays, veneers and crowns) made of ceramic or composite
- Bonding of composite to ceramic, composite, metal and amalgam
- Dentin sealing

Application techniques
a) Self-Etch: Apply ONE COAT 7 UNIVERSAL directly on enamel and dentin
b) Selective-Etch: Etch Enamel first. Then apply ONE COAT 7 UNIVERSAL on enamel and dentin.
c) Total-Etch: Etch Enamel and dentin first. Then apply ONE COAT 7 UNIVERSAL on enamel and dentin.

Contraindications
Known hypersensitivities to ingredients of ONE COAT 7 UNIVERSAL.
Insufficient dry working.
If it is impossible to keep the working area dry.

Side effects
Ingredients of ONE COAT 7 UNIVERSAL can cause sensitivity in predisposed persons.

Interactions
Avoid contact of ONE COAT 7 UNIVERSAL with phenolic containing substances and other substances which inhibit polymerization (e.g. zinc oxide eugenol).

Applications

1. Direct filling technique
1.1 Preparation of the cavity
According to the principle of the adhesive restoration technique.
For pulp protection, a small amount of hard setting calcium hydroxide cements should be placed close to the pulp. Finish the enamel again with a diamond bur.
In case of uncut enamel, etch with phosphoric acid gel or clean with pumice. Then thoroughly rinse with water and gently dry the enamel.
For isolation, the use of a rubber dam (e.g. from HYGENIC/ RÓEKO) is recommended.

1.2 Conditioning
a) No conditioning required in self-etch technique.

b) For Selective-Etch technique: etch enamel only. And for
c) Total-Etch technique: etch enamel and dentin.
Only supplied to dentists and dental technicians or upon their instructions.
Keep out of the reach of children!
Only dispose of completely empty packages.

Shelf life and labelling
Expiry date and ☒☒ number can be found on the containers and on the product.

Storage
4 - 8 °C / 39 - 46 °F
Do not expose to direct sunlight or other heat sources.

Hygiene tips
For optimal hygienic conditions using ONE COAT 7 UNIVERSAL, the following instructions are recommended:
Use single dose or disposable brush for one patient only.
Close the bottle immediately after use. Clean with standard disinfectants.

Caution
Avoid raw restricts this device to sale by or on the order of a dentist.

Date of issue of instructions for use
03-2018

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