Lunos®

Lunos® Fissurenversiegelung
Opaque mit Fluorid
Lunos® Fissure Sealing
Opaque with Fluoride

CE 0483

Made in Germany by
DETAX GmbH
Carl-Zeiss-Str. 4
76275 Ettlingen
www.detax.com
Distributed by
DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen, DE
Tel.: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

DÜRR DENTAL

Z.REO.2020.2020.0

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EN Operating Instructions

Lunos® Fissure Sealing Opaque with fluoride

The fissure sealant is available in the following variant:
Opaque (Order No: CP021A0050), contains fluoride

Intended use
Fissure sealant for dental prophylaxis

Indication
The fissure sealant is used for prophylactic sealing of fissures, pits and carious lesions of deciduous and permanent teeth. Also for sealing work after completed orthodontic treatment, and for sealing damaged enamel surfaces.

Patient Target Group
Persons being treated in the context of a dental procedure.

Intended Users
Dentist

Product Description
Fissure sealant, light-curing, composite-based

Procedure
Etching process
Clean the enamel surface and completely remove any cleaning pastes (fluoride-free). If necessary, delimit the treatment area, and dry it with oil-free air. Apply etching gel to the enamel only. Allow the etching gel to act for 20 – 30 seconds. Then rinse off the etching gel with water spray for approx. 20 seconds. Blow off excess water with air (oil-free), or absorb it with cotton pellet. The enamel surface should not be finally dried, but remain moist. Avoid contact of the etched surface with saliva.

Syringe application
Remove the sealing cap of the application syringe. Place the black application needle on the application syringe and lock them in place. The needle must be slightly bent. Hold the application syringe vertically, then vent the cannula channel by carefully pressing the plunger, and fill it with material up to the application tip. The fissure sealant flows evenly under slight to medium pressure. The fissure sealant is preferably applied onto the moist enamel surface, but can be applied onto slightly dried enamel as well. Apply a thin layer of fissure sealant into the fissure and allow it to penetrate for about 20 seconds. Before light curing, check for any bubbles and remove them with a probe or cannula.

Curing
(Halogen > 500 mW/cm² / LED > 300 mW/cm²)
Cure with a halogen or LED lamp for 20 seconds per surface. Then control the fissure sealing and marginal areas using a probe. Fill up any missing areas and carefully remove excess material, e.g. with a finishing diamond. Polish the fissure sealant with a rubber polisher.

Safety information
- Only for the specified intended use by trained specialists.
- Wear personal protective equipment (protective gloves, goggles) when handling.
- If the etched surface becomes contaminated, e.g. with saliva, it must be re-etched.
- Expose to the light of the operating lamp only for a short time, as this will start polymerization and increase viscosity.
- Avoid contact with mucous membranes, eyes and skin; wear personal protective equipment (protective gloves, goggles, face mask). In case of accidental contact, rinse with sufficient running water, seek medical advice if necessary.
- Uncured material may cause irritation and sensitization upon contact with skin and mucous membranes.
- Immediately close again after use. Remove the cannula and use the original closure. Use a new application needle for the next application. Protect from light.
- Wear suitable personal protective equipment during the finishing of the cured material.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes
- DETAX shall not be held liable for any damage caused by misuse.
- Read and understand the safety data sheet!

For users and/or patients:
Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the country in which the user and/or patient resides.

Contraindication
Contains (meth)acrylates.
Ingredients may cause allergic reactions in predisposed persons. In such a case, do not use the product any further.

Disposal
Disposal of the contents/container must be carried out in accordance with the local/ regional/national and international regulations.

UK-REP: DÜRR DENTAL (Products) UK Ltd.,
14 Linnell Way - Telford Way Industrial Estate, Kettering, Northants NN16 8PS, United Kingdom, info@duerruk.com

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