Heliobond

[en] Instructions for Use

Light-curing bonding agent

Rx ONLY

CE 0123

Date information prepared:
2023-09-12 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose
- Enamel bonding in conjunction with the enamel-etch technique
- Repair of light-curing composite restorations and ceramic restorations

Patient target group
Patients with permanent teeth

Intended users / Special training
- Dentists
- No special training required.

Application
For dental use only.

Description
Light-curing, single-component bonding resin for optimizing the enamel-etch technique in combination with light-curing restorative materials.

Indications
Missing tooth structure in the enamel area of anterior and posterior teeth

Areas of application:
- As a bonding agent in the adhesive restoration of teeth
- Repair of light-cured composite restorations and ceramic restorations

Contraindications
The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use
- The product must not be used if the stipulated working technique cannot be applied.
- The product must not be used in combination with purely self-curing composite materials.
- Dual-curing materials always have to be light-cured.
- When using the product, the 3s PowerCure curing program must not be employed.
- The applicators are for single use only and must not be reused due to hygienic reasons.

Side effects
In rare cases, components of the product may lead to sensitization. In these cases, the material should no longer be used.

Interactions
- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used ≤ 14 days prior to the application of the product.
- In order to achieve sufficient adhesion to enamel, the acid-etch technique must be employed. Contamination of the etched and dried enamel surface with e.g. saliva or blood will hamper this effect.
- Phenolic substances (e.g. eugenol) may inhibit the polymerization of the product.
- Before any treatment with the product, completely rinse off all iron- and aluminium-based styptic agents with water.

Clinical benefit
- Reconstruction of chewing function
- Restoration of esthetics

Composition

– Bis-GMA	60 wt%
– Triethylene glycol dimethacrylate	40 wt%

2 Application

2.1 Isolation
Ensure adequate relative or absolute isolation.

2.2 Preparation of the enamel
Please observe the Instructions for Use of the phosphoric acid gel.
- Depending on the case, bevel the enamel surfaces to be treated using a diamond or clean them carefully with an abrasive (e.g. pumice powder and water). Do not use any oily or greasy abrasive paste.
- Then rinse and dry the enamel surfaces.
- Apply the phosphoric etching gel on the enamel surfaces to be treated and allow it to react for 15–30 s.
- Then carefully rinse with water spray and dry the enamel surfaces with a stream of dry and oil-free air.
- The etched and dried enamel surface must not be contaminated prior to the application of the product (if the surface comes into contact with saliva or blood, the etching and drying procedures must be repeated).

2.3 Dosage
Note: Protect the product from light. Exposure to light may lead to premature polymerization.
- Dosage depends on the area of application.
- Dispense the required amount of material into a mixing well.
- Pick up and apply the product using a customary disposable applicator.
- Seal the bottle tightly after each use.
- When used as a bonding agent, the product should be applied in a very thin layer. If necessary, use an air syringe to disperse the bonding agent into an optimally thin layer.

2.4 Application of the product
2.4.1 As a bonding agent in the adhesive restoration of teeth
Please observe the instructions for use of the curing light employed.
- Pre-treat the respective enamel surfaces if required and etch them (see preparation of the enamel).
- Apply a thin layer of the product onto the etched enamel surface using a disposable applicator. An optimally thin layer can be achieved by dispersing the product with a stream of air.
- If used in combination with light-curing luting composites, the product does not require separate curing. For direct restorations, separate curing of the product is required: Light-cure it for 10 s at a light intensity of between 500 to 1300 mW/cm^2.
- Apply the composite.
- Light-cure the composite.
- Finish the restoration.

2.4.2 Preparation of resin repair work
- Roughen the smooth resin surfaces.
- If necessary, apply a universal primer (see Instructions for Use).
- Apply a thin layer of the product. An optimally thin layer can be achieved by dispersing the product with a stream of air.
- Light-cure the product: Use an exposure time of 10 s at a light intensity of between 500 to 1300 mW/cm^2.

Polymerization

Please observe the Instructions for Use of the curing light employed.

- An exposure curing time of 10 s at a light intensity of between 500 to 1300 mW/cm² is sufficient for thin layers.
- The light probe should be held less than 5 mm above the product surface and should not touch unpolymerized material.

Application notes

- Do not expose the product to intensive light while applying it as this will considerably reduce the working time.
- Do not use the bottle intraorally. For hygienic reasons, suitable application aids must be used.
- Prevent any contamination of the applicator with blood, saliva, sulcus fluid or water during application. In the event of contamination, the cavity must be rinsed again and the procedure must be repeated with a new applicator.

3 Safety information

- The product is an irritant. Contact with the skin, mucous membrane or eyes should be avoided.
- In case of skin contact, immediately wash with copious amounts of water.
- If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician or ophthalmologist.
- In rare cases, contact with the skin may lead to sensitization towards the ingredients.
- Customary medical gloves do not provide protection against the sensitizing effects of methacrylates.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at https://ec.europa.eu/tools/eudamed.
- Basic UDI-DI: 76152082AB0ND006HU

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Warnings for phosphoric acid gel: Phosphoric acid gel is corrosive. Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and operators). If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician or ophthalmologist. Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.
The following clinical residual risks are known:

- Failure of the adhesive bond
- Postoperative sensitivity

4 Shelf life and storage

- Close bottles immediately after use; exposure to light causes premature polymerization.
- Storage temperature 2–28 °C (36–82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
- Expiry date: See note on packaging.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!
Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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