Te-Econom Plus®

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

Rx ONLY

CE 0123

Date information prepared:
2024-01-08 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose
Direct restorations of anterior and posterior teeth

Patient target group
Patients with permanent teeth

Intended users / Special training
- Dentists
- No special training required.

Use
For dental use only.

Description
Te-Econom Plus® is a light-curing, radiopaque composite (290% Al) for the direct restorative treatment of anterior and posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1). Te-Econom Plus is also suitable for restoring occlusal surfaces.
100% aluminum has a radiopacity equivalent to that of dentin and 200% aluminum is equivalent to enamel.
Te-Econom Plus cures with light in the wavelength range of 400-500 nm.

Indications
Missing tooth structure in anterior and posterior teeth

Types of restorations:
Restorations of permanent teeth

Contraindications
The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use
- If a dry working field cannot be established, the product should not be used.
- If the stipulated working technique cannot be applied, the product should not be used.
- The product should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.

Side effects
In rare cases, components of Te-Econom Plus may lead to sensitization. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions
Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylatebased materials. Consequently, the application of such materials together with Te-Econom Plus must be avoided.
Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.
Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit
- Reconstruction of chewing function
- Restoration of esthetics

Composition
Barium glass, ytterbiumtrifluoride, Bis-GMA, UDMA, barium aluminium fluorosilicate glass, TEGDMA, Silica oxide, Si-Zr mixed oxide

Total content of inorganic fillers: 59 – 61 vol%
Particle size of inorganic fillers: between 0.03 µm and 16.3 µm.

2 Application

I. Shade determination
Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the Te-Econom Plus shade guide).

II. Isolation
Adequate relative or absolute isolation is required.

III. Cavity preparation
The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp, internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel the enamel margins in the anterior region. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, 25-40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

IV. Pulp protection / Base
Do not apply a base material when using an enamel/dentin bonding agent. In very deep cavities only, areas close to the pulp should be selectively coated with a calcium hydroxide liner. Subsequently cover with a pressure resistant cement. Do not cover the remaining cavity walls so that they can be used to generate a bond with an enameldentin adhesive.

V. Placement of matrix / interdental wedge
Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix band and secure it with wedges.

VI. Conditioning / Application of the bonding agent
Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

VII. Application of Te-Econom Plus
- For optimum results, apply Te-Econom Plus in layers of max. 2 mm thickness and adapt it with a suitable instrument.
- Adapt the material correctly in order to ensure an intimate contact of the composite resin with the cavity walls.
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time (Exposure time) per increment and light intensity (Light intensity) see table 1 (Table 1).
- **The instructions for use of the curing light must be observed.**
- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.

- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.
- Optionally, a flowable composite can be applied as an initial layer.
- Cure this layer separately according to the respective instructions for use.

VIII. Finishing / Checking the occlusion / Polishing

Remove any excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- Te-Econom Plus should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
- Syringes should not be disinfected with oxidizing disinfection agents.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/ Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at https://ec.europa.eu/tools/eudamed.
- Basic UDI-DI: 76152082AFILL004JG

Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG www.ivoclar.com).
- Unpolymerized Te-Econom Plus should not come in contact with skin, mucous membrane and eyes.
- Unpolymerized Te-Econom Plus may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.
The following known residual clinical risks exist:

- Failure of the adhesive bond (loss of the restoration)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2 - 28 °C (36 - 82 °F)
- Close syringes immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Date of expiration: see information on syringes and packages.

Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

Table 1

Light intensity	Exposure time
500 – 900 mW/cm²	20 s
1,000 – 1,300 mW/cm²	10 s

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- A se respecta Fișa cu date de securitate (SDS) (disponibilă la Ivoclar Vivadent AG www.ivoclar.com).