Kerr™

# Herculite™ XRV Ultra™

nanohybrid composite restorative

**Composite Material**

# Instructions For Use

# HERCULITE XRV ULTRA

Nanohybrid Composite Restorative

Nanohybrid composite Herculite XRV ULTRA is a light-cured, resin-based, dental restorative designed for direct placement. It is indicated for all cavity classes in anterior and posterior teeth.
The new generation of Herculite XRV Ultra employs state-of-the-art technology to offer improved aesthetics, handling and mechanical properties for you and your patient's convenience. Herculite XRV Ultra offers superior aesthetics as it is easy to polish, retains clinical luster and has a lifelike appearance replicating fluorescence and opalescence of the natural tooth. We have taken extraordinary steps to maintain the mechanical properties, fracture toughness and wear resistance, giving you the durability and reliability you have come to expect from Herculite brand.
Storage Conditions and Shelf Life
Herculite XRV Ultra, should be stored at ambient temperature.
The expiration date is marked on packaging. Do not use after expiration date.
Every effort should be made to protect the composite from visible light prior to use.

**PLEASE NOTE THE FOLLOWING:**
You have purchased this composite in either a bulk syringe or Unidose delivery. Either delivery has certain handling and storage requirements that if followed will provide excellent performance for the user over the life of the material. Please follow the appropriate recommendations indicated below for the delivery choice that you have made.
Syringe Delivery: To ensure optimal performance of syringes during use: After dispensing the composite, relieve pressure by turning the syringe handle back two turns. This will prevent material build-up in the cap and provide reliable composite use after storage.
Unidose Delivery: Unidose™ tips are designed for single patient use only, to prevent cross contamination between patients. Do not re-cap and/or re-use the Unidose tip once material has been dispensed for that patient.

**English**

**PRIOR TO PLACEMENT**
**RECOMMENDATIONS ON PROPER BONDING**

- Isolate tooth from adjacent teeth. Use clear, thin strips. Rubber dam is highly recommended.
- Moist dentin surfaces (rather than dry, desiccated dentin) are shown to actually enhance the bond strengths of today's modern bonding agents (such as OptiBond Solo Plus, OptiBond FL, or OptiBond All-In-One.) NOTE: "Moist" does not infer blood or saliva contamination.
- Placement sequences such as those outlined in this instruction booklet should be followed closely.
- Take a moment to identify the surface to which you are bonding.
- Etching dentin is recommended.
- Please take care to ensure that your air line is free of oil and other contaminants.

**OPTIBOND SOLO PLUS and OPTIBOND ALL-IN-ONE**

If you have purchased Herculite XRV Ultra in kit form, you will notice that OptiBond Solo Plus is included. This is recommended for use with all Kerr composites. OptiBond Solo Plus is 15% filled with the same 0.4-micron barium glass found in the Kerr XRV Ultra composite. This unique technology provides a high level of protection against microleakage, while sustaining high-bond strengths to a variety of surfaces. The filler not only reinforces the hybrid zone, but penetrates the dentin tubules as well, creating a true "structural bond." With OptiBond Solo Plus, adhesion promoters are carried in an ethanol solvent, thus greatly reducing the tedious need for multiple coats and constant re-application commonly found with acetone adhesives. OptiBond Solo Plus is indicated for both direct and indirect (in conjunction with a luting agent) bonding applications.

OptiBond Solo Plus is available in two delivery systems, Bottle and Unidose. Kerr's Unidose delivery adhesive system was the first of its kind, designed as a progressive step toward ensuring fresh material for each bonding procedure as well as the ultimate in infection control. It is protected under U.S. Patent 5,860,806.

OptiBond Solo Plus is a "total-etch" adhesive. OptiBond All-In-One is a single-component, "self-etch" adhesive and may also be used with Herculite XRV Ultra.

**DIRECT BONDING TO ENAMEL/DENTIN WITH OPTIBOND SOLO PLUS TOTAL-ETCH TECHNIQUE**

1. Isolate and dry.
2. Place Kerr Gel Etchant* on enamel and dentin for 15 seconds.
3. Rinse thoroughly, ensuring that all etch is removed.
4. Dry lightly, do not desiccate.

5. If using Solo Plus in bottle delivery, shake bottle briefly before dispensing. If using Solo Plus in Unidose delivery, twist capsule at middle, and discard stem side.
6. Apply OptiBond Solo Plus to enamel/dentin surface with light brushing motion for 15 seconds.
7. LIGHTLY air thin for maximum 3 seconds.
8. Light cure for 20 seconds.**
9. Proceed with placement of Herculite XRV Ultra composite.

Note: Kerr Disposable Mixing Wells, Gel Etchant tips and Applicator Tips are designed for single use only, to prevent cross contamination bet ween patients.

**DIRECT BONDING TO ENAMEL/DENTIN WITH ALL-IN-ONE – SELF-ETCH TECHNIQUE**

1. Prepare cavity. Wash thoroughly with water spray and air dry (Do not desiccate).

   Note: For optimal adhesion, always prep (bevel) all enamel and/or margins prior to application of OptiBond All-In-One adhesive.

2. Shake OptiBond All-In-One adhesive bottle briefly.

   **Note:** Upon first use, shake OptiBond All-In-One adhesive bottle vigorously for 10 seconds.

   - Bottle Delivery: Dispense 2-3 drops of OptiBond All-In-One adhesive into a clean well. Replace cap immediately after dispensing.
   - Unidose Delivery: Unidose™ tips are designed for single patient use only, to prevent cross contamination between patients. Open the Unidose device and insert the applicator brush into the container to saturate the tip.

3. Using the disposable applicator brush, apply a generous amount of OptiBond All-In-One adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.
4. Apply a second application of OptiBond All-In-One adhesive with a brushing motion for 20 seconds.
5. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air.
6. Light cure for 10 seconds.*
7. Apply composite according to manufacturer's instructions for use.

**PLACEMENT OF Herculite XRV Ultra COMPOSITE**

Herculite XRV Ultra is available in a wide variety of VITA Shades. For details on all available shades and their various uses, please refer to "Explanation of Herculite XRV Ultra Shading System" below for more information.

1. Select the desired shade.
   - Syringe Delivery: After extruding material out on pad, make sure to back off the syringe at least two full turns to relieve pressure on the material.
   - Unidose Delivery: Insert tip and properly align within dispenser. Place Unidose tip at the deepest portion of the preparation.

   Increments should be no more than 2mm at a time.
2. After placing an increment, stroke the composite to ensure marginal adaptation.
3. Light Cure each increment for a minimum of 40 seconds**. When multiple surfaces are available, cure each surface for this recommended time.

**EXPLANATION OF Herculite XRV Ultra SHADING SYSTEMS**

Herculite XRV Ultra is available in a variety of VITA Shades. Due to the extensive process involved in matching VITA shades, many composite manufacturers choose to include their own shade guide with approximations of the actual VITA shade. Since Kerr composites exhibit the ideal match to VITA, simply use the VITA shade guide in your office for selection.

(Note: VITA is a registered trademark of VITA Zahnfabrik Company, Germany)

The following is a listing of the available shades:

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The enamel and dentin shades have been developed to match the opacity of enamel and dentin and should be placed to correspond with the layers of tooth structure that have been removed during preparation process.

Incisal shade offers the dentist a high level of translucency for building incisal edges, or building translucent anatomy. It is an excellent overlay material on a variety of composite restorations.

XL Shade can be useful to match bleached as well as deciduous teeth, it is lighter than the lightest VITA Shade.

Unidose Composite Gun Sterilization Guidelines – Steam Autoclave: 30 Minutes @ 270F.

* **CAUTION:** Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water. Get medical attention for eyes. Do not take internally.

** **Recommended Cure Times:** The Demi Plus, Demi, 5 seconds, L.E.Demetron II, 5 seconds; L.E.Demetron I, 10 seconds; or Optilux 501 in Boost mode, 10 seconds, Ramp Mode, 20 seconds, or Regular Mode, 20 seconds. For all other lights, see manufacturer's recommendation.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

"Optibond" is trademarks owned by Kerr Corporation. "Optibond" is registered whit the U.S. PTO and other countries.

**Limited Warranty – Limitation of Kerr's Liability**
Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.
EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Distributor: KerrHawe S.A.