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SYMBOL DESCRIPTION
Full explanation of symbols used on Kerr packaging is located at:
http://www.kerrdental.com/symbols-glossary

Caution
Manufacturer
European Conformity
Authorized representative in the European Community
Consult instructions for use
Do not re-use
Use-by-date
Batch code
Catalogue number
Medical device
Temperature storage limitations: low 2°C/36°F and high 28°C/82°F
Keep away from sunlight
Importer

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GENERAL INFORMATION
Nanohybrid composite Herculite XRV Ultra is a intra-oral, light cured, resin-based, dental restorative designed for direct placement. The new generation of Herculite XRV Ultra employs state-of-the-art technology to offer improved aesthetics, handling and mechanical properties for you and your patient's convenience. Herculite XRV Ultra offers superior aesthetics as it is easy to polish, retains clinical luster and has a lifelike appearance replicating fluorescence and opalescence of the natural tooth. We have taken extraordinary steps to maintain the mechanical properties, fracture toughness and wear resistance, giving you the durability and reliability, you have come to expect from the Herculite brand. Product is available in a variety of VITA shades. Herculite XRV Ultra is available in Unidose™ compules and syringe delivery system suitable for occlusal surfaces.

PRODUCT COMPOSITION
The term nanohybrid is used to describe the filler system with composition of both nano- sized (particles within the range of 1 nm – 1.000 nm atoms) with other micron-sized fillers. 100% aluminium has a radio-opacity equivalent to that of dentine and 200% aluminium is equivalent to enamel.
All shades except bleach shade are radiopaque with about 250% aluminium.
The fillers comprise pre-polymerized (PPF)fillers, barium glass filler, silica, and ytterbium trifluoride(about 78.5% by weight or 60.5% by volume for the dentin and enamel shades, about 72.0% by weight or 52.0% by volume for the bleach shade). The total inorganic filler is up to 78% by weight or 59% by volume. All shades have a primary particle size range of 40 nm to 30 microns.
Product contains Bis-GMA, bis-GMA, TEGDMA, and/or UDMA resins.
None of the ingredients in the product contain pharmacologically active components.

INDICATIONS FOR USE
Herculite XRV Ultra is a dental composite restorative material intended to be used in all classes of cavities.

CONTRAINDICATIONS
None known.

WARNINGS
1. Do not use this product on patients with known methacrylates allergies.
2. Avoid product application in direct contact with pulp tissue since exposure to this product may injure the pulp.
3. Avoid skin, eye, and soft tissue exposure to this product since uncured methacrylate may cause an irritation reaction. Wash thoroughly with water after contact. Get medical attention for eyes. Do not take internally.
4. Opti*Bond Gel Etchant contains 37.5% phosphoric acid. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact. Get medical attention for eyes. Do not take internally.
5. Do not use for direct pulp capping procedure.

PRECAUTIONS
1. This product is intended to be used by qualified dental professionals on a general patient population.
2. Wear protective protective equipment when using and handling this product.
3. Use of protective gloves and a no-touch clinical technique is strongly recommended.
4. If the composite is touched with commonly used gloves, acrylate may penetrate the gloves. Dispose of gloves and replace with fresh gloves.
5. Suitable protective eyewear for patients and medical staff is recommended when light curing the material.
6. Use a rubber dam to isolate the treated tooth and avoid fluid contamination during material placement.
7. Protect the Herculite XRV Ultra material from ambient light that may initiate polymerization.
8. Once the Opti*Bond unit dose is mixed, opaque dispensing tips, applicators or protect from ambient light.
9. Replace the syringe cap after each use to prevent the resin from hardening in the syringe.
10. Opti*Bond Solo Plus bottle should be tightly closed immediately after use.
11. Unidose tips, Opti*Bond syringes tips, disposable mixing wells, and Kerr Applicators are intended for single patient use only, to prevent cross contamination between patients.
12. Do not use this product after the expiration date.

Precautions regarding use of Opti*Bond Gel Etchant:
1. Never use excessive force to syringe gel into dispenser to dispense etchant which may lead to over dispensing of the material. If resistance is felt replace current dispensing tip with a new tip and check flow outside the oral cavity.
2. Replace the syringe cap after each use to prevent the product from drying in the syringe.

ADVERSE EVENTS
If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS
1. Clean the tooth surfaces and water to remove any debris or stain that may interfere with the shade selection process.
2. Select the desired shade. A composite mock-up or VITA-matched shade guide may be used to check the color of the tooth to be restored and any adjacent teeth. In many cases, a single shade of Herculite Ultra is required for the full restoration or to restore the tooth. The shade of the restorative material is affected by the thickness of the composite. Shade selection and mock-ups should match the thickness of the tooth to be restored.
3. After shade selection, isolate the tooth with your preferred isolation method.
4. Prepare the tooth using your preferred rotary instrument.
5. Complete the bonding procedure. Herculite XRV Ultra is compatible with all methacrylate dental adhesive systems. However, Opti*Bond FL and Opti*Bond Solo Plus.
6. If desired, place a liner with a low viscosity flowable composite, such as Herculite Ultra Flow and light cure according to the manufacturer's instructions.
7. Herculite XRV Ultra features optimized, non-sticky handling. It permits the use of hand instruments to define most of the restoration anatomy, while helping to minimize the time and effort associated with rotary instrumentation.
8. For Unidose delivery: insert Unidose tip into any composite dispenser. Unit-dose should be firmly seated in a dispenser. Remove the cap and inject composite material from the Unidose tip directly into the cavity, starting with the deepest portion.
For Syringe delivery: After extruding material onto the pad, back off syringe at least 2 full turns to relieve pressure on the material. Recap the syringe immediately after use.
9. Apply composite layers up to 2 mm in thickness for enamel fracture(not cuspal) shades or 1.5 mm in thickness for dentin shades, taking care to properly adapt the material to the cavity surface. Avoid exposure of Herculite Ultra to intense operatory light during shaping of the material.
10. Cure each layer. Refer to the table at the end of these instructions for curing recommendations.
11. After curing, complete the finishing and polishing procedure. Adjust occlusion as needed.

***RECOMMENDED CURING TIMES**
Herculite XRV Ultra is intended to be cured by an LED light with a minimum light intensity of 1,000 mW/cm² and a light wavelength profile in the 400–520 nm range. Place refer to the table below for form curing recommendations. For other brands of curing lights, follow their manufacturing accommodations. Use polymer restoration, Kerr recommends additional 10 second cures on both the buccal and lingual surfaces of the tooth after removal of the matrix band.

Curing Time	Dentin Shades	Enamel/Translucent (Incisal) Shades
Incremental Depth	1.5 mm	2.0 mm
Demi™ Plus (or light with output > 1,000 mW/cm²)	10 seconds	10 seconds

CLEANING AND DISINFECTION
Unidose Tips: Do not re-use Unidose tips.
Unidose tip guns: Refer to the reprocessing instructions for the manufacturer of the tip gun you are using. Composite syringes: To prevent composite from getting into the plunger end, add a little or ether body fluids, use protective barriers. Always handle the syringe with clean gloves. Do not reuse if it becomes contaminated.

STORAGE AND DISPOSAL
Product should be stored between 2°C to 28°C (36°F to 82°F). Avoid excess heat. If product is stored in a refrigerator allow to warm to room temperature prior to use. Every effort should be made to protect the composite from visible light prior to use.
The expiration date is marked on packaging. Do not use restorative after expiration date.
For proper disposal always follow Federal or State laws applicable to that country.

LIMITED WARRANTY – LIMITATION OF KERR'S LIABILITY
Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect of material or workmanship, Kerr'sliability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages. EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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