**PRINT SPECIFICATION PAGE ONLY.**
**DO NOT PRINT THIS PAGE IN FINAL IFU.**

| Instructions for Use – Optibond FL | | |
|---|---|---|
| **PART NUMBER** | 41204 | |
| **REVISION** | 01-'22 | |
| **REVISION DATE** | 2022-01 | |
| **BINDING** | 1 glue point | |
| **CORNER** | N/A | |
| **COLORS** | 1/1 Black | |
| **SIZE** | Flat: 688 mm (w) x 400 mm (h) / Folded: 70 mm (w) x 45 mm (h)<br>Tolerance Open: +/- 10mm / Closed width: 45-47 mm Length 68-70 mm.<br>Thickness 8-10 mm (measured without any pressure) | |
| **PAGE COUNT** | 3 pages front and back / 2 pieces | |
| **MATERIAL** | Soramedico 40 gr/mg | |
| **ADHESIVE MATERIAL** | 2 Glue dots + glued together with 1 glue dot | |
| **COATING** | N/A | |
| **DIELINE** | Blueline Does Not Print | |
| **LANGUAGES: 27** | en-ENGLISH | ENGLISH |
| | es-SPANISH | ESPAÑOL |
| | fi-FINNISH | SUOMI |
| | fr-FRENCH | FRANÇAIS |
| | tr-TURKISH | TÜRKÇE |
| | de-GERMAN | DEUTSCH |
| | hu-HUNGARIAN | MAGYAR |
| | ro-ROMANIAN | ROMÂNÂ |
| | nl-DUTCH | NEDERLANDS |
| | no-NORWEGIAN | NORSK |
| | pl-POLISH | POLSKI |
| | pt-PORTUGUESE | PORTUGUÊS |
| | it-ITALIAN | ITALIANO |
| | cs-CZECH | ČESKÝ |
| | sv-SWEDISH | SVENSKA |
| | et-ESTONIAN | EESTI |
| | da-DANISH | DANSK |
| | bg-BULGARIAN | БЪЛГАРСКИ |
| | el-GREEK | ΕΛΛΗΝΙΚΑ |
| | sk-SLOVAK | SLOVENSKÝ |
| | sl-SLOVENIAN | SLOVENSKÍ |
| | ru-RUSSIAN | РУССКИЙ |
| | lt-LITHUANIAN | LIETUVIŲ |
| | sr-SERBIAN | SRPSKI |
| | is-ICELANDIC | ÍSLENSKA |
| | uk-UKRAINIAN | УКРАЇНСЬКА |
| | ar-ARABIC | العربية |
| **SUPPLIER** | RK Tisk | |

Printed by THOA.NGUYEN on 11 May 2023, at 10:28:18 a.m.

**SYMBOLS DESCRIPTION**
Full explanation of symbols used on/in packaging located at:
**http://www.kerrdental.com/symbols-glossary**

REF Catalogue number
Manufacturer
LOT Batch code
European Conformity
Consult instructions for use: http://www.kerrdental.com/ifu/ebooks-center
Caution
EC REP Authorised representative in the European Community
Rx only CAUTION: U.S. federal law restricts this device to sale by or on the order of a dentist.
Keep dry
Do not reuse
Keep away from sunlight
Temperature limit(s): Entre 2°C/35°F and High 25°C/77°F
Medical device
Importer

**CONTENTS**

1. Tube (5 g)
2. Dentine™ (3.5 ml)
3. 5.4 g Intact Syringe (3 g)
4. 1.4 g Delivery Tip
5. Mixing Well
6. Applicators

optiBond™
Two-Component Total-Etch Adhesive
Instructions For Use
Kerr Italia S.r.l.
Via M. Gioia 144
Scalenghe (TO) 10067 ITA
+39 011 9966 680
Kerr Dental S.r.l.
Via Piquet, 13
Verbania (VB) 28922 ITA
+39 0323 581 581
Rx only
CE 2797
Kerr
41204_REV01-20_OPTIBOND_FL_IFU_EU_PU_72.indd 1
12/28/21 10:30 AM

## en - ENGLISH

**OPTIBOND™ FL**
Two-Component Total-Etch Adhesive

**GENERAL INFORMATION**
OptiBond FL is a two-component adhesive system indicated for direct light-cured applications. This unique technology provides the highest level of protection against sensitivities, while maintaining high bond strengths to a variety of surfaces. The filled adhesive provides a true "one-step bond" with hands-on utilization of two systems. With OptiBond FL, ultimate preparations are carried out as a final cleaning step, eliminating the need for multiple acids.
OptiBond FL is available in both bottle and tube dispensers.

**PRODUCT COMPOSITION**
OptiBond FL is a two-part primer-adhesive system indicated for direct light-cured applications. OptiBond FL Adhesive is 45% filler, with 0.6 micron fluoride glass and co-polymers.

**INDICATIONS FOR USE**
OptiBond FL is a multi-purpose bonding agent designed to be used in the following situations: composite to tooth structure; composite repair; porcelain repair; porcelain repair; composite to metal; bonding composite over build-up materials.

**CONTRAINDICATIONS**
None known.

**WARNINGS**
1. Do not use this product on patients with known methacrylate allergies.
2. Avoid prolonged skin contact or direct contact with pulp tissue since exposure to this product may cause irritation.
3. Avoid contact with the eyes. In case of contact, immediately irrigate with copious amounts of water for 15 minutes. Consult a physician if irritation persists.
4. Use rubber dam to isolate the treated tooth and avoid fluid contamination during material placement.
5. Contaminated tooth should be isolated immediately after use.
6. OptiBond FL is sensitive to light. Protect from ambient light. Dispensing tips and tubes are light-cured by ambient light. Store in amber bottle.
7. Keep bottles, dispensing tips and tubes tightly capped when not in use. Do not use if bottle, dispensing tips and tubes are open for an extended period of time. Do not use if product has thickened.
8. Do not use polishers that create excessive heat.
9. Do not use applicators from the OptiBond FL PRIME kit as a separate applicators for the OptiBond FL kit.

**PRECAUTIONS**
1. This product is intended to be used by qualified dental professionals only as general purpose adhesive.
2. Wear universal protective equipment when using and handling this product.
3. Submit adverse reactions to this product and medical audit to recommended when light curing the material.
4. Use a rubber dam to isolate the treated tooth and avoid fluid contamination during material placement.
5. Contaminated tooth should be isolated immediately after use.
6. OptiBond FL is sensitive to light. Protect from ambient light. Dispensing tips and tubes are light-cured by ambient light. Store in amber bottle.
7. Keep bottles, dispensing tips and tubes tightly capped when not in use. Do not use if bottle, dispensing tips and tubes are open for an extended period of time. Do not use if product has thickened.
8. Do not use polishers that create excessive heat.
9. Do not use applicators from the OptiBond FL PRIME kit as a separate applicators for the OptiBond FL kit.

**Precautions in case of exposure of the fluid to the eyes:**
1. Never use excessive force or pressure when using a syringe or when a valve must be worn during the use of the material. If residue or fluid enters the eyes, immediately irrigate with water for 15 minutes. Consult a physician if irritation persists.
2. Replace the syringe cap after each use to prevent the product from drying in the syringe.

**ADVERSE EVENTS**
If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority in the country in which the user and/or patient is established.

**STEP-BY-STEP INSTRUCTIONS**
1. **Product Preparation**
   Bottle Delivery: Shake OptiBond FL Activator (B) bottle briskly (5 seconds) prior to dispensing. Dispense "X" drop of OptiBond FL Activator (B) and "X" drop of OptiBond FL Primer (A) into a clean, dry dish. Begin scrubbing immediately after dispensing. Two applications may be completed. Dip applicator tip into OptiBond FL Activator (B) bottle between applications. Discard dish after the second application.
   Tube Delivery: See Tube Use Directions for Tube Delivery System.
2. **Tooth Preparation and Etching**
   - Tooth Preparation: For a rubber dam, tie-in before the treated tooth. Premix gluma desensitizer and pumice to form a paste. Apply to tooth surface using gauze. Rinse thoroughly with water. Dry cavity.
   - Etching: Place Ivory Gel Etchant on enamel and dentine for 15 seconds. Rinse thoroughly with water. Blot dry with cotton pellet or cotton roll. Gently air dry with 5 sec. stream of air without overdrying the dentine.
3. **Product Application**
   1. Apply OptiBond FL PRIME bottle or tube brush tip applicator to tooth surface with firm pressure. Scrub for 15 seconds. Gently air thin with 5 sec. stream of air. Do not overdry. Overdrying may cause post-operative sensitivity. Re-wet with OptiBond FL PRIME if overdried.
   2. Apply OptiBond FL ADHESIVE (Optional): Dispense 1 drop each of OptiBond FL Activator (B) and OptiBond FL Adhesive (C) into a clean, dry dish. Mix with applicator for 5 seconds. Apply to tooth surface with firm pressure. Scrub for 15 seconds. Gently air thin with 5 sec. stream of air. Do not overdry. Overdrying may cause post-operative sensitivity. Re-wet with OptiBond FL PRIME if overdried.
   3. Light-cure for 10 seconds. Refer to the table at the end of this section for recommended curing times.
   4. Proceed with composite placement per manufacturer's instructions for use.

***Recommended Curing Times**
Actual curing times are dependent upon a variety of factors to ensure they are functioning and providing the expected irradiance. The OptiBond FL is intended to be cured by a halogen or LED light unit with a minimum output of 400 mW/cm² and a wavelength output within the 450–500 nm range. Please refer to the table below for curing recommendations.

| Light Type | Curing Time |
|---|---|
| Spektor™ High-Intensity LED Light (800 – 1000 mW/cm²) | 10 seconds |
| Demi™ Plus (as light with output > 2000 mW/cm²) | 5 Seconds |

For all other lights, see manufacturer's recommendations.

**CLEANING AND DISINFECTION**
To prevent cross-contamination between patients, use clean gloves when handling acids, bottles, dispensing tubes, mixing well and liner applicators.
Bottles: Do not reuse OptiBond FL bottles once contaminated.
Rubber: Do not reuse Dentine.
Tubes: Do not reuse tubes once contaminated with saliva or other body fluids, use protective barriers. Wipe tube with alcohol wipe, syringe cap and dispensing tip with clean gloves.

**STORAGE AND DISPOSAL**
Product should be stored between 2°C to 25°C (36°F to 77°F). If product is stored in a refrigerator allow the product to reach room temperature before use. Every effort should be made to protect the bonding agent from ambient light prior to use. The expiration date is marked on packaging. Do not use after expiration date. Dispose of product in accordance with local regulations. Mercury Free. The SDS for this product is available at kerrdental.com. It can also be obtained from your local reseller.

**LIMITED WARRANTY – LIMITATION OF KERR'S LIABILITY**
Kerr's entire liability and the customer's exclusive remedy, is designed to succeed only in using Kerr's product. The device's owners at this end liability for damages arising out of the improper use of Kerr's product, whether the cause is a defect in material or workmanship, or the failure to follow directions, warnings or precautions. Kerr is not responsible for any other expenses or damages. The duration of the warranty is limited to the shelf life of the product. To take advantage of this limited warranty, the defective product must be returned to Kerr or to a Kerr authorized reseller, accompanied by a written statement describing the defect. Kerr will, at its option, repair or replace the product or refund the purchase price. KERR MAKES NO WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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