SYMBOLS DESCRIPTION

	Warning or Caution
	Consult Instructions for Use
2°C (36°F) 25°C (77°F)	Temperature Storage Limitations: Low 2° C/35° F and High 8 ° C/46 ° F
	Manufacturer
	Use-by Date
	Do not re-use
	Date of manufacture
LOT	Manufacturing lot number
EC REP	European Authorized Representative
REF	Part Number

Kerr™

OptiBond™ Universal

Single-Component Universal Adhesive

Kerr Corporation
1717 West Collins Avenue, Orange, CA 92867 USA
1-800-KERR-123 | kerrdental.com

EC REP **Kerr Italia S.r.l. Via Passanti, 332, I-84018 Scafati (SA) Italy | +39-081-850-8311 | kerrdental.eu**

Distribution Facility:
Kerr Australia Pty. Ltd.
Unit 10, 112-118 Talavera Road
North Ryde NSW 2113, Australia

Rx ONLY

DMC3190_Rev.01

en - ENGLISH

OptiBond™ Universal
Single-Component Universal Adhesive

DEVICE DESCRIPTION

OptiBond™ Universal is formulated to provide reliable bond strengths in any etch technique, including self etch, selective enamel etch, and total etch. This flexibility allows the clinician to choose a technique based on either the indication or personal preference, all with a single adhesive system.

INDICATIONS FOR USE

Direct Applications
- Light-cured composite and compomer restorations.
- Composite/ceramic/metal repairs.
- Cavity sealing for amalgam restorations.
- Light-cured or dual-cured core build-ups (always light-cure the first incremental build-up of core material).

Indirect Applications
- Veneers.
- Porcelain, composite, and metal-based inlays, onlays, crowns and bridges.
- Endodontic posts.
- Cavity sealing as a pretreatment for indirect restorations.
- OptiBond Universal adhesive can be used to prime metal, ceramic, and composite substrates/restorations. This eliminates the need for other primers such as silane or metal primer.

IMPORTANT: For indirect applications, use NXG resin cement whenever there is no or limited light accessibility. NXG, either self-cured or dual-cured, is compatible with OptiBond Universal adhesive.

CONTRAINDICATIONS

1. Patients who have a history of severe allergic reaction to methacrylate resins.
2. Direct pulp capping.
3. In combination with a chemical-cure or dual-cure resin material (other than NXG) when there is insufficient light accessibility.

WARNINGS

Do not use on patients with known methacrylate allergies. Uncured methacrylate resin may cause contact dermatitis and damage to the pulp. OptiBond Universal adhesive contains acetone and is flammable. Keep away from sources of ignition. OptiBond Universal is a corrosive liquid. Gel Etchant is a corrosive acid - contains phosphoric acid. Avoid contact with skin, eyes and soft tissue. In case of contact with skin or eyes, flush immediately and thoroughly with water. Get medical attention for eyes. Do not take internally.

PRECAUTIONS

Refrigerate upon receipt between 2° and 8°C. OptiBond Universal adhesive bottle should be tightly closed immediately after use. OptiBond Universal adhesive is a light-cured material. Once the material is placed in a mixing well, use immediately or protect from ambient light.
Gel Etchant syringe: Never use excessive force on syringe plunger to dispense etchant which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.

ADVERSE REACTIONS

None Known.

INSTRUCTIONS FOR USE
Step-by-Step Instructions

Direct Restorations
Bonding Composites/Compomers to Dentin/Enamel
1. Prepare cavity. Wash thoroughly with water spray and air dry (do not desiccate).

For Self-Etch Technique, omit Steps 2-3 and proceed to Step 4.

For Selective Enamel Etch Technique
2. A - Etch enamel margins with a phosphoric acid gel etchant (e.g. Kerr's 37.5% Gel Etchant) for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to Step 3.

For Total Etch Technique
2. B - Etch enamel and dentin cavosurfaces with a phosphoric acid gel etchant (e.g. Kerr's 37.5% Gel Etchant) for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to Step 3.
3. Air dry lightly; do not desiccate.
4. Shake OptiBond Universal adhesive bottle briefly.

Note: Upon first use, shake OptiBond Universal adhesive bottle vigorously for 10 seconds.

5. **Bottle Delivery:** Dispense 2-3 drops of OptiBond Universal adhesive into a clean well. Replace cap immediately after dispensing.

Unidose™ Delivery
For single patient use only: Open the unidose device and insert the applicator brush into the container to saturate the tip.
6. Using the disposable applicator brush, apply a generous amount of OptiBond Universal adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.
7. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air. The surface should have a glossy uniform appearance. If not, repeat steps 6 and 7.
8. Light cure for 10 seconds.*

Indirect Restorations
(See indication-specific information in Section A-D below):

I. **Restoration Surface Preparation:**
Ceramic/Porcelain Restorations: Sandblast the surface with 50µ alumina with a pressure of about 30 psi (0.2 MPa) and etch the restoration with hydrofluoric acid for 1 minute. Apply a coat of OptiBond Universal adhesive† on the internal surface of the restoration, air thin with gentle air first and then strong air to avoid pooling of OptiBond Universal adhesive, and light-cure for 10 seconds* (light-curing of OptiBond Universal adhesive is optional if NXG resin cement is used for cementation).
†Alternatively, apply a coat of silane on the internal surface of the restoration, and lightly air thin. Protect the restorations from light as Kerr silane contains resin and is light-sensitive.

Metal-based, Zirconia-based, Alumina-based and Composite Restorations: Sandblast the internal surface with 50µ alumina with a pressure of about 60 psi (0.4 MPa) for metal-based, zirconia-based and alumina-based restorations, or a pressure of about 15 psi (0.1 MPa) for composite restorations. Apply a coat of OptiBond Universal adhesive on the internal surface of the restoration, air thin with gentle air first and then strong air to avoid pooling of OptiBond Universal adhesive, and light-cure for 10 seconds* (light-curing of OptiBond Universal adhesive is optional if NXG resin cement is used for cementation).

II. **Tooth Preparation and Conditioning**
1. Clean the preparation. Wash thoroughly with water spray and air dry (Do not desiccate).

For Self-Etch Technique, omit Steps 2-3 and proceed to Step 4.

For Selective Enamel Etch Technique
2. A - Etch enamel margins with a phosphoric acid gel etchant (e.g. Kerr's 37.5% Gel Etchant) for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to Step 3.

For Total Etch Technique
2. B - Etch enamel and dentin cavosurfaces with a phosphoric acid gel etchant (e.g. Kerr's 37.5% Gel Etchant) for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to Step 3.
3. Air dry lightly; do not desiccate.
4. Shake OptiBond Universal adhesive bottle briefly.

Note: Upon first use, shake OptiBond Universal adhesive bottle vigorously for 10 seconds.

5. **Bottle Delivery:** Dispense 2-3 drops of OptiBond Universal adhesive into a clean well. Replace cap immediately after dispensing.

Unidose™ Delivery
For single patient use only: Open the unidose device and insert the applicator brush into the container to saturate the tip.
6. Using the disposable applicator brush, apply a generous amount of OptiBond Universal adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.
7. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air. The surface should have a glossy uniform appearance. If not, repeat steps 6 and 7.
8. Light cure for 10 seconds.*
9. Apply composite according to manufacturer's instructions for use.

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