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PANASIL®

Vinyl polysiloxane impression material
ISO 4823

Panasil® binetics Putty Fast | **Panasil**® binetics Putty Soft
Panasil® tray Fast Heavy | **Panasil**® tray Soft Heavy
Panasil® monophase Medium

KETTENBACHDENTAL
Simply intelligent

INSTRUCTIONS FOR USE

English

Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation.

Indications / Techniques

Panasil® binetics Putty Fast and Panasil® binetics Putty Soft are to be used as preliminary materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step putty-wash impression technique using a foil (plastic putty spacer)
- One-step putty impression technique for forming functional peripheries

Panasil® tray Fast Heavy and Panasil® tray Soft Heavy are to be used as heavy bodied materials for:

- One-step impression technique (simultaneous technique) using single or dual viscosities
- Two-step impression technique using dual viscosities (Panasil® tray Fast Heavy)
- Functional impressions (one-step: Panasil® tray Soft Heavy)

Panasil® monophase Medium is to be used as a medium bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays
- Fabricating full or partial dentures
- Reline impressions
- Use in the simultaneous mixing technique as well as the putty-wash and Multi tray techniques
- Transferring root posts when fabricating posts and cores indirectly

Patient Target Group

Patients being treated in the context of a dental procedure.

Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

Panasil® binetics Putty impression materials are not suitable for detailed impressions when used alone.

Potential Allergenic Substances

Platinum compound, 4-methoxyphenol.

Cautions

Do not use after expiration date.

Normally, impression materials do not cause allergic reactions; however, an allergic reaction cannot be ruled out in sensitive individuals.

If in doubt, we recommend performing an allergy test before use.

Do not use the material in case of known allergies to any of the ingredients or in case of contact allergies.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

To ensure an optimal connection between the two materials, after the trial positioning, the preliminary impression should also be rinsed with water and dried before the correction material is applied. Both steps of the putty-wash impression technique (preliminary impression and wash impression) should also be carried out in quick succession.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth / jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23°C (74°F) before applying. Otherwise, working and setting times will be affected.

Impression tray: Preparation and Adhesives

In general, all conventional impression trays can be used if the necessary pressure built up can be generated. If it is not possible to achieve sufficient retention with the impression material, brush the impression tray with Panasil® adhesive / Kettenbach Universal Adhesive. The surface of the impression tray must be repeatedly and fully brushed so that the entire surface of the tray is visibly blue / red. Allow the adhesive to dry in accordance with the information in the direction for use and cover with sufficient impression material in the impression tray.

Instructions for the Safe Handling of the 380-mL cartridge system; 5:1

Prior to first use please unscrew the plug counter-clockwise and pull it out. Check the outlet opening. Then screw on the mixer.

Please use only the dispensing and mixing components compatible with the 380 mL cartridge system for dispensing and mixing impression material.

The 380 mL cartridge is heavy and could be damaged if dropped. Damaged cartridges must not be used under any circumstances.

When using the material in a Sympress II dispensing device, please activate the speed setting indicated on the product label. Failure to comply with this instruction may modify the curing characteristics of the impression materials.

If you are using a unit other than the Sympress mixing device, please follow the instructions for the automatic dispensing and mixing system you are using.

When changing to a different mixing device, always discard at least 3 cm of mixed material.

Using the 380-mL cartridge system; 5:1 in Dispensing Device Sympress / Sympress II

By pressing the plunger return knob, allow the plunger to return to its original position automatically (while the cover is closed. Next, open the cover.

Insert the 380 mL cartridge with the pre-assembled mixer in the Sympress / Sympress II dispensing and mixing unit. If you are working with a previously used cartridge, the filled mixer will already be located on the cartridge as a cap.

If you are using a Sympress II dispensing device, please activate the speed setting indicated on the product label. A light signal will confirm your selection.

Close the device cover.

Using the Cartridge for the First Time

Dispense the impression material (ensure that the mixer is wound into place fully) until the base material and catalyst are mixed homogeneously; discard at least the first 3 cm when dispensing.

At least 3 cm of the mixed material must be discarded each time the plunger is restarted or when using the cartridge in different mixing units.

Load the impression tray (or application syringe) with the required amount of material by pressing the appropriate feed key. Hold the impression at an angle when loading and press against it lightly. Leave the mixer in the material when dispensing to avoid air bubbles.

Release the feed key as soon as the required amount of material has been dispensed. The filled mixer should remain on the cartridge as a seal.

Remove the used mixer before the next application by turning it to the left, visually check the outlet opening for blockage and attach a new mixer.

Further Use of a Cartridge Already in Use

Wind a new mixer onto the same mixing unit and begin the new application (material no longer has to be discarded). Proceed with the application as described above.

Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection

Impressions can be disinfected using, for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials. See manufacturer instructions.

Model Preparation

Impressions with Panasil® may be poured immediately after disinfection of the impression up to several weeks with high quality dental stone (class III – V). Impressions can also be poured several times.

Composition

Vinyl polysiloxanes, Hydrogen polysiloxanes, Filler, Pigments, Additives, Platinum compound.

Disposal

Completely emptied cartridges / mixers and set material can be disposed in the household rubbish.

Trademarks

Panasil® is a registered trademark of
Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: July 15, 2025

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® binetics Putty Fast	Type 0, Putty	23	5:1 380 mL cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Mixer
Panasil® binetics Putty Soft	Type 0, Putty	23	5:1 380 mL cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Mixer
Panasil® tray Fast Heavy	Type 1, Heavy bodied	32	5:1 380 mL cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Mixer
Panasil® tray Soft Heavy	Type 1, Heavy bodied	32	5:1 380 mL cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Mixer
Panasil® monophase Medium	Type 2, Medium bodied	35	5:1 380 mL cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Mixer

Technical Data

Product name	Mixing ratio	Working time at 23°C/74° F ≤	Working time at 35°C/95°F ≤	Intraoral setting time at 35°C/95°F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® binetics Putty Fast	5:1	1:30 min.	Not applicable	2:30 min.	4:00 min.	A 63	– 0.20	99.5	2.5
Panasil® binetics Putty Soft	5:1	2:00 min.	Not applicable	3:00 min.	5:00 min.	A 56	– 0.20	99.5	3.5
Panasil® tray Fast Heavy	5:1	1:20 min.	Not applicable	2:00 min.	3:20 min.	A 62	– 0.20	99.7	2.5
Panasil® tray Soft Heavy	5:1	2:00 min.	Not applicable	2:00 min.	4:00 min.	A 55	– 0.20	99.7	3.0
Panasil® monophase Medium	5:1	2:00 min.	1:00 min.	2:00 min.	4:00 min.	A 60	– 0.20	99.7	3.5

* Total setting time (removal time from mouth) from start of mix

MD Medical Device CE 0123

Manufacturer: Kettenbach GmbH & Co. KG, Im Heerfeld 7, 35713 Eschenburg, Germany
Distributed in the U.S. by: Kettenbach LP, 62-64 Enter Lane, Islandia, NY 11749, USA
37424/0726 Made in Germany

PANASIL®

Vinyl polysiloxane impression material
ISO 4823

Panasil® binetics Putty Fast | **Panasil**® binetics Putty Soft
Panasil® tray Fast Heavy | **Panasil**® tray Soft Heavy
Panasil® monophase Medium

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