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PANASIL®

Vinyl polysiloxane impression material
ISO 4823

Panasil® tray Fast Heavy | **Panasil**® tray Soft Heavy | **Panasil**® tray Soft Heavy Fast
Panasil® monophase Medium | **Panasil**® initial contact Regular
Panasil® initial contact Light | **Panasil**® initial contact Light Fast
Panasil® initial contact X-Light | **Panasil**® initial contact X-Light Fast
Panasil® contact two in one Light | **Panasil**® contact plus X-Light

KETTENBACHDENTAL
Simply intelligent

INSTRUCTIONS FOR USE

English

Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation. Panasil® initial contact displays a high initial hydrophilicity.

Indications / Techniques

Panasil® tray Fast Heavy, Panasil® tray Soft Heavy and Panasil® tray Soft Heavy Fast are to be used as heavy bodied materials for:

- One-step impression technique (simultaneous technique) using single or dual viscosities
- Two-step impression technique using dual viscosities (Panasil® tray Fast Heavy)
- Functional impressions (one-step: Panasil® tray Soft Heavy)

Panasil® monophase Medium is to be used as a medium bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays
- Fabricating full or partial dentures
- Reline impressions
- Use in the simultaneous mixing technique as well as the putty-wash and Multi tray techniques
- Transferring root posts when fabricating posts and cores indirectly

Panasil® initial contact, Panasil® contact plus X-Light and Panasil® contact two in one Light are to be used as syringeable impression materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step impression technique using a spacer
- One-step impression technique (simultaneous technique) using dual viscosities
- Reline impressions
- Fabricating full or partial dentures

Patient Target Group

Patients being treated in the context of a dental procedure.

Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

Potential Allergenic Substances

Platinum compound, 4-methoxyphenol.

Cautions

Do not use after expiration date.

Normally, impression materials do not cause allergic reactions; however, an allergic reaction cannot be ruled out in sensitive individuals.

If in doubt, we recommend performing an allergy test before use.

Do not use the material in case of known allergies to any of the ingredients or in case of contact allergies.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

To ensure an optimal connection between the two materials, after the trial positioning, the preliminary impression should also be rinsed with water and dried before the correction material is applied. Both steps of the putty-wash impression technique (preliminary impression and wash impression) should also be carried out in quick succession.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth / jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23°C (74°F) before applying. Otherwise, working and setting times will be affected.

Impression tray: Preparation and Adhesives

In general, all conventional impression trays can be used if the necessary pressure built up can be generated. If it is not possible to achieve sufficient retention with the impression material, brush the impression tray with Panasil® adhesive / Kettenbach Universal Adhesive. The surface of the impression tray must be repeatedly and fully brushed so that the entire surface of the tray is visibly blue / red. Allow the adhesive to dry in accordance with the information in the direction for use and cover with sufficient impression material in the impression tray.

Instructions on the Use of Cartridge Material in the Applyfix 4 Dispensing Gun DS-50 1:1 / 2:1

Insert the cartridge into the Applyfix 4 dispensing gun DS-50 1:1/2:1. Ensure that the notches on the cartridge base are pointing down. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge and turn cap in a clockwise direction to lock in place.

Either dispense the necessary amount of impression material directly into the impression tray, on to the preliminary impression or into an impression syringe. Use the cartridge with the mixer attached and intraoral tip or the impression syringe to apply the material around the prepared teeth.

Replace the plug after use. Before using the cartridge again, remove the cartridge cap or the filled mixing tip and check the cartridge outlet openings to ensure that no polymerized materials have sealed the openings. Clear any clogged openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a new mixing tip and continue as described.

Cleaning and Disinfection of the Applyfix 4

Remove the cartridge.

Remove any uncured paste using a cloth soaked in alcohol.

The Applyfix 4 dispensing gun can be used several times if properly disinfected. Kettenbach recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Kettenbach recommends the use of glutaraldehyde-based solutions.

Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection

Impressions can be disinfected using, for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials. See manufacturer instructions.

Model Preparation

Impressions with Panasil® may be poured immediately after disinfection of the impression up to several weeks with high quality dental stone (class III – V). Impressions can also be poured several times.

Composition

Vinyl polysiloxanes, Hydrogen polysiloxanes, Filler, Pigments, Additives, Platinum compound.

Disposal

The empty cartridges / mixing tips / applicator tips and set material can be disposed in the household rubbish.

Trademarks

Panasil® is a registered trademark of
Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: July 15, 2025

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® tray Fast Heavy	Type 1, Heavy bodied	34	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® tray Soft Heavy	Type 1, Heavy bodied	34	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® tray Soft Heavy Fast	Type 1, Heavy bodied	34	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® monophase Medium	Type 2, Medium bodied	35	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® initial contact Regular	Type 2, Medium bodied	37	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® initial contact Light	Type 3, Light bodied	41	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® initial contact Light Fast	Type 3, Light bodied	41	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® initial contact X-Light	Type 3, Light bodied	44	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® initial contact X-Light Fast	Type 3, Light bodied	44	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® contact two in one Light	Type 3, Light bodied	38	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip
Panasil® contact plus X-Light	Type 3, Light bodied	42	1:1 50 ml cartridge	Applyfix 4 dispensing gun DS-50 1:1/2:1	Mixing tip

Technical Data

Product name	Mixing ratio	Working time at 23°C/74° F ≤	Working time at 35°C/95°F ≤	Intraoral setting time at 35°C/95°F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® tray Fast Heavy	1:1	1:00 min.	/	2:00 min.	3:00 min.	A 62	– 0.20	99.5	3.0
Panasil® tray Soft Heavy	1:1	2:00 min.	/	2:00 min.	4:00 min.	A 55	– 0.20	99.5	3.0
Panasil® tray Soft Heavy Fast	1:1	1:00 min.	/	2:00 min.	3:00 min.	A 55	– 0.20	99.5	3.0
Panasil® monophase Medium	1:1	2:00 min.	1:00 min.	2:00 min.	4:00 min.	A 60	– 0.20	99.7	3.5
Panasil® initial contact Regular	1:1	1:30 min.	1:00 min.	2:30 min.	4:00 min.	A 46	– 0.20	99.7	3.0
Panasil® initial contact Light	1:1	1:30 min.	1:00 min.	2:30 min.	4:00 min.	A 46	– 0.20	99.3	3.5
Panasil® initial contact Light Fast	1:1	1:00 min.	0:30 min.	2:00 min.	3:00 min.	A 46	– 0.20	99.3	3.5
Panasil® initial contact X-Light	1:1	1:30 min.	1:00 min.	2:30 min.	4:00 min.	A 46	– 0.20	99.3	3.5
Panasil® initial contact X-Light Fast	1:1	1:00 min.	0:30 min.	2:00 min.	3:00 min.	A 46	– 0.20	99.3	3.5
Panasil® contact two in one Light	1:1	2:00 min.	1:00 min.	2:00 min.		A 46	– 0.20	99.7	3.5
Panasil® contact plus X-Light	1:1	2:00 min.	1:00 min.	2:00 min.		A 46	– 0.20	99.7	5.0

* Total setting time (removal time from mouth) from start of mix

MD Medical Device CE 0123

Manufacturer: Kettenbach GmbH & Co. KG, Im Heerfeld 7, 35713 Eschenburg, Germany Distributed in the U.S. by: Kettenbach LP, 62-64 Enter Lane, Islandia, NY 11749, USA 37433/0726 Made in Germany

PANASIL®

Vinyl polysiloxane impression material
ISO 4823

Panasil® tray Fast Heavy | **Panasil**® tray Soft Heavy | **Panasil**® tray Soft Heavy Fast
Panasil® monophase Medium | **Panasil**® initial contact Regular
Panasil® initial contact Light | **Panasil**® initial contact Light Fast
Panasil® initial contact X-Light | **Panasil**® initial contact X-Light Fast
Panasil® contact two in one Light | **Panasil**® contact plus X-Light

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