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PANASIL®

Vinyl polysiloxane impression material
ISO 4823

Panasil® Putty | **Panasil®** Putty Fast | **Panasil®** Putty Soft

KETTENBACHDENTAL
Simply intelligent

INSTRUCTIONS FOR USE

English

Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation.

Indications / Techniques

Panasil® Putty (Putty, Putty Fast, Putty Soft) is to be used as a preliminary material for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step putty-wash impression technique using a foil (plastic putty spacer)
- One-step putty impression technique for forming functional peripheries

Patient Target Group

Patients being treated in the context of a dental procedure.

Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

Panasil® Putty impression materials are not suitable for detailed impressions when used alone.

Potential Allergenic Substances

Platinum compound, 4-methoxyphenol.

Cautions

Do not use after expiration date.

Normally, impression materials do not cause allergic reactions; however, an allergic reaction cannot be ruled out in sensitive individuals.

If in doubt, we recommend performing an allergy test before use.

Do not use the material in case of known allergies to any of the ingredients or in case of contact allergies.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

To ensure an optimal connection between the two materials, after the trial positioning, the preliminary impression should also be rinsed with water and dried before the correction material is applied. Both steps of the putty-wash impression technique (preliminary impression and wash impression) should also be carried out in quick succession.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth / jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23°C (74°F) before applying. Otherwise, working and setting times will be affected.

Impression tray: Preparation and Adhesives

In general, all conventional impression trays can be used if the necessary pressure built up can be generated. If it is not possible to achieve sufficient retention with the impression material, brush the impression tray with Panasil® adhesive / Kettenbach Universal Adhesive. The surface of the impression tray must be repeatedly and fully brushed so that the entire surface of the tray is visibly blue / red. Allow the adhesive to dry in accordance with the information in the direction for use and cover with sufficient impression material in the impression tray.

Instructions for Use: Jar Material for Hand-mix Putty Preparation

Only use materials with the same lot number.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

Use the scoops to dispense equal amounts of base paste and catalyst paste.

The different color scoops must only be used for the materials with corresponding colors.

Close the containers carefully after use, and ensure that the lids are not interchanged.

Contamination of base paste and catalyst paste in the container renders the material unusable.

Knead the base paste and catalyst paste for 45 seconds until the material is homogeneously colored.

Should gloves be worn, test them for compatibility prior to mixing with a sample of the material to be mixed. Certain types of gloves can prevent the polymerization (e.g., latex gloves) or impair the adhesive bond (e.g., nitrile gloves). Using vinyl gloves is recommended.

Prior to taking the wash impression, the preliminary impression must be cleaned with water, dried and trimmed using standard methods.

Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection

Impressions can be disinfected using, for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials. See manufacturer instructions.

Model Preparation

Impressions with Panasil® may be poured immediately after disinfection of the impression up to several weeks with high quality dental stone (class III – V). Impressions can also be poured several times.

Composition

Vinyl polysiloxanes, Hydrogen polysiloxanes, Filler, Pigments, Additives, Platinum compound.

Disposal

The set material can be disposed in the household rubbish.

Trademarks

Panasil® is a registered trademark of Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: July 15, 2025

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique
Panasil® Putty	Type 0, Putty	24	1:1, 900 ml in jars 1:1, 400 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)
Panasil® Putty Fast	Type 0, Putty	24	1:1, 900 ml in jars 1:1, 400 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)
Panasil® Putty Soft	Type 0, Putty	24	1:1, 900 ml in jars 1:1, 400 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)

Technical Data

Product name	Mixing ratio	Working time at 23°C / 74°F ≤	Intraoral setting time at 35°C / 95°F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® Putty	1:1	2:00 min.	2:00 min.	4:00 min.	A 66	– 0.20	99.0	2.7
Panasil® Putty Fast	1:1	1:30 min.	2:00 min.	3:30 min.	A 66	– 0.20	99.0	2.7
Panasil® Putty Soft	1:1	2:00 min.	2:00 min.	4:00 min.	A 60	– 0.20	99.0	2.7

* Total setting time (removal time from mouth) from start of mix

MD Medical Device CE 0123

Manufacturer: Kettenbach GmbH & Co. KG, Im Heerfeld 7, 35713 Eschenburg, Germany
Distributed in the U.S. by: Kettenbach LP, 62-64 Enter Lane, Islandia, NY 11749, USA
37480 / 0726 Made in Germany

PANASIL®

Vinyl polysiloxane impression material
ISO 4823

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