CHARISMA® Opal Flow

CHARISMA® Opal Flow

66040670/11

(GB) Instructions for use
(ES) Instrucciones de uso
(BR) Instruções de uso
(GR) Οδηγίες χρήσης
(CZ) Návod k použití
(HU) Használati utasítás
(LV) Lietošanas instrukcija
(LT) Naudojimo instrukcija
(PL) Instrukcja obsługi
(HR) Upute za uporabu
(RO) Instrucțiuni de utilizare
(TR) Kullanma talimatı
(RU) Инструкция по применению
(AE) تعليمات الاستخدام
(KR) 사용 설명서
(BG) Инструкции за употреба
(ID) Petunjuk penggunaan
(TH) คำแนะนำในการใช้

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Distributed in USA/Canada exclusively by:
Kulzer, LLC
300 Heraeus Way
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dentist.

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CE 0197

CHARISMA® Opal Flow — Instructions for use (GB)

Product description
Charisma Opal Flow is a flowable, light-curing, radiopaque nano-hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations.

Composition
Charisma Opal Flow is based on multifunctional methacrylate monomers (UDMA, EBADMA) and contains approximately 65% by weight or 41% by volume inorganic fillers such as Ba-Al-F silicate glass, YbF_3 and SiO_2. The filler particle size is between 0.02 µm and 5 µm.

Shades
A1, A2, A3, A4, B2, OM (opaque medium), Baseliner

Indications
• Enlarged fissure sealing • Cavity lining - as the first layer for Class I and II cavities • Class V fillings • Minimally invasive Class I and II fillings in areas not subjected to masticatory forces • Minimally invasive Class III fillings • Small repairs of direct and indirect restorations combined with a suitable bonding agent • Splinting of mobile teeth

Contraindications
Application to dentine close to the pulp (less than 1 mm). Do not use in cases of a known allergy to any of the constituents of Charisma Opal Flow

Side effects
This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Warnings / Precautions
May cause sensitisation with skin contact. May irritate eyes. Avoid contact with the skin and eyes. In case of contact with skin, wash skin with soap and water. In case of contact with eyes, rinse with plenty of water and contact a doctor.
The tips supplied for precise, hygienic application are intended for single use and must be discarded after use.

PLTs are intended for single use. To avoid cross-contamination between patients, PLTs must be discarded after a single use!

For deep cavities with only a minimal layer of dentine over the pulp or those which have already had a small area exposed, the pulp should be indirectly or directly capped locally using calcium hydroxide compound. The calcium hydroxide compound should be covered with a glass ionomer cement before the Charisma Opal Flow is inserted in the cavity.

Prepared enamel and dentine surfaces on which Charisma Opal Flow will be adhesively luted must not come in contact with preparations containing eugenol as these could inhibit optimal polymerisation of the composite.

The use of a dental dam is recommended.

Method of application

Shade selection should be done using the Charisma Opal shade guide before inserting the dental dam as the teeth change their colour noticeably once dry.

1. Remove any caries and protect any hard dental structures according to the general guidelines for adhesive fillings. Then clean and dry the cavity.
2. A quality adhesive system (e.g., GLUMA® Self Etch, GLUMA® 2Bond) should be applied to the enamel and dentine according to the manufacturer's instructions before inserting the filling material. Adhesive-treated cavities must not be contaminated with moisture, blood or saliva.
3. Apply Charisma Opal Flow in thin layers (max. 2 mm, baseliner max. 1 mm) to the cavity. When applying the Charisma Opal Flow, ensure that it flows slowly and evenly to prevent the formation of air bubbles. If necessary, remove any air bubbles using a probe before polymerisation.
4. Polymerise each layer according to the following table using a suitable curing light (e.g., Kulzer Translux®) for 20 sec (light output: > 550 mW/cm²).

 When using plasma arc curing lights, the polymerisation time can usually be reduced by 25%. When using LED (light emitting diodes) curing lights, adhere to the polymerisation times specified for curing with halogen units, provided direct contact between the light source and the composite can be ensured. If the light intensity falls below 550 mW/cm², it is recommended to increase the time by 50%. It must be noted that polymerisation lamps can vary immensely in intensity and emitted spectrum. All polymerisation times specified are therefore minimum times. In case of doubt contact the lamp manufacturer for more information and check the functionality of the curing light in vitro before using it on patients.

Recommended Charisma Opal Flow polymerisation times using Translux curing light or similar halogen curing lights:

Shade	Time
A1, A2, A3, A4, B2, OM, Baseliner	20 s

5. During polymerisation a dispersion layer forms on the surface. This layer should not be touched or removed if further composite layers are to be applied.
6. Charisma Opal Flow can be prepared and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers (e.g., Charisma Supra) and polishing brushes are suitable for preparation.

Instructions
During long restorative procedures the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerisation reaction.
The polymerisation times specified above are minimum times and should not be reduced so that an optimal polymerisation rate, and the associated physical material properties, can be reached. In case of repairs, the surface should be roughened with a rotating instrument. We recommend prior silanization using a Silane coupling agent (e.g. Silicoup) to improve adhesion. Afterward, apply and cure the flow composite.
After ejection of the material from the syringe, release the plunger to prevent continued flow. Replace the cap of the syringe after every use.
Visually inspect for damage before first use. Damaged products must not be used.

Further precautions
For use by dental personnel only and for its intended use.

Storage
Store below 25°C (77°F). Do not expose to direct sunlight.
Do not use the material after the expiration date.
Expiration date ⌛ and batch number LOT: See note on syringe/packaging.
Please quote batch number and expiry date in all correspondence about the product.

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