CHARISMA® TOPAZ

CHARISMA® TOPAZ

GB Instructions for use
ES Instrucciones de uso
GR Οδηγίες Χρήσης
CZ Návod k použití
HU Használati utasítás
LV Lietošanas instrukcija
LT Naudojimo instrukcija
PL Instrukcja obsługi
BR Instruções de uso
HR Upute za uporabu
RO Instrucțiuni de utilizare
TR Kullanma talimatı
RU Инструкция по применению
AE تعليمات الاستخدام
CN 使用说明
KR 사용 설명서
SK Návod na použitie
BG Инструкции за употреба
EE Kasutusjuhend
SI Navodila za uporabo

6606725 3/11

KULZER
MITSUI CHEMICALS GROUP

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

CE 0197

CHARISMA® TOPAZ — Instructions for use (GB)

Product Description
Charisma Topaz is a light curing, radio-opaque nano composite. It is indicated for adhesively bonded, tooth colored restorations. It is suitable for anterior and posterior applications. Due to the excellent material properties, a beautiful aesthetic result and a durable, high lustre polish is easily reached.

Composition
Charisma Topaz is based on TCD-DI-HEA and UDMA (approximately 59% filler by volume, with 58% anorganic filler by volume and a particle size of 5 nm – 5 µm). • Barium Aluminium Fluoride glass • Pre-polymerized filler • Highly discrete nanoparticles

Indications
• Direct restoration of Class I – V cavities (according to G.V. Black) • Direct composite veneers • Shape corrections of teeth (i.e. diastemas, congenital defects in teeth, etc.) • Splinting of teeth loosened by trauma or periodontal disease • Indirect restorations (inlays, veneers) • Restoration of primary teeth • Core build-up • Repairs of porcelain, composite (in combination with an adequate repair-system)

Contraindications
Application to dentine close to the pulp (less than 1 mm). Do not use in cases of a known allergy to any of the constituents of Charisma Topaz. The use of this product is contraindicated in case of known or suspected allergies against (meth)acrylate compounds.

Side effects
This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Precautions
• Avoid contact with skin. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. May cause sensitisation by skin contact. Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
• In the case of deep cavity preparations in close proximity to the pulp or a direct pulpal exposure, cover the dentine or the exposure locally with a calcium hydroxide liner covered with a glass ionomer cement before applying the adhesive.
• Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of Charisma Topaz.
• Visually inspect for damage before first use. Damaged products must not be used. Composites should only be used at room temperature.

Available Shades

Opaque		OB (Opaque Bleach)	OL (Opaque Light)		OM (Opaque Medium)		OD (Opaque Dark)	
Universal		BL (Bleach Light)	A 1	A 2		A 3	A 3.5	A 4
			B 1		B 2	B 3		
				C 2				
Incisal	translucent	CL (Clear)						
	opalescent	CO (Clear Opal)						

Application
Select an appropriate shade or shade combination from the Charisma Topaz Shade guide prior to placing the rubber dam. This is important because teeth will significantly change colour when desiccated. If the tooth surface being treated is not to be instrumented by a diamond stone, carbide bur, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoride free polishing paste before restoration is placed. A non-sticking plastic, Teflon coated or gold plated placement instrument (such as the Kulzer Plasmacoat Instrument) is recommended for use in clinically placing, sculpting and shaping the material.
1. For ideal results, use Ivory® rubber dam (or similar) to control contamination from moisture, blood and saliva.
2. Prepare the cavity according to generally accepted rules for the clinical composite application. Then clean and dry the cavity.
3. Prior to application of Charisma Topaz restorative material to the cavity preparation, pretreatment of the tooth structure with an adequate enamel/dentine adhesive is mandatory (e.g. GLUMA® Bondings by Kulzer). Make sure that no excess remains in the corners of the cavity. Once treated with the adhesive, avoid contamination with saliva or blood. The use of a flowable composite material (e.g. Charisma® Flowables) allows the creation of a smooth cavity floor and facilitates the adaptation of the filling material. When using a layer of a flowable material, apply it according to the specifications in the respective instructions for use and cure it separately.
4. Apply Charisma Topaz restorative material in thin layers (max. 2 mm) and adapt carefully to the tooth surface or cavity walls. In case of severe discoloration of the underlying tooth or surface being restored, it may be necessary to mask this by using a thin layer of one of the recommended Charisma Topaz opaque dentine shades first.
5. To avoid excess translucency and achieve an ideal aesthetic result in the case of large Class IV cavities and anterior tooth reconstruction, it is necessary to utilize a Charisma Topaz opaque dentine shade (max. 1 mm) in the deepest layers, which can then be covered by the universal Charisma Topaz shades if necessary.
6. Polymerise Charisma Topaz in thin layers (max. 2 mm) according to the table below using an acceptable QHT (Quartz-Halogen-Tungsten) or LED curing light such as the Kulzer Translux® light curing unit or any other blue-light unit (wavelength of 440 – 480 nm; light output of 1550 – 550 mW/cm²). During polymerisation, an inhibited-layer forms on the surface which should not be touched, contaminated or removed during the build-up phase of the restoration. It serves as chemical connection for any subsequent composite layers that are added.

Charisma Topaz Shades: recommended curing time and maximum layer thickness for Translux®-curing units or other blue-light-curing units.

A1, A2, A3, A3.5, B1, B2, B3, C2, CL, CO	20 s
A4, OB, OL, OM, OD, BL	40 s
Incisal- and Universal-Shades	2 mm
Opaque-Shades (OB, OL, OM, OD)	1 mm

The data given are correct to the best of our knowledge and are based on internal test results.
Note: Curing lights can have an immense variation in intensity and emitted spectral output. User should check with the light unit manufacturer and verify its efficacy non-clinically before using in a patient.
7. After final build-up and polymerisation is completed, Charisma Topaz may immediately be finished and polished. For gross finishing, use diamond burs. Fine finishing and shaping is best accomplished with fine-grained diamond stones and multifluted carbide burs. Proximal surfaces may be finished with abrasive strips. A high final polish can be achieved with flexible, silicone rubber polishing points (e.g. Charisma EasyShine), cups, discs, or similar. Composite polishing pastes (e.g. Signum® HP diamond only extraoral) can be applied with polishing brushes and small buffing wheels if desired.

Instructions
• The light curing times mentioned above are the minimum times recommended and should not be reduced in order to achieve optimal restoration quality and physical properties of the material. • When building up large restorations or those, which require lengthy time to complete, the overhead operating light should temporarily be turned off or aimed away from the working area in order to prevent premature curing of the composite. • To avoid waste, after the material has been applied from the syringe, relieve the internal pressure by unscrewing the plunger a quarter turn. • Close the cap of the syringe after every use.

Special advice
Only to be used by dentists and for its intended use.

Storage
Keep out of direct sunlight. Do not store above 25°C (77°F). Do not use the material after the expiration date. Please quote batch number and expiry date in all correspondence about the product.
Expiration date ⧖ and batch number LOT: See note on syringe/packaging.

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