CHARISMA® Opal

66037751/04 – D.O.G./03.11

CHARISMA® Opal

GB Instructions for use
ES Instrucciones de uso
GR Οδηγίες Χρήσης
HR Upute za uporabu
PL Instrukcja obsługi
RO Instrucțiuni de utilizare
CZ Návod k použití
HU Használati utasítás
TR Kullanma talimatı
SK Návod na použitie
BG Инструкции за употреба
RU Инструкция по применению
CN 使用说明
KR 사용 설명서
TH คำแนะนำในการใช้
AE تعليمات الاستخدام
ID Petunjuk penggunaan

Heraeus

CHARISMA® Opal

Instructions for use GB

Description:
CHARISMA Opal is a universal, light curing, radio-opaque, Microglass® composite
It is indicated for adhesive, tooth coloured restorations. It is suitable for anterior and posterior applications. The Microglass®-filler system is designed for aesthetic results and durable, high lustre polish restorations.

Composition:
CHARISMA Opal is based on a BIS-GMA matrix and contains 58% filler (by volume), which is:
• Barium aluminium glass (0.02–2 microns) • Highly dispersive silica (0.02–0.07 microns)

Range of application:
• Fillings of the classes I–V according to Black • Shape and shade corrections of teeth (i.e. diastema closure, congenital defects) • Splinting of teeth loosened by trauma or periodontal disease • Restorations of primary teeth • Temporary repairs of porcelain and composite (in combination with an adequate repair-system) • Direct composite veneers • Core build-up

Clinical treatment protocol:
Select an appropriate shade or shade combination with the CHARISMA Opal shade guide prior to placing the rubber dam. This is important because teeth will significantly change colour when desiccated. If the tooth surface being treated is not to be instrumented by a diamond stone, carbide bur, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoride free polishing paste before restoration is placed. A non-sticking plastic, Teflon coated or gold plated placement instrument (such as the Heraeus Kulzer Plasmacoat Instrument) is recommended for use in clinically placing, sculpting and shaping the material. CHARISMA Opal filling material should be at room temperature 23°C (73°F), when being processed.
1. The use of a rubber dam (Ivory® or similar) is recommended to control contamination from moisture, blood and saliva.
2. Prepare the cavity according to generally accepted rules for the clinical composite application. Then clean and dry the cavity.
3. Prior to application of CHARISMA Opal into the cavity, pre-treatment of the tooth structure with an adequate enamel/dentine adhesive is mandatory (e.g. GLUMA® adhesives by Heraeus Kulzer). Once treated with the adhesive, avoid contamination with saliva or blood.
4. Apply CHARISMA Opal to the cavity in thin layers (max. 2 mm) and adapt carefully to the tooth surface or cavity walls. In case of clear discoloration of the underlying tooth or surface being restored, it may be necessary to mask this by using a thin layer of one of the recommended CHARISMA Opal opaque shades first.
5. To avoid excess translucency and achieve an ideal aesthetic result in the case of large class III or IV cavities and anterior tooth reconstruction, it is necessary to utilise a CHARISMA Opal opaque shade in the deepest layers, which can then be covered by the enamel and clear CHARISMA Opal shades.
6. Polymerise CHARISMA Opal in thin layers (max. 2 mm) according to the table below using an acceptable curing light such as the Heraeus Kulzer Translux® curing unit or any unit with a similar spectrum and light-output (> 400 mW/cm², DIN-EN-ISO 10650-2). During polymerisation, an oxygen inhibited-layer is formed on the surface which should not be touched, contaminated or removed during the build-up phase of the restoration. It serves as chemical connection for any subsequent composite layer that is added.

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Distributed in USA/Canada exclusively by:
Heraeus Kulzer, LLC
300 Heraeus Way
South Bend, IN 46614-2517
1-800-431-1785

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CE 0197

CHARISMA Opal curing depths and times

	Curing time with QHT curing light	Maximum layer thickness
A1, A2, A3, A3,5, B1, B2, B3, C2, CO (clear opal)	20 s	2 mm
A4, OL (opaque light), OM (opaque medium), OD (opaque dark), BL (bleach)	40 s	2 mm

The specifications of Heraeus products are based on extensive measurements in production, research and applications technology. We reserve the right to make slight product changes in the interests of technical progress.
Recommended CHARISMA Opal curing times with Translux curing units or other halogen and LED light curing units.
For Plasma Arc light curing units, the curing time can usually be reduced by 25%.
Note: Curing lights vary greatly in intensity and spectral output. All curing times are minimum times. User should check with the light unit manufacturer and verify its efficacy non-clinically before using in a patient.
7. After final build-up and polymerisation is completed, CHARISMA Opal may immediately be finished and polished. For gross finishing, use diamond burs. Fine finishing and shaping is best accomplished with fine-grained diamond stones and multifluted carbide burs. Proximal surfaces may be finished with abrasive strips. A high final polish can be achieved with flexible, silicone rubber polishing points, cups, discs, or similar. Composite polishing pastes can be applied with polishing brushes and small buffing wheels if desired.

Further recommendations:
• Before using for the first time, please check for damage by visually inspecting. Do not use damaged products. • Close syringe carefully with cap immediately after use. • To avoid waste, after the material has been applied from the syringe, relieve the internal pressure by unscrewing the plunger a quarter turn. • When building up large restorations or those, which require lengthy time to complete, the overhead operating light should temporarily be turned off or aimed away from the working area in order to prevent premature curing of the composite. • The light curing times mentioned above are the minimum times recommended and should not be reduced in order to achieve optimal restoration quality and physical properties of the material. • In the case of deep cavity preparations in close proximity to the pulp or a direct pulpal exposure, cover the dentine or the exposure locally with a calcium hydroxide liner covered with a glass ionomer cement before restoring with CHARISMA Opal. • Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of CHARISMA Opal.
Storage:
• Do not store above 25°C (77°F). Keep out of direct sunlight. Do not use the material after the expiration date.
Precautions:
• Wear suitable gloves (e.g. latex or nitrile). • Avoid skin contact. • Avoid eye contact. In case of contact with eyes rinse with plenty of water with open eyelid and consult an eye specialist. • Do not use CHARISMA Opal in patients, which are allergic against the ingredients of CHARISMA Opal. There is no evidence of allergies but nevertheless allergies to the product cannot be ruled out. • Contains methylacrylate.
• Only to be used by qualified dental personal.

Disclaimer:
Our technical advice, suggestions and recommendations for the use of this product is given in good faith and is based on our own tests and the best available current knowledge. But it is not to be construed as a warranty for use. This also applies where other manufacturer products have been referred to or the proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to determine their final suitability for their intended purpose. The final application, use and processing of these products are beyond our control, and therefore, entirely the responsibility of the user. In the event that it can be shown that our product was delivered to you in a defective or damaged condition, our liability is limited to the replacement value of the goods delivered by us and used by you. We will, of course, provide consistent high quality products within the scope of our quality control standards and General Conditions of Sale and Delivery. Dated: 03.2011

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