oxasil®

GB Instructions for use
GR Οδηγίες Χρήσης
CZ Návod k použití
HU Használati utasítás
LV Lietošanas instrukcija
PL Instrukcja obsługi
RO Instrucțiuni de utilizare
TR Kullanma talimatı
RU Инструкция по применению
AE تعليمات الاستخدام
SK Návod na použitie
BG Инструкции за употреба
UA Інструкція по застосуванню
AL Udhëzime përdorimi
SI Navodila za uporabo
RS Uputstva za upotrebu

66078009/11

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Romania

CE

oxasil® Instructions for use (GB)

oxasil®	EN ISO 4823	Colours	Dosage[1]		Mixing-time	Working-time[2]	Time in mouth
			BASE	Activator paste	s	min : s	min : s
Soft Putty	Type 0	yellow	1 spoon (10.7 g/5.5 ml)	6–8 graduations (0.32 g/0.29 ml)	30	1:00	4:00
Putty	Type 0	yellow	1 spoon (14.4 g/5.5 ml)	4–6 graduations (0.21 g/0.19 ml)	30	1:00	4:00
Mucosa	Type 2	orange	12 graduations (8.0 g/7.6 ml)	12 graduations (0.65 g/0.59 ml)	30	1:30	4:00
Light Flow	Type 3	blue	12 graduations (8.0 g/7.6 ml)	12 graduations (0.65 g/0.59 ml)	30	1:30	4:00
Very Light Flow	Type 3	blue	12 graduations (8.0 g/7.6 ml)	12 graduations (0.65 g/0.59 ml)	30	1:30	4:00

*1) At room temperature of 23°C (73°F), 50% rel. humidity.
Higher temperatures shorten, lower temperatures prolong these times.
*2) 1 graduation = 1.3 cm (when using Kulzer mixing pad).
The stated dosage should be observed.

Special Advice
Medical Device. For use by health care professionals only.

Product description
Elastomeric, silicone-based precision impression material, condensation curing

Indication
oxasil impression materials: The materials are for use in all standard impression techniques according to the professional's experience. E.g. use for fully and partly toothed jaws as well as for edentulous ones to generate impressions for total and partial dentures, reparations, provisional prosthetics, bite splints or the situation of the opposing jaw.

Contraindication
Do not use in case of a known allergy against a component.
Do not use in pregnant or nursing patients or children under 6 years of age.
Do not use in patients with known allergies to tin or peppermint oil.

Interactions
For best results usage with other silicon materials and Activators is not recommended.
If there are substantial undercuts or highly divergent teeth in the impression area, perform appropriate blocking especially when using oxasil Putty and oxasil Soft Putty.

Storage conditions
Do not store above 25°C (77°F).
Do not use oxasil impression materials if the expiry date has elapsed.
The expiry date (EXP) and batch code (LOT) are to be found on the containers

Disposal
Small quantities can be polymerized with the matching system component(s) and the cured solid material can be disposed of with the regular garbage. Larger quantities must be disposed of following the regulations of the local authorities.

User Advice
Avoid surplus material on impression tray or in patient's mouth; this may provoke gag reflex, swallowing of material or aspiration.
Additionally, to reduce the aforementioned hazards, put the patient in an upright position during impression taking.
The oxasil assortment is a complete system. Components can be individually combined for all impression techniques.
Visually inspect for damage before first use. Damaged products must not be used. Close carefully after use.

Adhesion
We always recommend the use of our Universal Adhesive.

Pouring the model
Pouring can take place 30 minutes after the impression is removed from the mouth. Impressions remain dimensional stable up to 7 days.

Type 0 = putty consistency
Type 1 = heavy bodied consistency
Type 2 = medium bodied consistency
Type 3 = light bodied consistency

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Warning and precaution:
Flammable. Do not inhale. Inhalation can cause health problems. Eyes, respiratory organs and skin can be irritated. In an emergency, rinse off immediately with copious amounts of water and consult an eye specialist. Avoid skin contact. Wear protective gloves/protective clothing/eye protection/face protection (EN 374, Category III/Type B), e.g. Dermatril from KCL, Germany). Do not swallow and ingest. If health problems arise after swallowing impression material, seek medical attention immediately. Intestinal blockage may arise in rare cases. Hypersensitive persons may experience an allergic reaction to the product.
According to EU Medical Device Regulation users/ patients are obliged to report occurred serious incidents with a medical device to its manufacturer and to the competent authority of the country, the incident occurred.

Dated: 2019-02

Recommended model material
All commercially available modeling stones can be used.
Our suggestion: Class III Gips, MOLDANO™/MOLDAVIT®
Class VI Gips, MOLDASTONE®/MOLDASYN®

Disinfection
Impressions should generally be treated by immersion for 10 minutes® in 5,25 % Sodium Hypochlorite solution. Other silicon specific disinfectant solutions should be used as directed by the products manufacturer.
*Sidex® Lab-Aides, INDS, PhD. Etanol Zarenstift, Renata Cavo Barrera, and Caribe Abrefix Sentra Cruz, DDS.
*PhD Effect of disinfectant agents on dimensional stability of elastomeric impression materials. J Prosthet Dent, 61, 427-4,1989

Warning and precaution:
Flammable. Do not inhale. Inhalation can cause health problems. Eyes, respiratory organs and skin can be irritated. In an emergency, rinse off immediately with copious amounts of water and consult an eye specialist. Avoid skin contact. Wear protective gloves/protective clothing/eye protection/face protection (EN 374, Category III/Type B), e.g. Dermatril from KCL, Germany). Do not swallow and ingest. If health problems arise after swallowing impression material, seek medical attention immediately. Intestinal blockage may arise in rare cases. Hypersensitive persons may experience an allergic reaction to the product.
According to EU Medical Device Regulation users/ patients are obliged to report occurred serious incidents with a medical device to its manufacturer and to the competent authority of the country, the incident occurred.

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