CE

i-SHINE®

EN	Composite polishing paste
LT	Poliravimo pasta kompozitams
DE	Komposit-Polierpaste
BG	Композитна полираща паста
CS	Světlem tuhnoucí nano hybridní kompozit
DA	Komposit poleringspasta
ET	Komposiit poleerimispasta
FR	Pâte de polissage composite
EL	Σύνθετη στιλβωτική πάστα
HU	Kompozit polírozó paszta
IT	Pasta per lucidare il composito
LV	Kompozītmateriāla pulēšanas pasta
NO	Kompositтpoleringstannkrem
PL	Pasta polerska do wypełnień kompozytowych
PT	Pasta de polimento compósito
RO	Pastă de lustruit pentru compozit
SK	Kompozitná leštiaca pasta
SL	Kompozitna polirna pasta
ES	Pasta para pulir resina
SV	Komposit-poleringspasta
TR	Kompozit polisaj macunu
RU	Композитная полировочная паста
UA	Композитна полірувальна паста

INSTRUCTION FOR USE — EN

DESCRIPTION

i-SHINE is non-sterile thixotropic composite polishing paste, developed for professional use in dentistry.

COMPOSITION

Glycerol 50-70%, ultra-fine aluminium oxide particles 30-50%, excipients.
i-SHINE does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivates, of human origin; tissues or cells of animal origin, or their derivates, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

appearance	white paste

INTENDED PURPOSE AND CLINICAL BENEFITS

i-SHINE restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- **For polishing of composite restorations (to restore smooth surface and natural luster, soften margins, improve overall appearance, and ensure longevity of the restoration).**

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, i-SHINE may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

INTENDED USER

i-SHINE is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental materials. There is no need for specific training.

STERILITY

i-SHINE is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

i-SHINE is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of paste is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No consumable components and accessories are supplied with the device.

INSTRUCTION FOR USE

1. Remove syringe cap.
2. Delicate push on plunger and dispense a portion of polishing paste on a mixing pad.
3. Apply a small amount of paste on the surface needed to be polished and on the polishing brush or polishing cup.
4. Run the handpiece slowly and continue circular movements. Later, process the composite surface up to moderate speed and pressure. Follow manufacturer's instruction for use of the brush or cup.
5. After polishing, remove residual paste and rinse with water.

WARNINGS

After the desired amount of material extruded, immediately close the syringe cap. Do not use i-SHINE for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. i-SHINE does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHELF-LIFE

Shelf-life of i-SHINE is 2 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

i-SHINE is safe and performs as intended if it is used in accordance to manufacturer's instruction for use.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF **ISHP1**	4 g polishing paste in syringe

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Glycerol 50-70 %, extremt fina aluminiumoxidpartiklar 30-50 %, hjälpämnen. Product does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivates, of human origin; tissues or cells of animal origin, or their derivates, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

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