# SDI WAVE/MV/HV

**WAVE, WAVE MV, WAVE HV**

| | |
|---|---|
| **ENG** | FLUORIDE RELEASING FLOWABLE COMPOSITE |
| **BUL** | ФЛУРООТДЕЛЯЩ ТЕЧЕН КОМПОЗИТ |
| **CES** | TEKUTÝ KOMPOZITNÍ MATERIÁL UVOLŇUJÍCÍ FLUORIDY |
| **DAN** | FLUORID AFGIVENDE FLYDENDE KOMPOSIT |
| **DEU** | FLUORIDABGEBENDES FLIESSFÄHIGES SUBSTRAT |
| **ELL** | ΡΕΥΣΤΗ ΣΥΝΘΕΤΗ ΡΗΤΙΝΗ Η ΟΠΟΙΑ ΑΠΕΛΕΥΘΕΡΩΝΕΙ ΦΘΟΡΙΟ |
| **FIN** | FLUORIDIA VAPAUTTAVA KELLUVA-KOMPOSIITTI |
| **FRA** | COMPOSITE A FAIBLE VISCOSITÉ DÉGAGEANT DU FLUOR |
| **HUN** | FLUORID LEADÓ FOLYÉKONY KOMPOZIT |
| **ITA** | COMPOSITO FLUIDO CON RILASCIO DI IL FLUORO |
| **LIT** | FLUORIDĄ IŠSKIRIANTIS SKYSTAS KOMPOZITAS |
| **NOR** | FLUORIDFRIGIVENDE FLYTENDE KOMPOSITT |
| **POL** | PŁYNNY KOMPOZYT UWALNIAJĄCY FLUOR |
| **POR** | COMPÓSITO DE ALTA FLUIDEZ COM FLÚOR |
| **RON** | COMPOZIT FLUXABLE CU ELIBERARE DE FLUOR |
| **SLK** | TEKUTÝ KOMPOZIT UVOLŇUJÚCI FLUORID |
| **SPA** | COMPOSITE FLUIDO QUE LIBERA FLUOR |
| **SWE** | FLUORIDFRISÄTTANDE FLYTANDE KOMPOSIT |

Made in Australia by SDI Limited
3-15 Brunsdon Street, Bayswater, Victoria 3153, Australia

Australia 1 800 337 003
Australia 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 0800 0225 5734
UK 0800 0225 5734
USA & Canada 1 800 228 5166

EU REP
SDI Germany GmbH
Hansestrasse 85
51149 Cologne, Germany
T: +49 2203 9255 0

Last Revised: 06-2025

PAK441094 L

10°C 50°F 25°C 77°F Allergies & Precautions

0297 only

## ENG INSTRUCTIONS FOR USE

### INTENDED USE

Wave is designed to restore teeth and intended to be used by dental professionals.

### DESCRIPTION:

Wave/Wave MV/Wave HV is a versatile fluoride releasing, flowable, radiopaque, light cured composites. Wave/Wave MV/Wave HV is directly extruded into the cavity preparation for maximising adaptation to the preparation. The Wave shade range consists of 13 popular shades - enamel A1, A2, A3, A3.5, B1, B2, B3, C2, C3, D3, gingival, bleach; and OA2 dentin shade. Wave/Wave MV/Wave HV is a Class 2/Group 1 polymer-based restorative material according to ISO4049, actively releasing fluoride ions.

### COMPOSITION WAVE:

Dimethacrylates (43% wt) (63% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

Fillers (56%wt) (36% vol)
- Strontium Glass
- Silicon Dioxide (hydrophobic silica)
- Strontium Fluoride
- Initiators (< 1%)
- Stabilisers (< 1%)
- Pigments (< 1%)

### COMPOSITION WAVE MV:

Dimethacrylates (38% wt) (60% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

Fillers (61%wt) (39% vol)
- Strontium Glass
- Silicon Dioxide (hydrophobic silica)
- Strontium Fluoride
- Initiators (< 1%)
- Stabilisers (< 1%)
- Pigments (< 1%)

### COMPOSITION WAVE HV:

Dimethacrylates (36% wt) (57% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

Fillers (63%wt) (42% vol)
- Strontium Glass
- Silicon Dioxide (hydrophobic silica)
- Strontium Fluoride
- Initiators (< 1%)
- Stabilisers (< 1%)
- Pigments (< 1%)

### INDICATIONS FOR USE

- Class V restorations
- Tunnel / air abrasion preparations
- Pit and fissure sealants
- Minimal Class I restorations (excluding occlusal surfaces), II, III and IV restorations
- Radiopaque cavity liners
- Cervical abrasion/erosion lesions
- Abfraction lesions
- Implant seals
- Incisal abrasions
- Repair enamel defects
- Minor core build-ups
- Temporary crowns
- Cement porcelain/ceramic/ composite veneers
- Splint mobile teeth
- Attach fibre bridges
- Restorations in deciduous teeth
- Repair porcelain restorations
- Block out undercuts in inlay/ onlay/crown preparations
- Cover stains

### CONTRA-INDICATIONS:

- **Do not** use if dry field cannot be established.
- **Do not** use for pulp capping.
- **Do not** use in conjunction with any eugenol containing materials.
- **Do not** use in anyone who has an allergy to acrylics or who has an allergic reaction to the material of this product.

**NOTE:** This product may cause skin irritation to some people. In such cases, discontinue use and seek medical attention.

### SIDE EFFECTS:

- Areas around the pulp should be protected with a suitable pulp/dentin protector to prevent any potential irritation (selectively apply a calcium hydroxide-based preparation in areas near the pulp).
- The components of Wave / Wave MV / Wave HV could, on rare occasions, cause sensitisation. In these cases, Wave / Wave MV / Wave HV should not be used.

### INTERACTIONS:

Eugenol containing substances may inhibit the polymerisation process of composite materials. Therefore, using such materials in conjunction with Wave / Wave MV / Wave HV must be avoided. The use of chlorhexidine, cationic mouthwashes, and plaque-disclosing substances may result in discoloration.

### INSTRUCTIONS FOR USE:

**Recommendations:**

1. Clean and isolate tooth employing standard techniques and instruments. Rubber dam is the preferred mode of isolation.
2. A conservative cavity should be prepared, employing standard techniques and instruments, to form a slightly rounded internal form. Margins of the cavity preparation should end in sound and supported enamel with no bevels in stress bearing locations. If bevelling is desired in a low stress location it should extend no further than 1mm at an angle of no greater than 45°. Pre-wedging is also recommended to ensure that the restored tooth will have an adequate contact point.
3. Prophy all surfaces to be etched including surfaces adjacent to the cavity with a fine pressure fluoride containing paste or a slurry of pumice and water.
4. Rinse thoroughly with water.

### Shade Selection

1. Ensure the tooth surface is clean and dry before colour matching. Do not desiccate.
2. Use a VITA® shade guide to select your shade before placement of rubber dam.

**NOTE:** Do not use uncured paste for color matching as there is a slight color change upon curing.

**Recommendation:**
When the adjacent teeth are significantly different in shade, try to match the lighter shade. For visual confirmation of shade selected, place small amount of composite on a surrounding tooth surface and light cure before bonding.

### Isolation:

Isolation techniques must be used to prevent contamination. Rubber dam is the preferred mode of isolation. Do not use if a dry field cannot be established.

### Direction Restoration:

1. Thoroughly dry the surface to be etched with dry, oil-free air. Etch tooth surface using the 'selective etch' or 'total etch' technique with with 37% phosphoric acid according to the manufacturer's instructions.
2. Wash thoroughly with water.
3. Remove excess water. Keep moist. Avoid contamination e.g. saliva.
4. Apply bonding agent, such as Zipbond (SDI Limited), according to the manufacturer's instructions to saturate all internal surfaces.
5. Directly extrude the flowable composite into the prepared cavity in increments of 2mm or less

⚠ **CAUTION:** Ensure that the dispensing tip is firmly attached to the syringe to prevent the tip hub securely onto the syringe.

**NOTE:** For best results, apply restorative at room temperature (23 ± 1°C, 73 ± 2°F).

⚠ **CAUTION:** Air bubbles may be present when extruding the material. In these instances, use a probe to burst air bubbles before light curing.

6. Light cure the composite using a high power LED curing light (440-480nm wavelength with a minimum light intensity of 1,000 mW/ cm2), such as the Radii-Cal CX (SDI Limited), for a minimum of 20 seconds in increments of 2mm.
7. Finish the restoration using standard techniques. Instruct patient not to eat or drink anything for at least one hour after the procedure.

**NOTE:** Cooling with water spray and effective suction whenever possible when finishing the restoration is recommended.

② **CAUTION**
Wave / Wave MV / Wave HV syringes are for single use only. Reuse will result in cross contamination.

### Cementation of porcelain, ceramic and composite veneers:

Prior to cementation with Wave / Wave MV / Wave HV, the restoration must be pre-treated according to the manufacturer's instructions for use. Wave / Wave MV / Wave HV is only suitable for the cementation of restorations that have been etched and are of appropriate thickness to allow polymerisation light to pass through, to ensure that the light reaches Wave / Wave MV / Wave HV, so that it can cure completely.

### STORAGE AND HANDLING

Store at temperatures between 10° and 25°C (50° - 77°F). For best results, apply restorative at normal room temperature (23 ± 1°C, 73 ± 2°F). When not in use replace the cap tightly.

⚠ **PRECAUTIONS**
Avoid prolonged contact of etchant, adhesive and composite with the skin or oral tissue, as it may cause irritation and/or attack oral tissues or skin sensitisation. Patients with known material sensitivity to any of the materials to this product should not use it. For professional use only.
- Keep out of reach of children and pets.
- Do not take internally.
- Do not use after expiry date.
- Do not dispense material to direct light.
- Do not intravenously inject the material.
- Do not refrigerate.
- Polymerisation may be started by ambient lighting in the dental operating light.

### FIRST AID

**Eye (contact):** Wash thoroughly with water and consult medical attention.
**Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
**Ingestion:** Rinse mouth with water and discard. Drink plenty of water/milk. Do NOT induce vomiting. Seek medical attention.
**Inhalation:** If symptoms expected under normal conditions of use. If affected, move to well ventilated area.

If a serious incident has occurred due to this product, please report this incident to the manufacturer of this product and to the competent authority in your state or country. Dispose in accordance with local regulations.

SDS available at www.sdi.com.au or contact your regional representative.

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Диметакрилати (43% wt) (63% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

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Диметакрилати (38% wt) (60% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

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Диметакрилати (36% wt) (57% vol)
- Urethane Dimethacrylate
- Triethylene Glycol Dimethacrylate
- Ethoxylated Bisphenol A Dimethacrylate

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10. Aplikujte bond, ako je Zipbond (SDI Limited) podľa pokynov výrobcu, aby sa nasýtili všetky vnútorné povrchy.