Prior to use, carefully read the instructions for use. EN

G-CEM LinkForce™
ADHESIVE RESIN CEMENT
G-Premio BOND DCA
DUAL CURE ACTIVATOR FOR G-Premio BOND
G-CEM™ Try-In Paste
TRY-IN PASTE FOR G-CEM LinkForce™ AND G-CEM™ Veneer

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE
1. Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges.
2. Cementation of metal, ceramic, fiber posts, and cast post and cores.
3. Cementation of all ceramic and composite veneers.
4. Permanent cementation of crowns and bridges on implant abutments.

CONTRAINDICATIONS
1. Direct pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, acetone or ethanol.

PRODUCT DESCRIPTION
G-CEM LinkForce is a self adhesive resin cement system consisting of the cement paste and G-Premio BOND DCA and Try-In Paste, and used with a bonding agent (G-Premio BOND / G2-BOND Universal 1-PRIMER). The cement is a dual-cure, radio-opaque resin cement to be used intra-orally and classified as a Type 2 and Class 3 per ISO standard 4049.
G-CEM LinkForce has a refractive index equal to or greater than 2 mm aluminum (dentin = 1 mm, enamel = 2 mm).
The particle size of inorganic fillers range is 0.01-1.0 µm.
The total amount of inorganic filler is approximately 42% vol%.

COMPOSITION
G-CEM LinkForce
Base: Barium glass, dimethacrylate, initiator, pigment, silicon dioxide, stabilizer
Catalyst: Barium glass, dimethacrylate, silicon dioxide, initiator, stabilizer
G-Premio BOND DCA
Water, ethanol, initiator

DIRECTIONS FOR USE
A.CEMENTATION TECHNIQUE for all inlays, onlays, crowns, bridges and veneers
1. TRY-IT OF THE RESTORATION
 a) Remove the temporary restoration and clean the prepared tooth using a brush and slurry of pumice to remove any remaining temporary cement.
 b) Check the fit, adjustment and occlusion of the final restoration. As necessary, use G-CEM Try-In Paste for trial-fit.

HOW TO USE G-CEM Try-In Paste
1. Hold syringe upright and remove white cap by turning counterclockwise.
2. Promptly and securely attach dispensing tip to the syringe by turning clockwise.
3. Remove the temporary restoration and clean the prepared tooth using a brush and slurry of pumice to remove any remaining temporary cement. Rinse the prepared tooth with water spray.
4. Check fit and occlusion of the final restoration. Dispense G-CEM Try-In Paste into internal surface of the restoration and seat onto preparation.
5. Place wing cap back on promptly after use.
6. Check fit and aesthetics of the restoration. Remove restoration and thoroughly rinse G-CEM Try-In Paste with water. Make sure that G-CEM Try-In Paste is thoroughly removed before proceeding to bonding step as residual paste would cause decrease in bond strength.

2. PRE-TREATMENT OF THE RESTORATION
 a) Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.
 • For glass ceramics: etch with hydrofluoric acid. Rinse thoroughly with water and dry with an oil-free air syringe.
 • For metal, zirconia, alumina, hybrid ceramics (e.g. CERASMART) and composite: sandblast* and blow clean with an oil-free air syringe. Optionally, clean in an ultrasonic bath. Rinse thoroughly with water and dry with an oil-free air syringe. Follow manufacturer's instructions on the parameters of sandblasting and acid etching.
 *For CERASMART: If sandblaster is not available, etch with hydrofluoric acid for 60 seconds.
 b) Apply G-Multi PRIMER to the bonding surface and dry with an oil-free air syringe.
3. PRE-TREATMENT OF THE PREPARATION
 a) Clean the preparation with pumice and water. Dry by gently blowing with an oil-free air.
 b) For tooth preparation, select from the following three techniques.
 Self-etch technique
 Apply G-Premio BOND to bur-cut enamel and dentin, without a separate etching step.
 Note:
 Un-cut enamel should always be treated with a 35-40% phosphoric acid gel (e.g. GC ETCHANT) for 10-15 seconds, rinsed with water for 5 seconds and gently dried.
 Selective etching of enamel
 Prior to applying G-Premio BOND to enamel and dentin, etch enamel only with 35-40% phosphoric acid gel (e.g. GC ETCHANT) for 10-15 seconds, rinse for 5 seconds and gently dry.
 Note:
 In case of veneer cementation, select either selective etching or total etch.
 Clinical hint:
 G2-BOND Universal can be used in place of G-Premio BOND. If G2-BOND Universal is used, only the 1-PRIMER component is needed.
 c) For composite core and metal core, roughen the surface with a diamond bur and proceed directly to the application of G-Premio BOND.
 Clinical hint:
 G2-BOND Universal can be used in place of G-Premio BOND. If G2-BOND Universal is used, only the 1-PRIMER component is needed.
 d) The prepared surface (tooth, composite, metal abutment) should be dry. Apply G-Premio BOND onto the entire abutment surface. Leave undisturbed for 10 seconds. Using an air syringe dry the surface for 5 seconds with MAXIMUM AIR PRESSURE. Light cure using a visible light curing unit.
 Clinical hint:
 G2-BOND Universal can be used in place of G-Premio BOND. If G2-BOND Universal is used, only the 1-PRIMER component is needed.

Table 1. Irradiation time

	Distance from light guide tip <10 mm	Distance from light guide tip >10 mm
Halogen / LED (>700 mW/cm²)	10 seconds	20 seconds
High power LED (>1200 mW/cm²)	5 seconds	10 seconds

The effective wavelength range of each dental curing unit must cover 450-480 nm.

Notes:
1. Light cure completely for effective bond strength. Lower light intensity will cause insufficient adhesion.
2. In case you prefer not to light-cure G-Premio BOND before cementation, G-Premio BOND can be mixed with G-Premio BOND DCA (Dual Cure Activator) in a 1:1 ratio. Apply mixture, leave it undisturbed for 20 seconds then dry for 5 seconds with MAXIMUM AIR PRESSURE and expose surface to oxygen.
3. In case of veneer cementation, G-Premio BOND should always be used without DCA and light cured.

4. DISPENSING
 a) Before dispensing the first time from a new syringe or after a long interval in between use, bleed the pastes from the syringe.
 b) Remove cap of the syringe by rotating 1/4 turn counterclockwise. Attach a GC Automix Tip by aligning the V-shaped notch on the rim of the mixing tip with the V-shaped notch between the syringe barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise.
 c) Extrude material directly into restoration.
 Note:
 Used mixing tip can be left on the syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
5. CEMENTATION
 a) Coat internal surface of the restoration with sufficient cement and seat immediately. The working time is 3 minutes from start of mix at 23°C (73°F). Higher temperatures will shorten the working time.
 b) Maintain moderate pressure.
 Notes:
 1. Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.
 2. In case of veneer cementation it is recommended to cement a maximum of 2 veneers at a time due to relatively short working time.
6. EXCESS CEMENT REMOVAL
 a) Maintain moderate pressure and make sure the restoration remains in place.
 b) To remove excess cement, tack cure using a visible light curing unit for 1-2 seconds or allow the cement to self-cure until the excess cement feels rubbery)
 c) Remove excess cement with appropriate instruments, holding the restoration in position.
7. FINAL SET
 While maintaining moderate pressure, light cure each surface/margin for 20 seconds (Halogen/LED >700 mW/cm²). If the light application is difficult, allow to set for 4 minutes after seating the restoration. For veneer cementation, please make sure to use the light cure mode.
8. FINAL POLISHING AND ADJUSTMENTS
 If needed, finish and polish margins with appropriate instruments.
 Note:
 Be careful not to overheat.

B.CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and cores
1. PRE-TREATMENT OF THE POST
 Make sure that post is pretreated and handled according to manufacturer's instructions.
 a) Clean the bonding surface of posts with alcohol and dry with an air syringe.
 b) Apply G-Multi PRIMER to the bonding surface and dry with an air syringe. Refer to G-Multi PRIMER instructions.
2. TOOTH PREPARATION
 a) Prepare post space in usual manner.
 b) Clean post space and rinse with water. Sodium hypochlorite solution or EDTA is recommended to use for the chemical cleaning of the post space.
 c) Dry thoroughly using paper points. The post space should be dry.
 d) Mix G-Premio BOND and G-Premio BOND DCA in a 1:1 ratio and apply into the prepared post space. Leave for 20 seconds.
 Note:
 In case of post cementation, the use of G-Premio BOND DCA is essential for stable adhesion.
 Clinical hint:
 G2-BOND Universal can be used in place of G-Premio BOND. If G2-BOND Universal is used, only the 1-PRIMER component is needed.
 e) The prepared surface (tooth, composite, metal abutment) should be dry. Apply G-Premio BOND onto the entire abutment surface. Leave undisturbed for 10 seconds. Using an air syringe dry the surface for 5 seconds with MAXIMUM AIR PRESSURE. Light cure using a visible light curing unit.
 Clinical hint:
 G2-BOND Universal can be used in place of G-Premio BOND. If G2-BOND Universal is used, only the 1-PRIMER component is needed.
3. DISPENSING & CEMENTATION
 a) For dispensing, refer to A.4.a),b). In case of post cementation, GC Automix Tip for Endo can be used.
 b) Extrude material directly into the post space.
 c) Insert the post immediately into the post space within 1 minute after cement application.
 d) Maintain moderate pressure.
 Notes:
 1) The oral temperature accelerates the setting reaction of the cement.
 2) Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.
 3) When the cement is in contact with G-Premio BOND DCA, the setting reaction is accelerated.
 4) Keep the tip immersed in the cement while dispensing to avoid trapping air.
 5) DO NOT use lentulo spirals to place the cement in the post space as this may shorten the working time.
 6) The used mixing tip can be left on the syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
4. EXCESS CEMENT REMOVAL
 Keep moderate pressure making sure the post remains in place, and remove excess cement with appropriate instruments. In case of a cast post and core, tack cure each surface for 1-2 seconds with a visible light curing unit.
5. FINAL SET
 a) While maintaining moderate pressure, light cure all surfaces/margins for 20 seconds each (Halogen/LED >700 mW/cm²).
 b) Leave material undisturbed for 4 minutes after inserting post.

STORAGE
G-CEM LinkForce cement
Store in refrigerator (2-8°C / 36-46°F) when not in use.
G-Premio BOND DCA, G-CEM Try-In Paste
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

AVAILABLE SHADES FOR G-CEM LinkForce and G-CEM Try-In Paste
Translucent (clear translucent), A2 (A2 translucent), Bleach (bleach opaque), Opaque (universal opaque)
The shade guide of G-CEM LinkForce is available on the GC Website.

CAUTION
1. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
2. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid irritation, rather clean around GC COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. Do not use G-CEM LinkForce, G-Premio BOND and G2-BOND Universal 1-PRIMER in combination with eugenol containing materials as eugenol will hinder G-CEM LinkForce, G-Premio BOND and G2-BOND Universal 1-PRIMER from setting or bonding properly. Do not use G-CEM LinkForce, G-Premio BOND and G2-BOND Universal 1-PRIMER in combination with substances like hydrogen peroxide (H_2O_2) as H_2O_2 will hinder G-CEM LinkForce, G-Premio BOND and G2-BOND Universal 1-PRIMER from setting or bonding properly.
5. Use a protective light shield or similar protective eye wear during light curing.
6. When spilled on the table or floor, immediately wipe off with a dry cloth.
7. Avoid inhalation or digestion of material.
8. G-Premio BOND, G2-BOND Universal 1-PRIMER, G-Premio BOND DCA, G-Multi PRIMER are volatile. Use in a well ventilated place.
9. Dispose of all waste according to local regulations.
10. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
They can also be obtained from your supplier.

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients, device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants.
Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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