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Prior to use, carefully read the instructions for use. EN

G-CEM™ Try-In Paste

TRY-IN PASTE FOR G-CEM LinkForce™ AND G-CEM™ Veneer

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE
For selection of shade of cement for final cementation using G-CEM LinkForce and G-CEM Veneer.

CONTRAINDICATIONS
Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE
1. Hold the syringe upright and remove the wing cap by turning counterclockwise.
2. Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
 Note:
 Take care not to attach the dispensing tip too tightly. This may damage its thread. In order to ensure a tight connection make sure that the threads are free of remnant material.
3. Remove the temporary restoration and clean the prepared tooth using a brush and slurry of pumice to remove any remaining temporary cement. Rinse the prepared tooth with water spray.
4. Check the fit and occlusion of the final restoration. Dispense G-CEM Try-In Paste into the internal surface of the restoration and seat onto the preparation.
 Note:
 If material is not dispensed from the dispensing tip smoothly, remove the dispensing tip and extrude a small amount of material directly from the syringe.
5. Place the wing cap back on promptly after use.
6. Check the fit and aesthetics of the restoration. Remove the restoration and thoroughly rinse G-CEM Try-In Paste with water. Make sure that G-CEM Try-In Paste is thoroughly removed before proceeding to the bonding step as residual paste would cause the decrease in bonding strength.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

SHADE
A2, Translucent, Opaque, Bleach

PACKAGE
1 x 1.5 g (1.2 mL) syringe, 3 x DISPENSING TIP NEEDLE, 1 x Light protective cap

CAUTION
1. In case of contact with oral tissue or skin, remove and flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. To prevent cross contamination between patients, do not reuse the dispensing tip if used directly in the mouth. The dispensing tip cannot be sterilized in an autoclave or chemiclave.
4. Take care to avoid ingestion of the material.
5. Dispose of all wastes according to local regulations.
6. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
or for The Americas
https://www.gc.dental/america
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system:
vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire
MK16 8JS, United Kingdom

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