Prior to use, carefully read the instructions for use. EN

G-Multi PRIMER

PRIMER FOR GLASS CERAMICS, HYBRID CERAMICS, ZIRCONIA, ALUMINA, COMPOSITES AND METAL BONDING

G-Multi PRIMER is a priming agent designed to treat the adherent surface of glass ceramics, hybrid ceramics (e.g. CERASMART270), zirconia, alumina, composite and metal restorations. This is to promote a strong bond to resin-based materials as well as resin-based cements in case of luting procedures or direct composites in case of repairs.

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

1. Priming of adherent surface of indirect restorations made of glass ceramics, hybrid ceramics (e.g. CERASMART270), zirconia, alumina, composites and metal before luting procedures.
2. Intra-oral and extra-oral repair of fractured indirect restorations made of ceramics, hybrid ceramics, zirconia, alumina, and composites.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer or ethanol.

COMPOSITION

Ethanol, dimethacrylate, MDP, silane coupling agent, MDTP

DIRECTIONS FOR USE

A. For luting procedure
 1. Treat adherent surface of restoration according to the respective manufacturer's instructions and then clean it using a suitable cleaning method. We recommend cleaning in an ultrasonic bath or with alcohol, then rinse and dry with an oil-free air syringe.
 2. Dispense one drop of G-Multi PRIMER into a dispensing dish.
 3. Apply a thin layer to the adherent surface of the restoration using a micro-tip applicator and dry with an air syringe.
 4. Apply suitable resin-based luting cement to the primed surface of the restoration. For these steps, refer to the manufacturer's instructions for use of the corresponding luting cement.

B. For Intra-oral Repair
 1. Roughen adherent surface using a coarse diamond bur / a carborundum point or with sandblaster.
 2. Treat adherent surface of the restoration according to the respective manufacturer's instructions and then clean it using a suitable cleaning method. Rinse and dry with an oil-free air syringe.
 3. Dispense one drop of G-Multi PRIMER into a dispensing dish.
 4. Apply a thin layer to the fractured surface of the restoration using a micro-tip applicator and dry with an oil-free air syringe.
 5. Continue the repair by using a light-cured adhesive (e.g. G-Premio BOND) and a light-cured composite (e.g. G-ænial, Essentia), referring to their respective manufacturer's instructions for use.
 Note:
 If treated surface is contaminated with saliva or blood, rinse with water and dry. Re-treat the surface according to the manufacturer's instructions and re-apply G-Multi PRIMER according to the instructions above.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

CAUTION

1. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
2. G-Multi PRIMER is flammable and volatile. Do not use near open flame. Keep away from sources of ignition. Keep away from direct sunlight. Use in a well ventilated place.
3. In case of contact with oral tissues, immediately remove with a cotton pellet or a sponge. Rinse immediately after luting or repairing.
4. In case of contact with skin, remove immediately with a cotton pellet or a sponge. Rinse with soap and water.
5. In case of contact with eyes, flush immediately with water and seek medical attention.
6. Replace the cap of bottle immediately after use. If condensation (water) is present on the bottle of G-Multi PRIMER, immediately wipe it off.
7. Avoid inhalation or ingestion of material.
8. Do not mix with any other primers.
9. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
or for The Americas
https://www.gc.dental/america
They can also be obtained from your supplier.

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system:
vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire
MK16 8JS, United Kingdom

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Ethanol, Dimethacrylat, MDP, Silan-Haftvermittler, MDTP

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