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Prior to use, carefully read the instructions for use. EN

G2-BOND Universal

TWO COMPONENT LIGHT-CURED UNIVERSAL ADHESIVE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE
1. Direct bonding of light-cured composites and compomers to tooth structure
2. Bonding of dual-cured core build up composites to tooth structure as long as these materials are light-cured
3. Intraoral repair of composite, metal-based and zirconia / alumina-based restorations
4. Intraoral repair of ceramic restorations in combination with a silane coupling agent
5. Treatment of hypersensitive teeth
6. Sealing of tooth preparation (cavity or abutment) for indirect restorations
7. Treatment of exposed root surfaces
8. Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges and veneers
9. Post cementation

CONTRAINDICATIONS
1. Direct pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED IN COMBINATION WITH
1. Chemically cured composite resin.
2. Eugenol containing materials as they may hinder G2-BOND Universal from setting or bonding properly.
3. Desensitizers as they may hinder G2-BOND Universal from setting or bonding properly.

COMPOSITION
1-PRIMER: Acetone, water, dimethacrylate, 4-MET, initiator, MDP, silicon dioxide, MDTP, pigment, stabilizer
2-BOND: Dimethacrylate, initiator, silicon dioxide, stabilizer

DIRECTIONS FOR USE
I. BONDING PROCEDURE FOR DIRECT RESTORATIONS AND TREATMENT OF EXPOSED ROOT SURFACE
1. Cavity preparation
 Use rubber dam for protection. Make sure that any temporary material is thoroughly removed from the surface. Prepare tooth using standard techniques. The cavity should be sufficiently clean to ensure proper adhesion. Rinse and dry the prepared tooth surfaces by gently blowing with an air syringe.
 Note:
 For pulp capping, use calcium hydroxide.
2. Treatment of tooth surfaces
 Select one of the following 3 technique options:
 a) Self-etch technique
 Apply 1-PRIMER to bur-cut enamel and dentin, without a separate etching step.
 Note:
 Non bur-cut enamel should always be treated by selective etching of enamel technique.
 b) Selective etching of enamel
 Prior to applying 1-PRIMER to enamel and dentin, etch enamel only with 35-40% phosphoric acid gel for 10-15 seconds, rinse for 5 seconds and gently dry.
 c) Total-etch technique
 Prior to applying 1-PRIMER to enamel and dentin, etch enamel and dentin with 35-40% phosphoric acid gel for 10-15 seconds, rinse for 5 seconds and gently dry.
3. Application of 1-PRIMER
 a) Dispense
 Bottle type: Prior to dispensing, shake the bottle of 1-PRIMER thoroughly. Dispense a few drops into a clean dispensing dish. Replace bottle cap immediately after use.
 Note:
 Wipe any excess liquid from the nozzle before closing the cap.
 Unit Dose type: Before opening Unit Dose cap, ensure tip of the cap is facing up. Unscrew cap counterclockwise. Using the disposable applicator provided, mix the liquid thoroughly for 5 seconds, making sure to agitate material from the bottom and sides of the cap.
 Note:
 If the cap of the Unit Dose is too tight to open, use another Unit Dose. Do not try to open the Unit Dose cap with excessive force as this may cause the liquid inside to spill and cause accidental contact with skin and eyes.
 b) IMMEDIATELY apply to the prepared enamel and dentin surfaces using a disposable applicator.
 c) Leave undisturbed for 10 seconds after application.
 d) Dry thoroughly for 5 seconds with oil free air under MAXIMUM air pressure. Use vacuum suction to prevent splatter of the adhesive. The final result should be a thin adhesive film with a glossy surface which does not visibly move under further air pressure.
 Notes:
 1) When removing from refrigerator after long-term storage allow 1-PRIMER to stand at room temperature for several minutes prior to use.
 2) Apply 1-PRIMER immediately as the material contains a volatile solvent.
 3) 1-PRIMER is light-cured material. Protect the liquid from light after dispensing.
 4) Should the applied material be contaminated with water, blood or saliva prior to light curing, wash and dry the tooth and repeat the procedure by re-applying material.
4. Application of 2-BOND
 a) Dispense
 Bottle type: Dispense a few drops into a clean dispensing dish. Replace bottle cap immediately.
 Note:
 Take care not to apply 1-PRIMER and 2-BOND in the same well.
 Unit Dose type: Before opening Unit Dose cap, ensure tip of the cap is facing up. Unscrew cap counterclockwise.
 Note:
 If the cap of the Unit Dose is too tight to open, use another Unit Dose. Do not try to open the Unit Dose cap with excessive force as this may cause the liquid inside to spill and cause accidental contact with skin and eyes.
 b) Gently apply air to evenly disperse material into a uniform film thickness.
 Notes:
 1) 2-BOND is a light-cured material. Protect the liquid from light after dispensing.
 2) Remove excess material remaining on the tooth (other than the surfaces to be bonded) using a sponge or cotton pellet as residue material is difficult to remove after light curing.
 3) Should the applied material be contaminated with water, blood or saliva prior to light curing, wash and dry the tooth and repeat the procedure by re-applying material.
 4) Bring refrigerated 2-BOND to room temperature before use.
5. Light curing
 Light cure using a visible light curing unit.

Irradiation time		
	Distance from light guide tip <10 mm	Distance from light guide tip >10 mm
Halogen / LED (700-1200 mW/cm²)	10 seconds	20 seconds
High power LED (more than 1200 mW/cm²)	5 seconds	10 seconds

Notes:
1) Light cure completely for effective bonding strength. Lower light intensity may cause insufficient adhesion.
2) Use a protective light shield or similar protective eye wear during light curing.

6A. Placement of light-cured composites and compomers
 Follow manufacturer's instructions for use in order to place, contour and light-cure the composite resin material of choice.
6B. Placement of dual-cured composites
 To ensure sufficient bonding strength, it is recommended to light cure the dual-cure composite material (versus self-curing).
7. Finishing
 Adjust restoration, finish and polish using standard techniques.

II. REPAIR SITUATIONS
1. Preparation of restoration
 Roughen the adherent surfaces using a diamond point.
 Use rubber dam for protection. Rinse and dry the prepared surfaces by blowing with an air syringe. If the restoration surface contains glass ceramics or hybrid ceramics treat those surfaces using a silane coupling agent such as G-Multi PRIMER prior to application of G2-BOND Universal.
 Note:
 Do NOT apply G-Multi PRIMER to tooth surface.
2. Treatment of tooth surfaces
 Select one of the following three techniques described in I. BONDING PROCEDURE FOR DIRECT RESTORATIONS 2. Treatment of tooth surfaces.
 a) Self-etch technique
 b) Selective etching of enamel
 c) Total-etch technique
3. Application of G2-BOND Universal
 Refer to I. BONDING PROCEDURE FOR DIRECT RESTORATIONS 3-5.
4. Placement of light-cured composites / compomers
 Follow manufacturer's instructions for use in order to place, contour and light-cure the composite resin material of choice.
5. Finishing
 Adjust restoration, finish and polish using standard techniques.

III.TREATMENT OF HYPERSENSITIVE TEETH AND IMMEDIATE DENTIN SEALING
1. Application of 1-PRIMER
 Refer to I. BONDING PROCEDURE FOR DIRECT RESTORATIONS 2. Treatment of tooth surfaces, 3. Application of 1-PRIMER.
 Optionally, 2-BOND can also be applied after application of 1-PRIMER. Refer to I. BONDING PROCEDURE FOR DIRECT RESTORATIONS 4. Application of 2-BOND.
 Note:
 For the treatment of hypersensitive teeth, the self-etch technique is recommended.
2. Light curing
 Refer to I. BONDING PROCEDURE FOR DIRECT RESTORATIONS 5. Light curing.
3. Remove unpolymerized layer
 Remove the unpolymerized layer of G2-Bond Universal 1-PRIMER using a moist cotton pellet.

IV.CEMENTATION OF INDIRECT RESTORATIONS using G-CEM LinkForce with or without G-Premio BOND DCA
 Refer to the instruction for use of G-CEM LinkForce and follow its instruction.

V. POST CEMENTATION AND CORE BUILD-UP Using GRADIA CORE with G-Premio BOND DCA
 Refer to the instruction for use of GRADIA CORE and follow its instruction.

Bottle Covers
Use of G2- BOND Universal bottle covers is highly recommended. Bottle covers are designed to protect product from body heat, which might impact product working time. The bottle covers may be reused with subsequent refills. When bottle is empty, remove cover and clean with water or wipe with alcohol. Do not sterilize, autoclave, or clean with chemical disinfectants.
Note:
Bottle covers marked 1 and 2 are for 1-PRIMER and 2-BOND respectively. Be careful not to attach the bottle covers incorrectly.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F). If not in use for a prolonged period of time, store in refrigerator.

PACKAGE
G2-BOND Universal Bottle Kit: 1-PRIMER 5 mL liquid (1), 2-BOND 5 mL liquid (1), 1-PRIMER Bottle Cover (1), 2-BOND Bottle Cover (1), disposable dispensing dish (20), DISPOSABLE APPLICATOR II (50)
G2-BOND Universal 1-PRIMER Bottle Refill: 5 mL liquid (1)
G2-BOND Universal 2-BOND Bottle Refill: 5 mL liquid (1)
G2-BOND Universal Unit Dose Kit: 1-PRIMER 0.1 mL liquid (50), 2-BOND 0.1 mL liquid (50), DISPOSABLE APPLICATOR II (100)

CAUTION
1. In rare cases, the product may cause sensitivity in some people. If such reaction occurs, discontinue the use of the product and refer to a physician.
2. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
3. Avoid inhalation or ingestion of material.
4. 1-PRIMER is a highly flammable liquid and vapour. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
5. 1-PRIMER is volatile. Use in a well-ventilated place. Advise patients to breathe through nose.
6. G2-BOND Universal Unit Dose is for single use only. Do not attempt to recap or reuse. Unit Dose cannot be recapped properly once opened and the volatile solvent may evaporate.
7. In case of contact with eyes, flush immediately with water and seek medical attention.
8. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
9. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply COCOA BUTTER to areas which cannot be covered with a rubber dam.
10. In case of surface spills, immediately wipe off with a dry cloth.
11. Do not mix with non-designated products.
12. Dispose of contents / container in accordance with local regulations.
13. Do not mix 1-PRIMER with 2-BOND.
14. Do not mix 2-BOND with Dual Cure Activators.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
They can also be obtained from your supplier.

CLEANING AND DISINFECTION:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients, this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gge@gc.dental

Undesired effects - Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper country in your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised: 04/2025

CE 2797

GC

MANUFACTURED by:
GC DENTAL PRODUCTS CORP.
2-285 Terimatsu-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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