3M ESPE

RelyX™ Ultimate

(en) Adhesive Resin Cement

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3M Deutschland GmbH
Dental Products
Carl-Schurz-Str.1
41453 Neuss – Germany

3M, ESPE, Cavit, RelyX, Scotchbond and Sol•Vi are trademarks of 3M or 3M Deutschland GmbH.
© 2018 3M. All rights reserved.

CE 0123

78-8101-4683-0

ENGLISH

Product Description

RelyX™ Ultimate is a self-curing resin cement supplied in an automix syringe. It is used for the adhesive cementation of indirect restorations and available in various shades. RelyX™ Ultimate is used in combination with Single Bond Universal Adhesive.

The adhesive can be used either as a "Total Etch" or "Self-Etch" procedure. The additional application of a 37% phosphoric acid etching gel (Scotchbond™ Universal Etchant increases the adhesive strength of the adhesive even further. RelyX™ Ultimate contains silanated methacrylates. The proportion of inorganic fillers is about 43% by volume: the grain size (0.90%) is about 15 µm. The material is basic (pH value ca. 11 part base paste : 1 part catalyst).

ATTENTION: These Instructions for Use must be kept for reference for the duration of product use. Please see the pertinent information for use for details on all of the products mentioned below. The product may only be used when the product labeling is clearly readable.

Indications

Scotchbond Universal Adhesive:

- Etching of enamel and dentin using the acid-etching technique

RelyX™ Ultimate in combination with Single Bond Universal as adhesive for:

- Final cementation of all-ceramic, composite, or metal inlays, onlays, crowns and bridges; 2-3-unit Maryland bridges and 3-unit Inlay- Retain bridges
- Final cementation of fiber posts and screws
- Final cementation of all-ceramic or composite veneers
- Final cementation of all-ceramic, composite, or metal restorations on implant abutments

Indications such as Maryland and Inlay-only bridge (resin-bonded bridges) demand especially high adhesive bond strength. Regardless of the manufacturer of the cement and restoration, these indications may be exposed to a higher risk of debonding. Therefore, the use of a dentin adhesive (e.g., Scotchbond Ultimate) it is important to maximize bond strength. For the detailed bonding procedure, please refer to the instructions for use of the restoration material.

Contraindications

Maryland and Inlay-only bridges are contraindicated for patients with bruxism or parafunction.

Precautionary Measures

For Patients and Dental Personnel

- Single Bond Universal Adhesive: contact with eyes may cause severe eye damage. Wear eye protection to prevent injury. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- RelyX™ Ultimate base paste: contains sodium persulfate that may trigger an allergic reaction in certain individuals when in contact. This product may not be appropriate for use in those individuals with known sensitivity to sulfites since a cross-reaction may occur with sodium persulfate.
- Scotchbond Universal Adhesive: contact with eyes or skin may cause burns. Wear suitable protective clothing, gloves, and eye/face protection. In case of contact with eyes, rinse with plenty of water and seek medical advice.

For Patients

- RelyX™ Ultimate and Single Bond Universal Adhesive contain substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known sorbitic or peracrylic allergies. If undesired contact with unpolymerized material occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary, and discontinue future use.

For Dental Personnel

- RelyX™ Ultimate and Single Bond Universal Adhesive contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured paste. If skin contact occurs, wash immediately with soap and water.
- The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If RelyX Ultimate or Single Bond Universal Adhesive contacts the glove, remove and discard the gloves, wash hands immediately with soap and water and put on new gloves.
- If allergic reaction occurs, seek medical attention as needed.
- Single Bond Universal Adhesive: contains highly flammable solvents: electricity but5ss. Keep away from heat, sparks, open flame, pilot lights, and other ignition sources.

3M VOSS5s can be obtained from www.zmmm.com, or contact your local subsidiary.

Preparations during Processing

- Seal temporary restorations using a eugenol-free product (e.g., with RelyX™ Temp NE). The use of temporary cements containing eugenol can inhibit the polymerization process of RelyX™ Ultimate during the final cementation.
- Do not use hydrogen peroxide (H_2O_2) or sodium bicarbonate ($NaHCO_3$) as the residues can adversely affect the adhesive strength and setting reaction of RelyX Ultimate.
- For the best possible bond, do not treat the surface of Zirconia, aluminum oxide, or metal restorations with phosphoric acid.

Pulp Protection

To avoid cuff irritation, cover areas in close proximity to the pulp by applying small amounts of hard-setting calcium hydroxide material. To avoid tire impression during the cementing of the tire restoration, the pulp protection should be applied prior to taking an impression for the final restoration.

Tooth Preparation

Preparation of the Cavity/Tooth Stump

- Remove loose preparation debris by spraying with water.
- Clean the cavity or the cavity thoroughly; rinse with a water spray, and dry in with 2–3 bursts of water-free and -oil-free air, or use cotton pellets or similar materials to dry it off. Do not overdry.
 - The cavity should be just dry enough that the surface has a slightly glassy appearance. As is the case with dry permanent cement, overdrying can lead to post-operative sensitivity.

Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, mixing solutions containing EDTA, enzyme cleaners, etc., after the final cleaning. Their residues may have a detrimental effect on the bonding strength and setting reaction of the adhesive.

Preparation of Root Canals

- Treat the root canal endodontically as usual and seal canal filling with gutta-percha and removal of the root canal filling leaving 4 mm of gutta-percha apically.
- Clean the root canal with a 2.5%–5.25% solution of sodium hypochlorite ($NaOCl$).
- Rinse immediately with water, dry with paper tips, and continue with the application of the adhesive; see "Dosing and Applying Adhesive".
 - The use of a rubber dam during the cementation of posts is recommended.

Preparation for Veneers

- Prepare and clean the tooth as is usual for veneers.
- Subsequently, it must be etched, as described under "Etching of the Tooth Structure".

Preparation for Maryland and Inlay/Onlay Bridges

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as dimpler resin and/or approximal grooves.
- Subsequently, it must be etched, as described under "Etching of the Tooth Structure".
- Use a rubber dam and keep the cavity free of any contamination during cementation.

Preparation for Chairside CAD/CAM Restorations

Please follow the instructions for use of the restoration material. For Resin Nano Ceramic restorations (e.g., Lava Ultimate), please refer to the Lava Ultimate instructions for use.

Trial Fitting of the Restoration

- Remove the temporary restoration and thoroughly clean the tooth structure of any residues from the temporary cement.
- Make a trial insertion of the final restoration and check its fit, color, and contact points. Use RelyX™ Try-In Paste or similar material for this purpose.
 - If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely, e.g., with alcohol.
 - For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.

Pre-treatment of Restorations

Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and during the final cementation.

Pre-treatment of Metal Restorations

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use the disposable applicator to apply Single Bond Universal Adhesive to the entire surface to be luted. Let it react for 20 sec.
 - During treatment, rewet the disposable applicator as required.
- Subsequently, direct a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely.

For Maryland and Inlay-only bridges, please refer to "Pre-treatment of Maryland and Inlay-only Bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with water-free and -oil-free air.
- Use the disposable applicator to apply Single Bond Universal Adhesive to the entire surface to be luted and let react for 20 sec.
 - During treatment, rewet the disposable applicator as required.
- Subsequently direct a gentle stream of air over the liquid for about 5 sec until it no longer moves and the solvent has evaporated completely.

Alternatively, instead of Single Bond Universal Adhesive any silane can be applied in accordance with the instructions for use. After silanizing, RelyX™ Ceramic Primer allow to react for 5 sec, then air dry completely so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use the disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec or use a suitable silane in accordance with the instructions for use.

For Maryland and Inlay-only bridges, please refer to "Pre-treatment of Maryland and Inlay-only Bridges".

Pre-treatment of Composite and Resin Nano Ceramic (e.g., Lava Ultimate) Restorations

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance. Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Glass Fibers-reinforced Posts

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance. Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.
- Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Glass Fibers-reinforced Composite Posts

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance. Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.
- Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Glass Fibers-reinforced Composite Posts

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance. Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.
- Clean the blasted surface with alcohol and dry it with water-free and -oil-free air.

Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) 5% solution is recommended), rinsed with water, and dried.

- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Glass Fibers-reinforced Posts

Please follow the instructions for use of the restoration material. In the absence of cleaning instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and -oil-free air.
- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Metal Restorations or Inlay Bridges

In the absence of bleaching instructions, we recommend the following procedure for metal, Zirconia, and aluminum oxide ceramics:

- Coat (diffuse) the surfaces to be bleached with the micro-etching device (Cavi Prep) and the blast-coating agent Cojet Sand from a distance of 2–10 mm and blast for 10 sec, see instructions for use for Cojet Prep and Cojet Sand.
- Blow away the residues with the blasting agent with water-free and -oil-free air.
- Use a disposable applicator to apply Single Bond Universal Adhesive over the entire surface to be luted. Let it react for 20 sec before blowing a gentle stream of air over the liquid for about 5 sec. As an alternative, a suitable silane can be used in accordance with the instructions for use.

Pre-treatment of Ceramic Inlays, Onlays, Veneers, Crowns and Bridges

- Carefully check the surface of the treated teeth and the surrounding areas: remove any cement residues still remaining. This is especially important when cementing to implant abutments.

Hygiene and Disinfection

- We recommend the use of commonly available hygienic protective covering to avoid any contamination of the RelyX Ultimate automix syringes during the intraoral application.
- Clean any contaminated RelyX Ultimate automix syringes, Single Bond Universal Adhesive bottles, and Scotchbond Universal Etchant syringes with commonly used agents and disinfect them in accordance with the commonly used hygienic regulations.

Protection

- Scotchbond Universal Etchant syringes must be completely emptied before being disposed of for recycling. Small amounts of Etching Gel can be disposed of with normal waste water when diluted with water or when neutralized.
- Dispose of empty RelyX Ultimate automix syringes and Single Bond Universal Adhesive bottles with the normal household waste.

Notes

- RelyX Ultimate and Single Bond Universal Adhesive can prematurely polymerize when exposed to natural or artificial lighting. Avoid intensive light exposure during application.

Storage and Stability

Store RelyX Ultimate in the lid pouch at 15–25°C/59–77°F. After removal from the lid pouch, use RelyX Ultimate within 6 months and before the lapse of the expiration date. Constantly high humidity influences the setting and must be avoided.

Store Single Bond Universal Adhesive at Scotchbond Universal Etchant at 02–25°C/36–77°F.

Do not use after the expiration date.

Customer Information

No person is authorized to provide any information which deviates from the information given in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Customer is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

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3M, ESPE, Cubit, RelyX, Scotchbond and Sol-X are trademarks of 3M or 3M Deutschland GmbH.
© 2009. 3M. All rights reserved.
3M, ESPE, Cclet, RelyX, Scotchbond and Sol•Lve are trademarks of 3M or 3M Deutschland GmbH,