3M ESPE

RelyX™ U200

(en) **Self-adhesive Resin Cement**

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3M Deutschland GmbH
Dental Products
Carl-Schurz-Str. 1
41453 Neuss – Germany

3M, ESPE, ESPE, Lava, RelyX and Sol-Lex are trademarks of 3M or 3M Deutschland GmbH.
© 2014, 3M. All rights reserved.

44/2007/8075/02

(en) ENGLISH

Product Description

RelyX™ U200 is a dual-curing, self-adhesive resin cement for manual mixing. It is used for the adhesive cementation of indirect all-ceramic, composite, or metal restorations and for posts and screws. When RelyX U200 is used, bonding and conditioning of the tooth structure are not necessary.
The cement is available in various shades.
RelyX U200 contains bi-functional (meth)acrylate. The proportion of inorganic fillers is about 43% by volume; the grain size (D 90%) is about 12.5 µm. The mixing ratio, based on volume, is 1 part base paste : 1 part catalyst.

☞ These Instructions for Use must be kept for reference for the duration of product use. Please be see the pertinent information for use for details on all of the products mentioned below.

Indications

- Final cementing of all-ceramic, composite, or metal inlays, onlays, crowns, and bridges; 2–3-unit Maryland bridges and 3-unit inlay/onlay bridges (contraindicated for patients with bruxism or periodontitis)
- Final cementing of posts and screws
- Final cementation of all-ceramic, composite, or metal restorations on implant abutments
- Final cementation of Lava™ zirconia build-ups for two-piece abutments – only in accordance with the Instructions for Use for Lava™ Frame or Lava™ Plus respectively

Indications such as Maryland and inlay/onlay bridges (resin bonded bridges) demand especially high adhesive bond strength. Regardless of the manufacturer of the cement and restoration, these indications may be exposed to a higher risk of decementation. To achieve an optimal result when using RelyX U200, please refer to the sections "Tooth Preparation" and "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Precautionary Measures

For Patients and Dental Personnel

- Base paste: contains sodium persulfate that may trigger an allergic reaction in certain individuals when inhaled. This product may not be appropriate for use in those individuals with known sensitivity to sulfites since a cross-reaction may occur with sodium persulfate.

For Patients

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate and/or peroxide (especially peroxide) allergies.
- If prolonged contact with moist tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use.

For Dental Personnel

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water.
- The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If the product contacts the glove, remove and discard the glove, wash the hands immediately with soap and water and put on new gloves.
- If allergic reaction occurs, seek medical attention as needed.

3M MSDSs can be obtained from www.mmm.com, or contact your local subsidiary.

Precautions during Processing

- Seat temporary restorations using a eugenol-free product (e. g., with RelyX™ Temp NE). The use of temporary cements containing eugenol can inhibit the polymerization process of RelyX U200 during the final cementing.
- Do not use hydrogen peroxide (H_2O_2) as the residues can adversely affect the adhesive strength and setting reaction of RelyX U200.

Pulp Protection

To avoid pulpal irritation, cover areas in close proximity to the pulp by applying small amounts of hard-setting calcium hydroxide material. To avoid bite increases during the cementing of the later restoration, the pulp protection should be applied prior to taking an impression for the final restoration.

Tooth Preparation

Preparation of the Cavity/Tooth Stump

- Prior to final cementation, clean the prepared stump or the cavity thoroughly with pumice slurry, rinse with a water spray, and lightly air dry in only 2–3 bursts of water-free and oil-free air, or use cotton pellets to dry it off. Do not overdry!
 - The cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any enamel or dentin etching material, overdrying can lead to post-operative sensitivity.

Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions containing EDTA, etc., after the final cleaning with pumice slurry and water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.

Preparation of Root Canals

- Treat the root canal endodontically as usual (root canal filling with gutta-percha and removal of the root canal filling, leaving 4 mm of gutta-percha apically).
- Clean the root canal with a 2.5%–5.25% solution of sodium hypochlorite (NaOCl).
- Rinse immediately with water and dry with paper points.
- We recommend the use of a rubber dam during the cementation of posts.

Preparation for Maryland and Inlay/Onlay Bridges

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal condition should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as cingular rests and/or approximal grooves.
- Use a rubber dam and keep the cavity free of any contamination during cementation.
- Etch the enamel surface of the cavity with 37% phosphoric acid for 15–20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to selectively etch the enamel to avoid post-operative sensitivity.

Trial Fitting of the Restoration

- Remove the temporary restoration and thoroughly clean the tooth surface of any residue from the temporary cement.
- Make a trial insertion of the final restoration and check its fit and contact points.
 - If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.
 - For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.
- Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.

Pre-treatment of Restorations

Pre-treatment of Metal Restorations

Please follow the instructions for use of the metal/ceramic material. In the absence of devoting instructions, we recommend the following procedure:

- Blast the restoration surface to be luted with aluminum oxide ≤ 40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided devating instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use, e.g.: RelyX™ Ceramic Primer to react for 5 sec; then air dry so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

Alternative 1:

- Blast the restoration surface to be luted with aluminum oxide ≤ 40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Alternative 2:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet Prep and the blast-coating agent CoJet Sand from a distance of 2–10 mm and vertically to the surface for 15 sec; see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".

Pre-treatment of Composite Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Blast the composite restoration surface to be luted with aluminum oxide ≤ 40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and oil-free air.
- When using RelyX™ Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.

Pre-treatment of Implant Abutments

Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).

Pre-treatment of Maryland and Inlay/Onlay Bridges

Metal, Zirconia and Aluminium Oxide Surfaces:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet Prep and the blast-coating agent CoJet Sand from a distance of 2–10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

Dosing

When dosing the first time:

- Remove the clicker dispenser from the foil pouch and note the removal date on the clicker label. Discard the pouch.
- Remove the cap from the clicker dispenser.
- Dispense a small quantity of paste on a mixing pad by pressing down the clicker lever. This levels the filling quantity of base paste and catalyst paste in the clicker. Discard the paste which has been dispensed.

When dosing, always:

- Remove the cap.
- Dispense the paste on a mixing pad:
 - Press the clicker lever and hold it down.
 - When the paste stops flowing, wipe off the clicker openings on the mixing pad. Doing so dispenses the correct volume of the base paste and catalyst paste.
 - Dispense only as much paste on the mixing pad (number of clicks) as required for cementing.
- Cleaning the clicker openings:
 - Wipe the clicker openings with a clean cloth to prevent blockage from cross-contamination of the pastes.
- Replace the cap by pressing until it locks into place noticeably.

Times

The processing and setting times depend on the ambient and oral temperature. The times shown are valid for conditions relevant in practice. As is the case with every composite cement, the setting of RelyX U200 cement slows down significantly at room temperature.
RelyX U200 is a dual-curing cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!

	min:sec
Working time from start of mixing	02:30
Light-curing in a lightwave range of 400–500 nm – Single surface, from occlusal – Any other surface additional – RelyX Fiber Post posts, from occlusal	00:20 00:20 00:40
Self-curing – Setting time after start of mixing	06:00

Mixing

- Mix base paste and catalyst paste into a homogenous paste within 20 seconds using a spatula. Avoid incorporating air bubbles.
- Then cover the cement with a light protection cover to protect it from intensive light.

Application

- Protect the working area from contamination with water, blood, saliva and sulcal fluid during the application and setting phase.
- Apply the mixed cement evenly to the entire cavity and as appropriate to the bottom side of the inlay/onlay, or fill the crown with material.
- Seat the restoration firmly and stabilize long enough for the cement to fully set.

Application in the Root Canal

- Do not use Lentulo-Spirals to insert the cement in the root canal as this can exceed the working setting.
- Spread cement to the post.
 - As an alternative, the cement can be inserted into the root canal with the aid of an application tip, e. g., Centrix™ Needle Tubes (Centrix is a registered trademark of the U.S.-based company Centrix); this will minimize the entrapment of air bubbles.
- Immediately after the placement of the pretreated root canal; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX U200 is subject to oxygen inhibition, i. e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of the uncured layer during shaping/polishing without leaving behind a deficit.

Using light-curing:

- After brief curing (about 2 sec using a standard polymerization device), use a suitable instrument (e.g., scaler) to remove the excess. Larger volume excess can be removed more easily.
- Immediately light-cure after removal (see Polymerization and Shaping).

Using self-curing:

- Use a suitable instrument (e.g., scaler) to remove the excess during the self-curing phase (about 3 min after beginning the mixing in the "gel phase").
- Immediately cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.

Alternatively:

- Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.
- Cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.
- Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

Polymerization and Shaping

- We recommend light-curing the cement through the restoration when doing so entails no aesthetic restrictions. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts. The exposure time depends on the translucency of the post being used; for RelyX Fiber Post, it is 40 sec.
- Shape any remaining uneven areas on the edges of restorations and polish the marginal area with diamond devices, aluminum oxide coated discs (e.g., Sof-Lex™), and diamond polishing paste.
- Then check the occlusion.

After Shaping

- Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing on implant abutments.

Hygiene and Disinfection

- Clean a contaminated RelyX U200 with cleaning agents commonly used in the dental practice.
- Use a cloth soaked in a disinfectant commonly used in the dental practice to disinfect the clicker.

Notes

- RelyX U200 can prematurely polymerize when exposed to natural or artificial lighting. Therefore dispense the paste on a mixing pad only just before mixing and apply immediately after mixing; avoid intensive light exposure during application.
- Deviations from the specified mixing ratio (1:1) adversely affect the material properties.

Storage and Stability

Store RelyX U200 in the pouch at 15°–25 °C / 59°–77 °F.
After removal from the pouch, use the product within 6 months and before the lapse of the expiration date. Constantly high humidity accelerates the setting and must be avoided.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Date of the information: April 2013

RUS: Authorised representative of the manufacturer in RF:
ZAO "3M Russia", Russia 121614, Moscow, Krylatskaya str. bldg. 17/3
T. +7/495/7847474

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3M, ESPE, CoJet, Lava, RelyX and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
3M, ESPE, ClipJet, Lava, RelyX and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.