solventum

Photac™ Fil Quick Aplicap™

(en) **Light-Curing Glass Ionomer Dental Restorative – Radiopaque**

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Solventum Germany GmbH
Edisonstrasse 6
59174 Kamen, Germany
Solventum.com

MD CE 0123

Solventum Customer Care/Safety Data Sheet/Safety Information Sheet Information: U.S.A.1-800-634-2249

Vita is a trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

© Solventum 2026. Solventum, the S logo, Aplicap, Ketac and Photac are trademarks of Solventum or its affiliates.

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0240431.02
34-2002-4655-0

(en) ENGLISH

Intended Use
Glass ionomer based dental restorative.

Product Description
Solventum™ Photac™ Fil Quick Aplicap™ is a light curing resin modified glass ionomer cement in a capsule according to ISO 9917-2, Class 2. Since it bonds chemically to enamel and dentin, it permits gentle preparation. The product is available in A1, A2, A3, A3.5, B2, B3, and C4 shade based on the VITA Classical shade guide; it can be applied without cavity lining. It releases fluoride ions and is radiopaque. The applicable quantity of a capsule amounts to at least 0.13 ml.

Aplicap capsules are single-use devices.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use. The Instructions for Use are available at: eIFU.Solventum.com.

Indications for Use
- Fillings in Black's class V
- Wedge-shaped defects
- Fillings in deciduous teeth
- Temporary restorations in Black's class I and II

Intended Patient Population
All patients requiring dental treatment, unless the patient's condition limits the use.

Contraindications
The product is contraindicated for use with patients who have a known allergy to acrylate-based materials or contained allergens as listed in section "Precautions and Warnings".

Intended Users
Educated dental professionals who have theoretical and practical knowledge on usage of dental products.

Clinical Benefit
Restoration of initial and cavitated tooth structure lesions.

Precautions and Warnings

For Patients
Photac Fil Quick Aplicap liquid contains acrylates that may cause allergic reactions by skin contact in certain individuals. If prolonged contact with oral tissue occurs, flush with large amounts of water. If an allergic reaction occurs, remove the product, and discontinue future use; seek medical attention as needed.

For Dental Personnel
Photac Fil Quick Aplicap liquid contains acrylates that may cause allergic reactions by skin contact in certain individuals. To reduce the risk of allergic response, avoid exposure to these materials; in particular, avoid exposure to uncured materials. If skin contact occurs, wash skin with soap and water. The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If the liquid contacts the glove, remove and discard the gloves. Wash the hands immediately with soap and water and re-glove.
If an allergic reaction occurs, seek medical attention as needed.

The safety and performance of the product has not been clinically evaluated in pregnant or lactating women.

Potential Side-Effects
Potential undesirable side effects are local irritations, or local or systemic allergic reactions.

Residual Risks
Restoration debonding or fracture may lead in rare cases to swallowing or aspiration of the restoration or parts of it. In situations with deep carious lesions, inadequate restoration contouring, poor patient oral hygiene, or if the IFU were not followed, the following risks may occur:
Tooth sensitivity, pulpal inflammation, enamel decalcification, secondary caries, discoloration, gingivitis, tissue abrasion, changes in occlusion that could lead to temporomandibular disorder.

Instructions for Use

Preparation
- Remove carious dentin only. Undercuts are not required.
- Use water spray to remove debris.
- Blow the cavity dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry! **The cavity should be just so dry that the surface has a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid any contamination by saliva.

Pulp Protection
- To prevent pulpal irritation, the areas around the pulp should be covered by dabbing them with a hardsetting calcium hydroxide preparation.

Conditioning
- No further preparatory procedures are required. In case of an unfavorable cavity geometry, the use of Solventum™ Ketac™ Conditioner can further improve the adhesive stability obtainable.

Activation/Mixing
- Place the Solventum™ Aplicap™ Activator on a sturdy surface and insert the Aplicap capsule in the activator.
- Using the ball of your hands, depress the activator lever firmly and completely to the stop and hold it down for 2 to 4 sec.
 - Firmly and completely depressing the lever to the stop and holding it in this position is the only way to ensure that all of the liquid is pressed into the powder. If inadequate pressure is applied or if the lever is not depressed completely, it is possible that too little liquid will be mixed with the powder, increasing the viscosity and changing the product characteristics.
- Capsules must be mixed immediately after activation. Pauses between mixing and application must be avoided; otherwise, as soon as the material starts to set, removal from the capsule will be rendered difficult or even impossible.
- Mix with a high frequency mixing device (approx. 4300 oscillations per min), such as CapMix/Capsule Mixing Device or with the RotoMix/Rotating Capsule Mixing Device.
- Do not reuse. Reuse may result in compromising product integrity or lead to device failure.
- Apply with the Solventum™ Aplicap™ Applier.
 - The Aplicap Activator and the Aplicap Applier are color-coded in red.
 - **Caution: Avoid any delays between activation, mixing, and application as the initiated setting will render the application difficult or prevent it entirely.**

Application
- After mixing, insert the capsule into the applier and secure it in place by pressing the lever once. Before applying the paste intraorally, open the nozzle to the stop and press the applier several times until the mixed paste becomes visible in the intra-oral tip.
- Use Aplicap Applier **to apply layers of up to 2 mm maximum**. Avoid all types of contamination during application.
- Use halogen or LED light for curing.
 - Place the light source as closely as possible to the filling.
 - The maximum layer thickness (without intermediate curing) may not exceed 2 mm.
 - In case of cavities of more than 2 mm in depth, several layers have to be applied.
 - The self-curing mechanism ensures that the material also cures in areas sheltered from the light. In order to attain its optimum physical properties and to ensure immediate preparation, the material has to be light-cured carefully.
 - It is possible to add layers on material already cured if the surface is clean and dry.
- The material adheres to metal instruments and thus has to be washed off in cold water before it sets.

Times
The following times apply at a room temperature of 23° C/73° F:

	min:sec
Activation	0:02
Mixing in RotoMix/Rotating Capsule Mixing Device	0:10
Mixing in a high-frequency mixer, such as CapMix/Capsule Mixing Device	0:15
Manipulation from the beginning of the mixture	2:00
Light curing	0:20

At room temperatures of over 23° C/73° F the curing is accelerated. Shorter mixing times prolong the manipulation and curing time, whereas longer mixing times shorten the manipulation and curing time. The mixing time in high-frequency mixers can be selected at between 8 and 15 sec, and the mixing time in RotoMix/Rotating Capsule Mixing Device can be set to between 7 and 12 sec. If the manipulation time is exceeded, bonding to the enamel and dentin is lost.

Finishing
Finishing can be started directly after light curing.
- Use a sharp instrument to remove excess material.
- Use rotary instruments (such as polishing discs or silicone polishers) to finish under cold water.

Storage Conditions/Shelf Life
Store the product at 15–25° C/59–77° F.
Do not use after the expiration date.

Handling Before and/or After Use
- Please keep the single-use capsule in the original blister to guarantee the traceability of the batch (LOT).
- Humidity in the cavity prevents bonding significantly.
- Photac Fil Quick Aplicap polymerises in daylight or artificial lighting. Therefore, the material should not be dispensed from the capsule until immediately prior to manipulation. **The manipulation time is reduced drastically in case of application under the operating light!**

Inspection, Maintenance, Testing
Examine the Aplicap Activator and Aplicap Applier for damage, discoloration, and contamination before every use. Do not use damaged discolored, or contaminated Aplicap Activator and Aplicap Applier under any circumstances.

Cleaning the Aplicap and Maxicap Activator and Applier
Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Aplicap/Maxicap Activator and Applier.

Ingredients
The ingredients of each component of the product are indicated in decreasing order of their concentration.
Powder: Silane treated glass powder, aromatic amine.
Liquid: Polymeric acid, Methacrylate (HEMA), water, phosphorylated methacrylate magnesium salt, urethane dimethacrylate (UDMA), stabilizer.

SDS and Disposal
See the Solventum Safety Data Sheet (SDS)/Safety Information Sheet (SIS), available at www.Solventum.com or contact your local subsidiary, for disposal and ingredient information.

Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste, waste with sharps to avoid health risks because of improper handling.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Please report an undesirable outcome (e.g. adverse event, serious incident (EUI) occurring in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.

A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. The EUDAMED public website is ec.europa.eu/tools/eudamed. The product Basic UDI-DI is 01979982761020000000026R5.

The summary of safety and clinical performance (SSCP) is also available at: eIFU.Solventum.com.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code	LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number	REF	Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one single use only.
ISO 15223-1 5.4.3 Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use.
CE Mark	CE 0123	Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
ISO 15223-1 5.7.7 Medical Device	MD	Indicates the item is a medical device.
ISO 15223-1 5.7.10 Unique device identifier	UDI	Indicates a carrier that contains unique device identifier information.
Rx Only	Rx Only	Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Non-corrugated fiberboard	21 PAP	Indicates product packaging is made of non-corrugated fiberboard. Official Journal of the EC; Commission Decision (97/129/EC)
Low-density polyethylene	04 PE-LD	Indicates product packaging is made of low-density polyethylene. Official Journal of the EC; Commission Decision (97/129/EC)

For more information, see eIFU.Solventum.com

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