Prior to use, carefully read the instructions for use. EN

DIAPOLISHER PASTE

DIAMOND POLISHING PASTE FOR COMPOSITE AND PORCELAIN

DESCRIPTION
DIAPOLISHER PASTE is a diamond polishing paste for direct and indirect use.

INTENDED PURPOSE
Final high luster polishing of composite resin and porcelain restorations.

INDICATIONS FOR USE
Natural teeth or restoration lacking smoothness and aesthetics.
Areas of application:
- Extraoral and intra-oral polishing of direct and indirect restorations.

CLINICAL BENEFITS
Restoration of aesthetics and gloss

CONTRAINDICATIONS
Patients with known allergies to benzoate preservatives.

INTENDED USERS/SPECIAL TRAINING
- Dentists
- No special training required.

PATIENT TARGET GROUP
- Patient with direct or indirect restoration.

PERFORMANCE CHARACTERISTICS
The product's performance characteristics satisfy the requirements of the intended use.

DIRECTION FOR USE
I. Preparation
1) Finish and polish restoration using standard techniques.
2) Hold the syringe upright and remove the wing cap by turning counterclockwise.
3) Attach the dispensing tip correctly to the syringe by turning clockwise.
Note:
Take care not to screw the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnant material.
4) After attaching the dispensing tip, protect it with the protective cap until ready for use in order to avoid contamination or drying-out.

II. Dispensing
Slowly push the plunger and dispense sufficient material into a clean dispensing dish or on a mixing pad. The syringe can be stored either with the original cap or with the protective cap on the dispensing tip.

III. Polishing
Apply a small amount of paste onto an appropriate composite / porcelain polishing brush or cup. Polish all surfaces using moderate speed and pressure (up to 10,000 rpm). Refer to the respective manufacturer's instructions for use of the brush or cup. After polishing, rinse with water to remove the paste.
Note:
In case of contact with adjacent natural tooth structure, use moderate speed and pressure up to 5,000 rpm.

DISPOSAL
Dispose of all waste according to local regulations.

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

SIDE EFFECTS
There are no known side effects to date.

RESIDUAL RISKS
There are no known residual risks to date, other than those specified as part of contraindications, limitations, precautions, warnings and side effects.

PRECAUTION, WARNINGS AND LIMITATION OF USE
Do not use if the prescribed application technique cannot be applied.
1. To prevent splatter of paste and / or excessive heat generation always use rotary instruments at speeds below 10,000 rpm for direct / indirect restoration, and below 5,000 rpm for natural tooth. Wear protective eye glasses during use.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Wear plastic or rubber gloves during use to prevent possible sensitivity.
4. When pulling the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of dust particles.
5. The dispensing tip cannot be sterilized in an autoclave or chemiclave.
6. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
8. In case of inhalation, take the affected person into fresh air and keep them calm.
9. In case of swallowing, rinse out mouth and then drink plenty of water. If symptoms persist, consult a doctor.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
or for The Americas:
https://www.gc.dental/america
They can also be obtained from your supplier.

PACKAGE
Kit contains: DIAPOLISHER PASTE 2 g x 1,
Dispensing tip (needle type) x 1, Protective cap x 1
Refill: Dispensing tip (needle type) 20 pieces

REPORTING OBLIGATION ACCORDING TO EU MDR
In the case of serious incidents related to the device, please report it to GC Europe N.V. (vigilance@gc.dental) and to the competent authority of the Member State in which the user/patient resides.

Note:
The current Instructions for Use are available in the download section of the https://www.gc.dental/ website. To access the appropriate version, please select your region from the footer.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell Buckinghamshire
MK16 8JS United Kingdom

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