Prior to use, carefully read the instructions for use. EN

GC Mi Paste Plus™

TOPICAL CREME WITH BIO-AVAILABLE CALCIUM, PHOSPHATE AND FLUORIDE

DESCRIPTION
GC Mi Paste Plus is a creme enriched with fluoride (900 ppm F) and RECALDENT™ (CPP-ACP: Casein Phosphopeptide - Amorphous Calcium Phosphate), a special milk-derived protein. RECALDENT delivers calcium and phosphate – the essential minerals that make up tooth enamel – in a soluble form. This allows GC Mi Paste Plus to offer additional protection for teeth, neutralise acids produced by bacteria in dental plaque, reduce tooth sensitivity and protect teeth from the effects of acidic foods and drinks. Contains non-cariogenic Xylitol (sugar free, fluoride content <0.05% after swallowing (<0.5 mg)). "RECALDENT and RECALDENT logo are trademarks of Mondelēz International Group used under license."

INTENDED FOR
GC Mi Paste Plus is intended to be used as protection for teeth by:
- neutralising acid attacks from dental plaque and intrinsic sources (e. g. bacteria in dental plaque, gastric reflux) and/or extrinsic sources (e. g. foods and drinks)
- reducing tooth hypersensitivity
- preventing tooth demineralization (mineral loss) and promoting tooth remineralization (mineral gain), thereby supporting the growth of white spot lesions
- supporting the growth of beneficial oral commensal bacteria, maintaining a healthy oral microbiological balance

GC Mi Paste Plus may be used as an alternative method for topical fluoride application in children aged 6 years and above.

INDICATIONS FOR USE
For natural teeth desiring enhanced protection, prevention, and stronger enamel
Particularly for these teeth:
1. Teeth prone to erosion
2. Teeth subjected to acid challenges from acid-producing bacteria in plaque.
3. Teeth subjected to acid challenges from gastric reflux and external acid sources.
4. Hypersensitive teeth
5. Teeth with developmental defects, such as molar (incisor) hypomineralization (‚chalky teeth, cheese molars')
6. Teeth with exposed root surfaces
7. Maintaining teeth and gum health

CLINICAL BENEFITS
1. Reduction of tooth hypersensitivity
2. Prevention of initial lesions (initial caries) and tooth lesions from various acid sources (erosion)
3. Reduction of white spot lesions
4. Inhibits and suppresses the growth of harmful bacteria
5. Helps maintain a healthier oral microbiome, reducing overall harmful bacteria levels, thereby supporting a healthy oral balance for teeth and gums

CONTRAINDICATIONS
Do not use this product if you have a proven or suspected milk protein (casein) allergy, and/or with a sensitivity to any of the ingredients, soybean products or flavoring agents.

INTENDED USERS/SPECIAL TRAINING
Adult users and professionals. This product may be used on children; however, it must be applied by an adult or under adult supervision. No special training is required for the use of this product.

TARGET GROUP
All ages and groups starting from 6 years and older, particularly in these situations:
- Prior to and following tooth bleaching at home or in the dental practice.
- After professional tooth cleaning (scaling and root planing).
- During and after orthodontic treatment.
- After tooth whitening.
- If you suffer from dry mouth (e.g. Xerostomia/Sjögren's syndrome).
- If you are at risk of dental erosion (e.g. due to frequent acidic foods (erosion).
- If you have developmental defects of the tooth enamel such as molar (incisor) hypomineralisation (‚chalky teeth, cheese molars').
- If you have exposed root surfaces
- If you suffer from tooth sensitivity
- If you are prone to gum disease

COMPOSITION
RECALDENT™ (Casein phosphopeptide - amorphous calcium phosphate, 5%-15%), water, glycerine, d-sorbitol, carboxymethyl cellulose, propylene glycol, silicon dioxide, titanium dioxide, hydroxyethyl cellulose, flavor, paraben, phosphoric acid, sodium fluoride, sodium saccharin

PERFORMANCE CHARACTERISTICS
The product's performance characteristics satisfy the requirements of the intended use.

DIRECTION FOR USE
Talk to your dentist/hygienist how often and how long you should use GC Mi Paste Plus. GC Mi Paste Plus is typically recommended for use twice daily, ideally in the morning and at night, after brushing your teeth. For more effective results:

Application without tray (e.g. at home maintenance):
1. Wash hands, if possible, dry them.
2. Apply a sufficient amount of GC Mi Paste Plus to upper and lower teeth. A pea size amount for each quadrant is the minimum recommended. The material should be applied to tooth surfaces using a clean dry finger or cotton tip. If applied by a medical/dental professional, a brush may be used to apply a thicker layer of material between (between the teeth) use an interproximal tooth cleaning brush or dental floss coated with GC Mi Paste Plus.
3. Leave GC Mi Paste Plus on the tooth/teeth for a minimum of 3 minutes.
4. Then use your tongue to distribute GC Mi Paste Plus. GC Mi Paste Plus is also suitable for use in toothpaste dispensers. Apply as a ribbon (between 1 - 2 minutes) avoiding spitting out and delaying swallowing. The longer GC Mi Paste Plus remains in the mouth, the better the results. Do not rinse.
5. Spit out and do not rinse. Any GC Mi Paste Plus remaining on the surface can be left to gradually dissipate. Do not eat or drink for 30 minutes following application.

Note:
Duration of application:
Saliva will enhance the effectiveness of RECALDENT (CPP-ACP) and the flavor will help stimulate saliva flow. The longer RECALDENT™ (CPP-ACP) and saliva are in contact with the tooth surface, the more beneficial the result. For nighttime application, GC Mi Paste Plus is the after brushing your teeth and avoid eating or drinking after application and before going to bed. GC Mi Paste Plus is also suitable for use in custom fit trays. "GC Tech Mousse™ Instead for overnight use in children under 12."

Custom tray application (e.g. at the dental office / by a bleaching tray can also be used):
1. Rinse the custom tray under running water before use.
2. Extrude a sufficient amount of GC Mi Paste Plus into the tray and apply to the upper and / or lower teeth.
3. Leave the tray undisturbed in the mouth for a minimum of 3 minutes.
4. Remove the tray.
5. The product should be spread using the tongue to spread the remaining GC Mi Paste Plus evenly over the teeth. Alternatively, leave the tray in the mouth for 1 - 2 minutes avoiding expectoration and delaying swallowing. The longer GC Mi Paste Plus remains in the mouth, the better the results. Do not rinse after the treatment.
6. The product is spit out and rinsing is avoided afterwards. Any GC Mi Paste Plus remaining on the surface can be left to gradually dissipate. Do not eat or drink for 30 minutes after the patient has left the practice. Do not eat or drink for 30 minutes following application.
7. Any material left in the GC Mi Paste Plus tray should be rinsed or brushed off under running water immediately after use.

DISPOSAL
Dispose of all waste according to local regulations.

STORAGE
Store at room temperature, store in a cool and dark place (4-25°C / 39-77°F).

SIDE EFFECTS: There are no known side effects to date.

RESIDUAL RISKS
There are no known residual risks to date, other than those specified as part of contraindications, warnings, precautions, and undesirable effects.

PRECAUTIONS, WARNINGS AND LIMITATIONS OF USE
Do not use if the packaging or application technique cannot be applied.
- Always replace cap after use and remove any residual paste or moisture around the tube neck to avoid the cap sticking.
- It is not recommended to use GC Mi Paste Plus on patients under the age of 6.
- The overnight application is not recommended for children under 12.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Contains Xylitol: Can cause GI distress, avoid in pets.
- If the product's odor, taste or appearance changes after opening, discontinue the use of the product and contact the distributor from whom it was purchased.
- GC Mi Paste Plus must always be used under adult supervision.
- If any unexpected side effects occur during product application, this may indicate sensitivity or allergy to the paraben, soybean products, flavoring agents, casein, or to some other ingredient of the product. In this event, discontinue the use of the product and refer to a physician.

Some products referenced in this Instructions for Use may be classified hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at: http:// www.gcdental.de/sdb. They can also be obtained from your supplier.

FLAVORS
Strawberry, Melon, Vanilla, Mint, tutti-frutti

REPORTING OBLIGATION ACCORDING TO EU MDR
In the case of serious incidents related to the device, please report it to GC Europe N.V. (vigilance@gc.dental) and to the competent authority of the Member State in which the user/patient resides.

NOTE: The current Instructions for Use are available in the download section of the https://www.gc.dental/ website. To access the appropriate version, please enter the GC Mi Paste Plus.
Last revised 05/2026

16
INSTRUCTIONS FOR USE (IFU)
GC Mi Paste Plus™
30002596
Rev 05/2026

MANUFACTURED BY
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
DISTRIBUTED BY
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

EU: GC EUROPE N.V.
Interleuvenlaan 33, 3001 Leuven, Belgium
TEL: +32 16 74 01 00

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