GC RELINE™

Prior to use, carefully read the instructions for use. EN

CHAIRSIDE HARD DENTURE LINER

DESCRIPTION
GC RELINE is a chairside hard denture liner to be used for hard denture linings and/or extension of denture bases.

INTENDED PURPOSE
GC RELINE is intended to be used mainly for lining of acrylic dentures when their fitting onto the dental arches is no longer acceptable from the clinical point of view.

INDICATIONS FOR USE
Improper fit of acrylic dentures.
Areas of application:
1. Chairside hard denture linings
2. Extension of denture borders

CLINICAL BENEFITS
Used to help minimize irritation of oral tissue due to ill-fitting dentures.


CONTRAINDICATIONS
Do not use on patients with a known sensitivity to methacrylate. In case of allergy, consult a physician.

INTENDED USERS/SPECIAL TRAINING
- Dental professionals, such as general dentists, dental assistants, and hygienists, who possess both theoretical and practical knowledge regarding the use of dental products.
- No special training required.

THE PATIENT TARGET GROUPS
Target populations are edentulous or partially edentulous patients in need of relining for acrylic dentures or extension of denture borders.

PERFORMANCE CHARACTERISTICS
The product's performance characteristics satisfy the requirements of the intended use.

DIRECTION FOR USE
1. Roughen the area of denture to be relined. Clean and dry the denture thoroughly.
2. Coat the artificial teeth, labial and buccal surface of the denture with a blocking barrier using a cotton swab. Do not apply coating within 3mm (1/8 inch) of the peripheral border. For dentures with plastic teeth, apply a coating of a bonding barrier to the teeth to prevent GC RELINE material adhering to the tooth surface.

3. Coat the denture to be relined with BONDING AGENT using the enclosed brush. Air dry the solvent for approximately 10 seconds.
Note: BONDING AGENT is volatile, be sure to close cap tightly after using.
4. Standard consistency is obtained by mixing the ratio of 1.8 g powder and 1.0 mL liquid. First measure 1.8g of powder and transfer the powder into a mixing cup. Then measure 1.0mL of liquid. Dispense liquid over the powder in the mixing cup and stir gently for 15 seconds.
5. Using a plastic spatula, apply the resin mix evenly onto the tissue contact surface of the denture.
6. Insert the denture into the mouth at approximately 1 minute from start of the mix. Instruct the patient to close lightly into occlusion. Muscle trim the borders as necessary.
7. For a full denture reline without undercuts: Remove the excess portion of the resin intraorally. Keep the denture in the mouth until curing is completed (approximately 5 minutes and 30 seconds from the start of the mix).
8. For a full denture reline with undercuts or a partial denture: Remove the denture when the material reaches a rubbery state (approximately 2' to 2'30" from the start of the mix). Trim excess and reset denture in the mouth before it begins to generate heat. If denture is to remain out of the mouth, place it in warm water (approximately 40° to 50°C = 105° to 120°F) for 3 minutes until completely cured.
Note: When undercuts or retentive areas are present, the denture should be removed from the mouth while the material is in a rubber-like state and before it begins to generate heat. This will minimize the difficulty to disengage and remove the denture.
9. For denture border extensions: Mix material slightly thicker than standard powder/liquid ratio.
10. Remove the denture from the mouth at approximately 5 minutes and 30 seconds from the start of the mix. Finish and polish the denture in a usual manner.

DISPOSAL
Dispose of all waste according to local regulations.

STORAGE
Recommended for optimal performance, store in a cool and dark place (15-25°C / 59-77°F).

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS:
(1) Wipe the device (and accessories, if included) with a lint-free disposable wipe wet with a cleaning agent/alkaline and quaternary ammonium based until it is completely wet and allow to dry for the entire exposure time in accordance with the instructions of the manufacturer of the cleaning and disinfectant solution.
(2) Repeat the process a second time.
(3) Visually inspect the cleanliness.
(4) Repeat steps 1 to 3, if residues are still visible on the medical devices.
Immediately after use inspect device and label for deterioration. Discard device if damaged.

SIDE EFFECTS
There are no known side effects to date.

RESIDUAL RISKS
Inflammatory reactions of the oral tissues.

PRECAUTIONS, WARNINGS AND LIMITATIONS OF USE
Do not use if the prescribed application technique cannot be applied.
1. Liquid and BONDING AGENT liquid are volatile and flammable, avoid sources of ignition.
2. Wash the brush with GC RELINE liquid or ethanol after applying the BONDING AGENT.
3. Shake the powder bottle lightly before use.
4. In case of contact with skin, wash thoroughly with water.
5. In case of contact with eyes, flush immediately with water and seek medical attention.
6. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at: https://www.gc.dental/europe
They can also be obtained from your supplier.

PACKAGE
Powder (80g) x 1
Liquid (50mL; 49.7g) x 1
Bonding Agent (12g; 13.2mL) x 1
Mixing Cups x 4
Vial x 1
Dropper x 1
Brush x 1
Spatula x 1

REPORTING OBLIGATION ACCORDING TO EU MDR:
In the case of serious incidents related to the device, please report it to GC Europe N.V. (vigilance@gc.dental) and to the competent authority of the Member State in which the user/patient resides.

Note: The current Instructions for Use are available in the download section of the https://www.gc.dental/ website. To access the appropriate version, please select your region from the footer.

Last revised: 11/2025

MANUFACTURED BY
GC AMERICA INC.
Alsip, IL 60803, U.S.A.
TEL: +1-708-597-0900
www.gc.dental/america

DISTRIBUTED BY
EU : GC EUROPE N.V.
Interleuvenlaan 33, 3001 Leuven, Belgium
TEL: +32 16 74 10 00

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