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Prior to use, carefully read the instructions for use. EN

TEMPSMART™ DC

DUAL-CURED COMPOSITE FOR TEMPORARY CROWN AND BRIDGE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE
1. Fabrication of temporary crowns, bridges, inlays, onlays and veneers
2. Fabrication of long-term temporary restorations

CONTRAINDICATIONS
Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

COMPOSITION
Paste A: Dimethacrylate, silicon dioxide, initiator, pigment, stabilizer
Paste B: Dimethacrylate, silicon dioxide, initiator, stabilizer

DIRECTIONS FOR USE
Fabricating temporary restorations
1. Taking impression
 1) Take the impression using alginate or silicone impression material. Alternatively, temporary can be made using a laboratory-prepared thermoforming matrix.
 2) Block out preparation undercuts.
2. Tooth preparation
 1) Block out preparation undercuts.
 2) Clean tooth preparation with water spray and lightly dry.
 3) Lightly lubricate with separating medium, such as petroleum jelly on prepared teeth or resin core.
3. Dispensing
 NOTE:
 Prior to first use of a new syringe or after a long interval between uses, extrude paste from the syringe onto a mixing pad.
 1) Remove cap of syringe by rotating 1/4 turn counter-clockwise. Attach a mixing tip by aligning V shaped notch on rim of the mixing tip with V shaped notch between syringe barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of mixing tip 1/4 turn clockwise.
 2) Discard first centimeter of extruded material before use.
 3) Extrude material directly into impression.
 NOTE:
 Used mixing tip can be left on syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
4. Fabricating a temporary restoration
 Fabricate a temporary or provisional bridge material according to timing guidelines. The indicated starts after mixing paste:

 Option 1: Fabricating in the mouth

0'00"–0'45"	Place TEMPSMART DC into impression and reposition in mouth.
0'45"–2'00"	Setting in mouth.
2'00"–2'30"	Remove impression from teeth. If necessary, tack cure each tooth unit with curing light for 1–2 seconds.
Final cure	Remove temporary restoration from mouth and cure light for each tooth unit with light curing unit (High power LED (more than 1200 mW/cm²), Halogen / LED (700-1200 mW/cm²)) for 5 seconds. Irradiate occlusal surface. If chemical cure only, keep temporary restoration in 23°C/74°F until 5'00".
After final cure	Completely remove oxygen inhibition layer from surface of temporary restoration with alcohol.

 Option 2: Fabricating on the model

0'00"–0'45"	Place TEMPSMART DC into impression and reposition on model.
0'45"–3'00"	Setting on model.
Final cure	Remove impression from the model and final cure with an extra oral light curing oven for 1 minute. Remove temporary restoration from model. If chemical cure only, keep temporary restoration in 23°C/74°F until 5'00".
After final cure	Completely remove oxygen inhibition layer from surface of temporary restoration with alcohol.

 CAUTION
 Due to high polymerization temperature created, do not irradiate temporary restoration for more than 2 seconds. Do not final cure on tooth.
 NOTE:
 1) Above processing times are applicable for products stored and processed at a temperature of 23°C/74°F and at 50% relative humidity. Higher temperatures shorten these times, while lower temperatures extend them. The size of the restoration also has impact on setting times: as size of restoration increases, setting takes longer.
 2) When final curing, it is important to cure from the occlusal surface to prevent deformation. In case of bridge, irradiate both ends of the bridge alternately.
5. Finishing
 1) Remove excess material using fine carbide bur.
 2) Check occlusion and if necessary correct shape by using a resin material such as TEMPSMART DC, G-aenial Flo X, G-aenial Universal Injectable or UNIFAST Trad for custom-shaping and shape correction, usually there is no need for a bonding agent. Completely remove saliva contaminated part using a carbide bur. If necessary, apply a ceramic primer such as CERAMIC PRIMER II or G-Multi PRIMER to surface of temporary restoration.
 (2) When using TEMPSMART DC for custom-shaping and shape correction, light cure mixing paste with curing light for 10 seconds.
6. Cementing
 Cement using temporary cement and follow manufacturer's instructions.
 NOTE:
 Do not use eugenol cement as it will prevent polymerization of resin.
7. Repairing / Modifying temporary restorations
 1) Abrade surface of temporary restorations using carbide bur.
 2) If necessary, apply a ceramic primer such as CERAMIC PRIMER II or G-Multi PRIMER and dry with air.
 3) Place a resin material such as TEMPSMART DC, G-aenial Flo X, G-aenial Universal Injectable or UNIFAST Trad and cure according to manufacturer's instructions.
 NOTE:
 Completely remove saliva contaminated part using a carbide bur.

SHADES
A1, A2, A3, A3.5, B1 and BW

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

PACKAGES
Refills:
1. Cartridge 48 mL (60 g) (1), MIXING TIP II SSS (16)
2. Cartridge 10 mL (12.5 g) (1), GC Automix Tip Regular (10)

Mixing tip package:
1. MIXING TIP II SSS (60)
2. GC Automix Tip Regular (20)

Dispenser package:
1. CARTRIDGE DISPENSER II (1)
2. GRADIA CORE DISPENSER GUN (1)

CAUTION
1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate operating field from skin or oral tissue.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Take care to avoid ingesting material.
4. Do not mix with components of other products.
5. Mixing tips are for single use only. To prevent cross contamination between patients, do not reuse the mixing tip. The mixing tip cannot be sterilized in an autoclave or chemiclave.
6. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Wearing gloves during use avoids direct contact with air inhibited resin in order to prevent possible sensitivity.
7. Wear protective eye glasses during light curing.
8. When polishing polymerized material, use dust collector and wear a dust mask to avoid inhalation of dust.
9. Avoid getting material on clothing.
10. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
11. Do not use TEMPSMART DC in combination with eugenol containing materials as eugenol may delay or prevent TEMPSMART DC from setting.
12. Do not wipe mixing tip with excessive alcohol as this may cause contamination of the paste.
13. Dispose of all wastes according to local regulations.
14. In rare cases product may cause sensitivity in some people. If any such reactions are experienced, discontinue use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
https://www.gc.dental/europe
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last Revised: 12/2024

CE 2797

MANUFACTURED by:
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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