Prior to use, carefully read the instructions for use. EN

REVISED SIZE FOR TEST - 4x1½

GC Tooth Mousse™

TOPICAL CREME WITH BIO-AVAILABLE CALCIUM AND PHOSPHATE

For use by or under the supervision of a dental professional in the recommended indications.

Description:
GC Tooth Mousse is a water based creme containing RECALDENT™ (CPP-ACP: Casein Phosphopeptide – Amorphous Calcium Phosphate). When CPP-ACP is applied in the oral environment, it will bind to biofilms, plaque, bacteria, hydroxyapatite and soft tissue localising bio-available calcium and phosphate ions.
Important:
Xylitol will enhance the effectiveness of CPP-ACP and the flavour will help stimulate saliva flow. The longer CPP-ACP and active are maintained in the mouth, the more effective the result.
A topical paste containing bio-available calcium and phosphate
1. Provides extra protection to the teeth.
2. Helps neutralize acid challenges from acidogenic bacteria in plaque.
3. Helps neutralize acid challenges from other internal and external acid sources.

CONTRAINDICATIONS
RECALDENT™ (CPP-ACP) is derived from milk casein. Do not use this material on patients with a proven or suspected milk protein allergy and/or with a sensitivity or allergy to benzoate preservatives.
DIRECTION FOR USE
INDICATIONS
a. Prior to and following in-office bleaching.
b. After ulcer therapy, hand scaling or root planing.
c. Following professional tooth cleaning (P.T.C.).
d. After application of topical fluorides.
e. To provide a topical coating for patients suffering from erosion, xerostomia or Sjögrens syndrome.
f. During orthodontic treatment.
g. For high risk caries patients.
h. For special needs patients.

In office application:
I. Custom tray application
1. Before use, make a custom tray thoroughly under running water.
2. Extrude a generous layer of GC Tooth Mousse into the tray and apply to the upper and / or lower teeth.
3. Leave the tray undisturbed in the mouth for a minimum of 3 minutes.
4. Remove tray.
5. Instruct the patient to use the tongue to spread the remaining GC Tooth Mousse throughout the mouth. Instruct the patient to retain for as long as possible (1+ 2minutes) avoiding expectoration and deglug swallowing. The longer GC Tooth Mousse and saliva are maintained in the mouth, the more effective the result.
6. Ask the patient to expectorate and if possible avoid rinsing. Any GC Tooth Mousse remaining on the surface can be left to gradually dissipate. Advise the patient not to eat or drink for 30 minutes following application.
7. Any excess GC Tooth Mousse in the tray should be rinsed or brushed off under running water immediately after use.

II. Non tray application
1. If necessary, remove any excess saliva on the tooth surface with a cotton roll. However, it is not necessary to dry the teeth with compressed air.
2. Using sufficient amount of GC Tooth Mousse to tooth surfaces using an application swab, gloved finger or in difficult interproximal areas using an Interproximal Tooth Cleaning Brush.
3. Leave GC Tooth Mousse undisturbed for a minimum of 3 minutes.
4. Then instruct the patient to use the tongue to spread the remaining GC Tooth Mousse throughout the mouth. Instruct the patient to hold in the mouth for as long as possible (an additional 1 – 2 minutes) avoiding expectoration (spitting out) and deglug swallowing. The longer GC Tooth Mousse and saliva are maintained in the mouth, the more effective the result.
5. Ask the patient to expectorate and if possible avoid rinsing. Any GC Tooth Mousse remaining on the surface of tooth should be left to gradually dissipate. Advise the patient not to eat or drink for 30 minutes following application.

IV. Night time application after tooth brushing as recommended by a dental professional
1. Apply a sufficient amount of GC Tooth Mousse to the upper and lower teeth. A pea size amount for each arch is the minimum required. The material should be applied to the tooth surfaces using a clean dry finger or cotton tip. For difficult areas (between the teeth) use an Interproximal Tooth Cleaning Brush or dental floss coated with GC Tooth Mousse.
2. Leave GC Tooth Mousse on the teeth undisturbed for a minimum of 3 minutes.
3. Then use your tongue to spread the remaining GC Tooth Mousse throughout the mouth. Hold in the mouth for as long as possible (a further 1 – 2 minutes) avoiding expectoration (spitting out) and deglug swallowing. The longer GC Tooth Mousse and saliva are maintained in the mouth, the more effective the result.
4. Expectorate thoroughly and if possible avoid rinsing. Any GC Tooth Mousse remaining on the surface can be left to gradually dissipate and/or rubbed in. Do not eat or drink for 30 minutes following application.

FLAVORS
Strawberry, Melon, Vanilla, Mint, Tutti-Frutti
STORAGE
Recommended for optimal performance, store in a cool dry low humidity area away from direct sunlight (8-25°C; 46,4-77,0°F).
PACKAGES
Tube containing 40g (35mL)
Assorted Package
- 10 x 40g (2 each of Strawberry, Melon, Vanilla, Mint, Tutti-Frutti)

2. Single Flavor Package
- 10 x 40g in each flavor
Vanilla, Strawberry, Mellon, Tutti-Frutti, Mint
3. Refill Package
- 10 x 40g tube in each flavor

NOTE
GC Tooth Mousse with its unique CPP-ACP formulation was developed at The School of Dental Science, The University of Melbourne, Victoria, Australia. *RECALDENT is a trademark used under licence. CPP-ACP technology has related patents or patents pending in Australia, NZ, Europe, Canada, and USA.

CAUTION
1. Always replace cap after use and remove any residual paste or moisture around the tube neck or inside the cap.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. In case of spillage on clothes, wash off with water.
4. If any ampicillin symptoms are experienced, this may indicate sensitivity or allergy to the benzoate preservatives, or to some other components of the product such as a flavoring agent. In this event, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at:
http://www.gceurope.com/
They can be obtained from your supplier.
Last revised : 06/2025

MANUFACTURED BY:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
DISTRIBUTED BY:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
GC CORPORATION
76-1 Hasuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
GC EUROPE N.V.
Researchpark Haasrode 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00
GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01, Tampines Plaza 2, Singapore 529541
TEL: +65-6546 7588
GC AUSTRALASIA DENTAL PTY. LTD.
17/12 Rodborough Rd, Frenchs Forest, NSW 2019, Australia
TEL: +61 2 9301 8200
PRINTED IN THE USA

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