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EN Calcium hydroxide paste for temporary root canal filling

BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY
TO BE USED BY DENTIST ONLY

COMPOSITION: Calcium hydroxide, hydroxyapatite, barium sulphate, excipients.

INTENDED USE
CALCIPAST is intended for use in dental treatment for temporary root canal filling.
Additionally, CALCIPAST may be used for temporary filling of deep cavities during indirect or direct pulp capping.
It has remineralizing effect on dental tissues thanks to its high content of calcium hydroxide. The product supports restoration of damaged periapical tissues.
pH 12.5–13. Radiopaque, non-hardening material.

INSTRUCTIONS FOR USE
Before treatment, dip the provided applicator in ethanol, allow to dry and then attach it onto the syringe.
Before use, check if the applicator is correctly attached onto the syringe and then expel a small amount from the syringe to check the paste flow through the applicator.
The syringe grip provides a comfortable grip for easy dispensing management, reducing the risk of pushing the product through the root apex and minimising the difficulties associated with expelling the product.
For MEGA PACK, manually put the grip onto the syringe before use. The grip plate should be put onto the capped syringe, as shown on the product packaging (a distinct „click" is heard when correctly applied).
When used directly in the tooth canal: measure the length of the canal, adjust the endostop on the applicator to a length 2-3 mm shorter than the measured length of the canal, dispense the paste in the canal, remove excess paste, close the cavity tightly. After the paste has been applied, an X-ray should be taken to check that the paste has adequately filled the canal and has not gone beyond the apex.
When used in a cavity: dispense the paste at the bottom of the cavity, close the cavity tightly.
Leave the filling in the cavity or root canal for 1 week to 30 days.
The product is intended for multiple use, while the applicator is for single use only. Reusing the applicator may pose a risk of secondary infection .
The grip plate is for multiple use.

CONTRAINDICATIONS
Do not use for patients who are hypersensitive to the product components.

WARNINGS AND PRECAUTIONS
Irritant to oral mucosa and skin. Causes serious eye damage.
In case of contact, rinse immediately with plenty of water and seek medical help.
If ingested, do not induce vomiting. Drink plenty of water. Seek medical help.
Use dental dam, disposable gloves, safety goggles/visors and protective clothing when handling the product.
Endodontic treatment may involve temporary pain and postoperative tenderness.
The paste should be carefully removed from the root canal before the final filling, as the presence of the formulation may affect the bond strength of dentin to the tooth restorative materials. The use of temporary materials may be associated with bacterial microleakage. CALCIPAST has cytotoxic potential.
Do not expel the paste if the applicator tip is wedged inside the canal. The applicator must be loosely positioned in a canal.
Before using CALCIPAST, the working length of the root canal should be determined. Excessive amounts of temporary material used to fill the root canal may lead to it being pushed beyond the canal apex. Particular care should be taken in anatomically sensitive areas such as the mandibular canal, mandibular opening or maxillary sinuses.
Keep out of reach of children.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Keep in original packaging below 25°C.
Expiry date is placed on the main container. Once container is open, the expiry date does not change if the container is securely closed after each use.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
CALCIPAST: Syringe containing 2.1 g of preparation, disposable applicators (2 pieces), grip plate.
CALCIPAST MEGA PACK: 4 syringes of preparation (each containing 2.1g of paste), disposable applicators (4 pieces), grip plate.

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