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EN Temporary filling

BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY
TO BE USED BY DENTIST ONLY

COMPOSITION: zinc oxide, zinc sulfate, calcium sulfate, excipients.

INTENDED USE
CAVIPACK is intended for temporary filling of cavities and for temporary sealing of root canal orifices during endodontic procedures.
CAVIPACK is a ready-to-use paste. No preliminary preparation is required. **It does not contain eugenol.** Binds chemically when exposed to moisture (saliva). Use for class I and II cavities.

INSTRUCTIONS FOR USE
Spray the cavity with water to moisten. Place the material in the cavity with an appropriate instrument, and condense it by pressing gently against the cavity walls. Shape the filling to fit the occlusion.
When sealing the mouths of root canals: measure the length of the canal, introduce the paste into the canal, remove excess paste from the canal, seal the cavity.
The surface of the material hardens within 20-30 minutes when exposed to moisture (saliva). Avoid biting for 2-3 hours. Remove from the cavity with a drill or scaler, rinse the remains with water spray.

CONTRAINDICATIONS
Do not use the product in patients with hypersensitivity to the preparation ingredients. Do not use as a temporary filling for multi-surface defects.

PRECAUTIONS
The use of temporary materials may be associated with bacterial microleakage. Overfilling of the root canal can cause the tooth to crack or break. Causes serious eye damage. In case of contact, rinse immediately with plenty of water and seek medical attention. It can cause allergic skin reaction.
If swallowed, drink plenty of water, do not induce vomiting. Contact a physician for any symptoms. Use a cofferdam, disposable gloves, safety goggles/visors and protective clothing when
handling the product. Keep out of reach of children.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Store in the original package below 25°C. Protect against moisture. After each use, close the packaging tightly. Expiry date is placed on the main container.
After opening the direct packaging, the shelf life does not change, provided that the packaging is resealed after each use.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
Packaging containing
3x10 g of the preparation.

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