**PL** Gutaperka

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**EN** Gutta percha points

## BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY
TO BE USED BY DENTIST ONLY

COMPOSITION
Gutta Percha, zinc oxide, barium sulfate, coloring agent.

INTENDED USE
Hand-rolled endodontic points. Used to fill a prepared root canal prior to tooth restoration.

INSTRUCTIONS FOR USE
Check the product if it is damaged, presence of impurities, and shape is suitable for use. Check if the product is within the expiration date. After treating the root canal of the damaged tooth, use the proper size gutta percha to fill. This product is disposable and discarded after use.

CONTRAINDICATION
Do not use in patients allergic to this product.

PRECAUTIONS
Product for dental use only. Do not use products that have expired. Check the product packaging abnormality before use. Before use, check for infection inside the root canal. Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
The expiry date is on the immediate packaging. Store at room temperature away from direct sunlight. Protect from moisture. Keep out of reach of children.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
GUTTA PERCHA POINTS 2% 4%, 6%, ProTaper - marked, unmarked.
Packaging available: slide box; endo-aide.

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