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EN Disposable dental cuspidor liners

BEFORE USE PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY

TO BE USED BY DENTIST ONLY

INTENDED USE
Disposable dental cuspidor liners SPIT SPOT are intended for use as an accessory to a medical device - a dental unit - during dental procedures.

The device is placed in the cuspidor bowl to collect contaminants generated during dental procedures, including saliva, liquids, organic debris, and small waste materials. The device supports the safe and hygienic use of the dental unit in accordance with its intended purpose and limits direct contact of personnel with potentially biologically contaminated material during emptying and cleaning of the cuspidor.
SPIT SPOT is intended for single use in dental practices by dentists. The device is not intended for direct contact with the patient's body and is supplied non-sterile.

INSTRUCTIONS FOR USE
1. Remove one liner from the packaging.
2. Place the liner in the cuspidor bowl, ensuring that it does not block the liquid drain.
3. After the procedure, grasp the liner using the holes designed to facilitate removal and carefully remove it from the cuspidor.
4. Dispose of the used liner. Do not reuse.

PRECAUTIONS
Before use, check the condition of the packaging and the device. Do not use the device if it is damaged, soiled, or deformed. The device is intended for single use only. Do not reuse.
After use, the device may contain potentially biologically contaminated material. When removing the liner, use personal protective equipment in accordance with the procedures applicable in the dental practice, in particular protective gloves and, where necessary, a protective mask, safety glasses, or a face shield.
Do not leave the used liner in the cuspidor after completion of the dental procedure. Do not use the device if its size or shape prevents correct placement in the cuspidor bowl or interferes with the drainage of liquids from the dental unit.

Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Protect from moisture.

DISPOSAL
After use, the device may be potentially contaminated with biological material or other contaminants generated during the dental procedure. The used device should be disposed of in accordance with the procedures for handling medical or potentially infectious waste adopted by the healthcare facility and in compliance with local regulations.

PACKAGING
The packaging contains 50 pcs.

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