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**EN** Instrument for placing the retraction cord in gingival fissure

## BEFORE USE, PLEASE READ THE FOLLOWING INSTRUCTION CAREFULLY

FOR DENTAL USE ONLY

NON-STERILE - STERILIZE DIRECTLY BEFORE USE

INTENDED USE:
Instrument for placing the retraction cord in gingival fissure.

HANDLING DENTAL TOOLS MADE OF STAINLESS STEEL
The instrument is made of stainless steel. Special processes have been undertaken during the manufacture of the dental instruments to ensure corrosion-resistance. In spite of this, the durability of these dental instruments depends on their proper use and correct maintenance. Before each use, the instrument should be inspected for contamination, damage to the material structure (cracks, bends, splits, scratches – deeper scratches). Failure to adhere to the above stated instructions decreases instrument durability and reduces instrument lifetime.

CLEANING, RINSING AND WASHING
The packaged instruments are non – sterile. Instruments that are brand new should be thoroughly cleaned with a non-scratching soft bristle brush and an instrument cleaning solution (detergent: ANIOSYME X3; solution prepared according to the manufacturer's instructions). Ultrasonic cleaning is also recommended. After cleaning, all components should be rinsed with purified water and dried, as moisture can adversely affect the effects of sterilisation. Instruments must be cleaned, rinsed and dried again and sterilised before each use on a subsequent patient.

STERILIZATION
All (reusable) dental instruments must be disinfected, then washed out and then dried before sterilisation.
Instruments should be sterilized in accordance with national requirements.

Recommended sterilisation methods:

| Method | Minimal parameters: |
|---|---|
| Steam autoclave | 134°C, for 5 minutes |
| | 121°C, for 20 minutes |

**PRECAUTIONS**
Adhere to the guidelines for cleaning and sterilization of the dental instruments.

**STORAGE**
Dental instruments should be stored in temperatures ranging between 5-30°C and relative humidity below 70%.
Failure to adhere to the above stated instructions decreases instrument durability and reduces instrument lifetime.

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