Dentsply Sirona

Unimetric™

EN	Instructions For Use	2
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	22
DE	Gebrauchsanweisung	32
ES	Instrucciones de uso	42
PT	Instruções de utilização	52
NL	Gebruiksinstructies	62
SV	Användarinstruktioner	72
DA	Brugsanvisning	82
NO	Brukerinstruksjoner	92
LT	Naudojimo instrukcija	101
LV	Lietošanas norādījumi	111
ET	Kasutusjuhised	121
CS	Návod k použití	131
SK	Návod na použitie	141
HU	Használati utasítás	151
RO	Instrucţiuni de utilizare	161
SR	Uputstva za upotrebu	171
BS	Uputstva za upotrebu	181
CNR	Uputstva za upotrebu	191
SQ	Udhëzimet për përdorim	201
EL	Οδηγίες χρήσης	211
TR	Kullanım Kılavuzu	221
PT-BR	Instruções de uso	231
KO	사용 지침	241

Dentsply Sirona

EN

Unimetric™

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
NON-STERILE – MANDATORY PROCESSING BEFORE FIRST USE
SINGLE USE – REPROCESSING PROHIBITED
REUSABLE – LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

POSTS: SINGLE USE NON-STERILE IMPLANTABLE ENDODONTIC DEVICES
REF. C215T (Unimetric™ X.X mm)
DRILLS FOR POSTS: REUSABLE NON-STERILE ENDODONTIC INSTRUMENTS
REF. C0212 (Unimetric™ X mm - Special Penetration Drill),
C0213 (Unimetric™ X mm - Special Calibrating Drill),
C0214 (Unimetric™ X mm - Special Precalibrating Drill)
POSTS SETS[1]
REF. C0321 (Unimetric™ T Demo Kit),
C226A (Unimetric™ X mm Integral Set),
C226T (Unimetric™ X mm Set)
[1.] Posts and Drills for Posts instruments can be sold in sets with other instruments and components. See sections [4] and [17] for details.

1 PRODUCT DESCRIPTION

- Product description:

 Unimetric™ Root Post
 Unimetric™ Root Post are single use non-sterile Titanium Post. Unimetric™ Root Post are conical and are available in different sizes (see section [4]).

 Unimetric™ - Drills for Posts
 Unimetric™ - Drills for Posts are Stainless steel instruments made of a working part for precision drilling, a shank to be connected to a contra-angle, and a color coding to identify the corresponding post.

- Clinical benefits:
 The clinical benefit for the patient being treated with Unimetric™ is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.

- Intended purposes:
 - Posts are devices intended to retain the core material in a tooth with extensive loss of coronal structure.
 - Drills for Posts instruments are devices used for the creation of hole for post placement.
 - Post sets are devices intended to place a post for post-endodontic restoration.
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

Insufficient residual dentin: at least 2 mm of tooth structure is required around the preparation.

4 DELIVERY FORM

The Unimetric™ brand comprises the following devices:

Posts

SKU	Product	Composition	Available sizes						Packaging and content
			N°X	Color	Anchor		Head		
					Length (mm)	Diameter (mm)	Length (mm)	Diameter (mm)	
C215TL00XXX00	Unimetric™ L X.X mm	Titanium Grade 5 (DIN 3.7165)	N°108L	White	7.5	0.8	2.5	1.51	Plastic bag with 25 devices
			N°110L	Blue	7.5	1	2.5	1.71	
			N°208L	Yellow	9.25	0.8	2.85	1.68	
			N°210L	Green	9.3	1	2.8	1.88	
			N°308L	Red	11.00	0.8	3.2	1.85	
			N°310L	Black	11.00	1	3.2	2.05	
C215TS00XXX00	Unimetric™ S X.X mm	Titanium Grade 5 (DIN 3.7165)	N°008	Purple	6.25	0.8	2	1.39	
			N°108S	White	7.5	0.8	2	1.51	
			N°110S	Blue	7.5	1	2	1.71	
			N°208S	Yellow	9.25	0.8	2.25	1.68	
			N°210S	Green	9.3	1	2.2	1.88	
			N°308S	Red	11.00	0.8	2.5	1.85	
			N°310S	Black	11.00	1	2.5	2.05	

Drills for Posts

SKU	Product	Available sizes				Packaging and content
		N°X	Color	Ø	Working part length	
C0212200XXX00	Unimetric™ X mm Special Penetration Drill	N°008	Purple	Ø 0.80 mm	7.00 mm	Box with 6 instruments
		N°108	White	Ø 0.80 mm	8.25 mm	
		N°208	Yellow	Ø 0.80 mm	10.00 mm	
		N°308	Red	Ø 0.80 mm	11.75 mm	
		N°110	Blue	Ø 1.00 mm	8.45 mm	
		N°210	Green	Ø 1.00 mm	10.20 mm	
		N°310	Black	Ø 1.00 mm	11.95 mm	
C0213200XXX00	Unimetric™ X mm Special Calibrating Drill	N°008	Purple	Ø 0.78 mm	6.25 mm	
		N°108	White	Ø 0.78 mm	7.50 mm	
		N°208	Yellow	Ø 0.78 mm	9.25 mm	
		N°308	Red	Ø 0.78 mm	11.00 mm	
		N°110	Blue	Ø 0.98 mm	7.50 mm	
		N°210	Green	Ø 0.98 mm	9.25 mm	
		N°310	Black	Ø 0.98 mm	11.00 mm	
C021420000000	Unimetric™ X.X mm Special Precalibrating Drill	/		Ø 0.78 mm	11.00 mm	

Posts Sets

- The Unimetric™ instruments can also be sold in sets with the Unimetric™ Wrench component REF. C0201. The different configurations of sets are detailed in section [17]).

5 COMPOSITION

Unimetric™ Root Post

Unimetric™ Root Post are made of Titanium Grade 5.

Unimetric™ - Drills for Posts and Wrench

Unimetric™ - Drills for Posts and Wrench are made of Stainless steel.

6 COMPATIBLE DEVICES

- Posts can be used in combination with primer adhesive combined with dual cure resin cement, and the devices below.
- REF. C0214, C0213 and C0212 shall be used in combination with a contra-angle.
 We recommend using REF. C0214, C0213, C0212 at a constant rotation speed within the 800 - 1200 rpm range.

7 PERFORMANCE CHARACTERISTICS

- Expected life of Posts: 10 years.
- Posts are radio opaque.

8 SAFETY NOTES

8.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the Unimetric™ devices, the patient and/or the user:

- Breakage of instrument.
- Contamination and cross-contamination.
- Swallowing of working part of the instrument.
- Post debonding.
- Toxic or allergic reactions caused by processing residues.

- The use of inappropriate rotation speed, pressure and contact duration would speed up the wear of drills for posts, increase the risk of breakage and heat creation.

8.2 Precautions

General

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Processing step is mandatory before use according to section [10] for all the products mentioned in this IFU.

Posts

- Check the Posts devices before use, do not use and discard the devices if damaged.
- Do not apply excessive pressure on the Posts devices when inserting in the canal using the Wrench component.
- Do not touch the Posts devices with your fingers after cleaning / sterilization.

Drills for Posts

- Check the Drills for Posts instruments before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Check the Drills for Posts instruments and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Before using any Drills for Posts instrument, make sure it is it is well connected to the contra-angle head.
- Always use minimal apical pressure. Never force the Drills for Posts instruments down the canal.

Posts Sets components

- Check the Wrench components before use, do not use and discard the products if damaged.

8.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

8.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8.5 MRI (Magnetic Resonance Imaging) safety information

Posts are MR Conditional. Unimetric™ Post devices which remain in the patient's body can be affected by MR energy. For more information refer to the Annex at the end of this IFU.

9 STEP BY STEP INSTRUCTIONS

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
2) Place the transparent drawing (with full sizes outlines of the different sized posts) over the radiography to select the appropriate Unimetric™ X.X mm size:
 - ⇨ At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
 - ⇨ Select the post size so that 2/3 of its length will be inserted in the root canal. Choose a long or a short head, according to the final level of the reconstructed tooth.
3) Select the appropriate Unimetric™ X mm Drills for Posts instruments colors / sizes according to the canal anatomy and the selected Unimetric™ X.X mm size as indicated in the following table:

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4) Remove the gutta percha and perform penetration drilling of the canal cavity with the selected Unimetric™ X mm - Penetration drill. The final position of the head is indicated by either of the two rings on the drill: one representing a long head, the other a short head post.
 Drill as deep as necessary for the post to be inserted in the root canal at 2/3 of its length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
5) If necessary, use the black universal Unimetric™ X mm - Precalibrating drill to precalibrate the root canal. While drilling, continue until resistance is felt.
6) Give the canal its final size by using the selected Unimetric™ X mm - Precision calibrating drill and the same method as for the Unimetric™ X mm - Penetration drill.
7) Rinse the canal with water or air to remove drilling residues.
8) Try each root post in its correctly sized canal, but without pushing completely in. Check that the heads of the posts cannot interfere with each other during insertion. If any post is too high, repeat operations n° 4 & n° 6.

9) For dentin conditioning prior to post cementation, refer to adhesive or cement manufacturer's instructions for use if applicable.

10) Insert the sealing adhesive or cement in the root canal and on the post itself following adhesive or cement manufacturer´s instructions for use[2].

11) Insert the selected Unimetric™ X.X mm post in the root canal using the appropriate Unimetric™ Wrench to 2/3 of its length.

12) Wait for the cement to harden, then remove any excess.

[2.] The adhesion between the Posts devices and the Prime&Bond Active with SCA and Core-x Flow has been tested.

10 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Unimetric™ posts are single use implantable devices, re-use is prohibited. Re-use can increase the risk of cross contamination, poor cementation, or breakage.
- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps. This IFU covers first use of Unimetric™ posts, and first use and re-use of Unimetric™ drills and wrench.

11 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

12 EXPIRATION DATE

None.

13 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

14 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

15 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information will be available at both the European database on medical devices (EUDAMED) and the IFU website, where it will be linked to the Basic UDI-DI.

EUDAMED

https://ec.europa.eu/tools/eudamed

IFU

dentsplysirona.com/ifu

The Basic UDI-DI numbers are:

- Posts made of Titanium: ++J00310059EE
- Post Sets: ++J00310129EA

16 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

17 POSTS SETS

The different configurations of Posts Sets are detailed in the following table:

SKU	Product	Available sizes XXX	Packaging and content
C0321000XXX00	Unimetric™ T Demo Kit	208 / 210	Box with: 1 Unimetric™ Wrench 1 Unimetric™ X mm - Special Penetration Drill N°XXX 1 Unimetric™ X mm - Special Calibrating Drill N°XXX 2 Unimetric™ X.X mm N°208 S 2 Unimetric™ X.X mm N°208 L
C226AT00XXX00	Unimetric™ X mm Integral Set	110 / 210	Box with: 1 Unimetric™ Wrench 1 Unimetric™ X mm - Special Penetration Drill N°XXX 1 Unimetric™ X mm - Special Calibrating Drill N°XXX 1 Unimetric™ X mm - Special Precalibrating Drill 5 Unimetric™ X.X mm N°XXX L 15 Unimetric™ X.X mm N°XXX S
C226ATL0XXX00 or C226ATS0XXX00	Unimetric™ X mm Integral Set	008 (only for S) / 108 / 110 / 208 / 210 (only for S) / 308 (only for L) / 310	Box with: 1 Unimetric™ Wrench 1 Unimetric™ X mm - Special Penetration Drill 1 Unimetric™ mm - Special Calibrating Drill 1 Unimetric™ mm - Special Precalibrating Drill (only for X10 sizes) 20 Unimetric™ X.X mm N°XXX L/S

SKU	Product	Available sizes XXX	Packaging and content
C226T000XXX00	Unimetric™ X mm Set	X10 / X08	Box with: 3 Unimetric™ Wrench 1 Unimetric™ X mm - Special Penetration Drill N°3XX 1 Unimetric™ X mm - Special Penetration Drill N°2XX 1 Unimetric™ X mm - Special Penetration Drill N°1XX 1 Unimetric™ X mm - Special Penetration Drill N°008 (only for X08 size) 1 Unimetric™ X mm - Special Calibrating Drill N°3XX 1 Unimetric™ X mm - Special Calibrating Drill N°2XX 1 Unimetric™ X mm - Special Calibrating Drill N°1XX 1 Unimetric™ X mm - Special Calibrating Drill N°008 (only for X08 size) 1 Unimetric™ X mm - Special Precalibrating Drill (only for X10) size) 15 Unimetric™ L X.X mm N°3XX S 15 Unimetric™ L X.X mm N°3XX L 20 Unimetric™ L X.X mm N°2XX S 20 Unimetric™ S X.X mm N°2XX L 25 Unimetric™ S X.X mm N°1XX S 25 Unimetric™ S X.X mm N°1XX L 25 Unimetric™ S X.X mm N°008 (only for X08 size)
C226TR00XXX00	Unimetric™ X mm Set	X10 / X08	Box with: 3 Unimetric™ Wrench 1 Unimetric™ X mm - Special Penetration Drill N°3XX 1 Unimetric™ X mm - Special Penetration Drill N°2XX 1 Unimetric™ X mm - Special Penetration Drill N°1XX 1 Unimetric™ X mm - Special Penetration Drill N°008 (only for X08 size) 1 Unimetric™ X mm - Special Calibrating Drill N°3XX 1 Unimetric™ X mm - Special Calibrating Drill N°2XX 1 Unimetric™ X mm - Special Calibrating Drill N°1XX 1 Unimetric™ X mm - Special Calibrating Drill N°008 (only for X08 size) 1 Unimetric™ X mm - Special Precalibrating Drill (only for X10 size) 5 Unimetric™ L X.X mm N°3XX S 5 Unimetric™ L X.X mm N°3XX L 5 Unimetric™ L X.X mm N°2XX S 5 Unimetric™ S X.X mm N°2XX L 5 Unimetric™ S X.X mm N°1XX S 5 Unimetric™ S X.X mm N°1XX L 5 Unimetric™ S X.X mm N°008 (only for X08 size)

CE 2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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ANNEX

MRI Safety Information for Metallic Dental Posts (SST and Ti)

Dentsply Sirona Metallic Posts made of Stainless Steel (SST) or Titanium (Ti) are MR Conditional.		**MR Conditional:** Indicate an item with demonstrated safety in the MR environment within defined conditions including conditions for the static magnetic field, the time-varying gradient magnetic fields and the radiofrequency fields.

A patient with Dentsply Sirona Stainless Steel or Titanium Dental Posts can be safely scanned under the following conditions.

Name/Identification of device	Dentsply Sirona Stainless Steel and Titanium Dental Posts
Nominal value(s) of Static Magnetic Field [T]	1.5 T or 3 T
Maximum Spatial Field Gradient [T/m and gauss/cm]	40 T/m (4000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil
Maximum Whole Body SAR [W/kg]	2.0 W/kg (Normal Operating Mode)
Maximum Head SAR [W/kg]	3.2 W/kg (Normal Operating Mode)
Limits on Scan Duration	2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact of 69 mm.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

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Dentsply Sirona Europe GmbH