*Radix-Anker® - Standard, Radix-Anker®, Unimetric™, Tenons RS Radix-Stifte™*

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# Summary of safety and clinical performance

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device.

The SSCP is not intended to replace the Instructions For Use (IFU) as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The following information is intended for users/healthcare professionals.

## 1. DEVICE IDENTIFICATION AND GENERAL INFORMATION

### 1.1 Device trade name including any model number or version

| | | |
|---|---|---|
| Trade Name(s): | SST | Radix-Anker® - Standard |
| | Ti | Radix-Anker®<br>Unimetric™<br>RS Radix-Stifte™ Posts |

**Note**: The RS Radix-Stifte™, Unimetric™ and Radix-Anker® are sold alone or as component of the following set:

RS Radix-Stifte™ Posts Integal Set
RS Radix-Stifte™ Posts Set
Unimetric™ T Demo Kit
Unimetric™ X mm Integral Set
Unimetric™ X mm Set
Radix-Anker® Long Set
Radix-Anker® Standard Demo Kit
Radix-Anker® Standard Integral Set
Radix-Anker® Standard Set
Radix-Anker® Long Integral Set

### 1.2 Manufacturer's name and address

| | |
|---|---|
| Name: | Maillefer Instruments Holding Sàrl |
| Address: | Chemin du Verger 3 |
| City / Postal Code: | CH-1338 Ballaigues |
| Country: | Switzerland |

### 1.3 Manufacturer single registration number (SRN)

| | |
|---|---|
| SRN No: | CH-MF-000016301 |

### 1.4 Basic UDI-DI

| | | |
|---|---|---|
| Basic UDI-DI: | SST | ++J00310058EC |
| | Ti | ++J00310059EE |

### 1.5 Medical device nomenclature

| | | |
|---|---|---|
| Description: | MD-code and terms: | 0402 – Dental materials |
| | GMDN and terms: | 38609 – Root anal post, preformed |
| | EMDN and terms: | Q010199 – Devices for conservative dentistry and endodontics - other |

### 1.6 Class of device

| | |
|---|---|
| Risk class: | According to Annex VIII of Regulation (EU) 2017/745 of 5. April 2017 on medical devices, Metallic posts are classified as a class IIa medical device (rule 8). |

### 1.7 Year when the first certificate (CE) was issued covering the device

| | |
|---|---|
| Year: | RS RADIX Root Post – 1989<br>Unimetric™ X.X mm – 1985<br>Radix-Anker® Standard (Ti), Radix-Anker® Long (Ti), Radix-Anker® Standard (SST)- 1975 |

### 1.8 Authorized representative if applicable; name and SRN

| Authorized representative: | **European Representative** | **United Kingdom Representative** |
|---|---|---|
| Name (and Address): | Dentsply DeTrey GmbH<br>De-Trey-Straße 1<br>78467 Konstanz<br>Germany | Dentsply IH Limited,<br>Brunel Way,<br>Stonehouse,<br>Gloucestershire,<br>GL10 3GB |
| SRN No: | DE-AR_000005904 | N/A |

### 1.9 Notified body (NB) involved in the conformity assessment of the device

| | |
|---|---|
| Notified Body (NB): | BSI The Netherlands B.V. |
| NB single identification number (SIN) | CE 2797 |

# 2. INTENDED USE OF THE DEVICE

## 2.1 Intended purpose

Device intended to retain the core material in a tooth with extensive loss of coronal structure.

## 2.2 Indication(s) and target population(s)

The products are indicated for treatment of endodontic disease.

The target population is all patients who need root canal treatment can be treated with these devices, without restrictions in terms of gender and origin.

## 2.3 Contraindications and/or limitations


Insufficient residual dentin: at least 2 mm of tooth structure is required around the preparation.

# 3. DEVICE DESCRIPTION

## 3.1 Description of the device

Description:

Posts are inserted in the pulpo-dentinal structure of the tooth to support and secure the coronal part for the restoration in case of insufficient residual tooth substance (< 4 mm). The practitioner shall preliminarily drill a space in the coronal part of the tooth with the corresponding drill for post instrument, then a precision drilling will be performed with the right drill. Then posts shall be inserted in the hole created and sealed with cement.

| Device name | Image of the device |
|---|---|
| **RS RADIX Root Post** |  |
| **Unimetric™ X.X mm** |  |
| **Radix-Anker® Standard (SST)** |  |
| **Radix-Anker® Standard (Ti)** |  |
| **Radix-Anker® Long (Ti)** |  |

Application:

| | | |
|---|---|---|
| ***Single Use*** | ☒ Yes | ☐ No |

Information about Sterilization:

| | | |
|---|---|---|
| ***Sterile (Yes/No)*** | ☐ Yes<br><br>**Method of sterilization:**<br>☐ EO Gas<br>☐ Radiation<br>☐ Steam or Dry Heat<br>☐ Aseptic Techniques | ☒ No |

Information about constituents:

| | | |
|---|---|---|
| ***Medicinal substance part of the device*** | ☐ Yes<br>Please list: | ☒ No |
| ***Contains tissue / cells of human origin or their derivates*** | ☐ Yes<br>Please list: | ☒ No |
| ***Contains tissue / cells of animal origin or their derivates*** | ☐ Yes<br>Please list: | ☒ No |
| ***Device is intended to be absorbed or locally dispersed in the human body*** | ☐ Yes | ☒ No |
| ***Device contain CMR or endocrine substance*** | ☐ Yes<br>Please list: | ☒ No |
| ***Device contains materials that could result in sensitization or an allergic reaction by the patient or user*** | ☐ Yes<br>Please list: | ☒ No |

## 3.2 A reference to previous generation(s) or variants if such exist, and a description of the differences

| Previous generation / variant | Description of differences to current generation/variant |
|---|---|
| Not applicable | Not applicable |

### 3.3 Description of any accessories, which are intended to be used in combination with the device

No accessories.

### 3.4 Description of other devices and products, which are intended to be used in combination with the device

| Other devices and products intended to be used in combination with the device |
|---|
| Posts can be used with CE mark medical devices (according to their intended purpose and manufacturer recommendations) intended to be used during the treatment such as Drills (e.g., Radix-Anker® Drills). See section **6) compatible device** and **Step by step instructions** in Radix Anker IFU, RS Radix-Stifte™ IFU, Unimetric™ IFU, for details about compatible devices and in which step of the treatment they must be used.<br>Post can also be sold as part of the set with the following corresponding trade names components:<br>- Gauge<br>- Wrench<br>- Mandrel<br>See IFU section 4) in IFU for more details. |

## 4. RISKS AND WARNINGS

### 4.1 Residual risks and undesirable effects

The residual risks are:

- Temporary post-operative discomfort / pain for the patient resulting from the overall endodontic treatment*[1]
- Post debonding / fracture, leading to the swallowing / inhalation of the post / restoration in case of use error
- Root fracture, leading to the swallowing / inhalation of the post / restoration in case of use error

*[1] *Note: it can be reasonably accepted that a temporary post-operative discomfort/pain is related to the endodontic process and that no matter how competent the dentist and how modern the equipment and techniques, some pain is unavoidable due to the inflammatory and infectious status of the tooth.*

### 4.2 Warnings and precautions

- Strictly follow this Instruction for Use to minimize following risks to devices, the patient and/or the user:
  - Breakage of instrument. ○ Contamination and cross-contamination.
  - Swallowing of working part of the instrument.
  - Post debonding.
  - Toxic or allergic reactions caused by processing residues.
- The use of inappropriate rotation speed, pressure and contact duration would speed up the wear of drills for posts, increase the risk of breakage and heat creation.
- *Please refer to the Lentulo® brand IFU for the Warnings related to these instruments available on Dentsply Sirona website at www.dentsplysirona.com/ifu.

* only applicable for devices sold as component of RS Radix-Stifte™ Posts Integral Set, RS Radix-Stifte™ Posts Set, Radix-Anker® Long Set, Radix-Anker® Long Integral Set, Radix-Anker® Standard Demo Kit, Radix-Anker® Standard Integral Set, and Radix-Anker® Standard Set post sets.

### 4.3 Other relevant aspects of safety, including a summary of any field safety corrective action (FSCA including FSN)

There was no field safety corrective action between April 2019 and September 2023.

## 5. SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP (PMCF)

### 5.1 Summary of clinical data related to equivalent device

Equivalence has been shown between **Radix-Anker® Long (Ti)** and **Radix-Anker® Standard (Ti)** with regard to their clinical, technical, and biological characteristics, with special attention to the performance and principles of operation. Based on the same intended use, the same material, the same manufacturing process, and the part of the treatment they are dedicated for, the devices result as equivalent. The sole difference is the length and diameter of the head and the total length of the posts.

No equivalence is claimed for **Radix-Anker® Standard (SST), RS RADIX Root Post and Unimetric™ X.X mm.**

Clinical evidence identified for the equivalent device **Radix-Anker® Standard (Ti)** (four published clinical studies, see section 5.3) should be deemed relevant for **Radix-Anker® Long (Ti)**.

| Product | Reference | Legal Manufacturer | Basic UDI-DI |
|---|---|---|---|
| **Radix-Anker® Long (Ti)** | C026100000X00 | Maillefer Instruments Holding Sàrl<br>Chemin du Verger 3<br>CH-1338 Ballaigues<br>Switzerland | ++J00310059EE |
| **Radix-Anker® Standard (Ti)** | C193T0000XX00 | | |

### 5.2 Summary of clinical data from conducted investigations of the device before the CE-marking

N/A, no pre-market clinical investigation has been performed for these products.

### 5.3 Summary of clinical data from other sources

The following data were collected by the manufacturer in order to provide the clinical safety and performance of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts:

- Clinical data from published literature retrieved for **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, and **Radix-Anker® Standard (Ti)** (applicable to **Radix-Anker® Long (Ti)** due to equivalence demonstration) posts
- Post Market Surveillance data
- Data from Advance Event databases

*Table 1 Comparison of the main parameters for evaluation of benefit/risk ratio for Metal Posts*

| Source of evidence | Clinical success rate | Clinical survival rate | Failure rate | Catastrophic failure (tooth loss) rate | Post-operative pain |
|---|---|---|---|---|---|
| **Acceptance criteria** | | | | | |
| State-of-the-Art | 76.24-97.56%, follow-up time 6 months up to 17 years | 71-100%, follow-up time 6 months up to 17 years | 2.44-23.76% follow-up time 6 months up to 17 years | 0-29% follow-up time 6 months up to 17 years | Rate not described in literature |
| **RS Radix Root Post** | | | | | |
| Scientific literature [5] | 87%<br>**Acceptance criteria is met.** | 81.5%<br>**Acceptance criteria is met.** | 13%<br>**Acceptance criteria is met.** | 18.5%<br>**Acceptance criteria is met.** | No mention of post-operative pain reported |
| PMS data / safety databases | N/A | N/A | N/A | No cases of tooth extraction reported | No complaints or adverse events reporting on pain |
| **Radix-Anker® Standard (Ti)** | | | | | |
| Scientific literature [4-7] | 80-93.9%<br>**Acceptance criteria is met.** | 81.5-100%<br>**Acceptance criteria is met.** | 6.1-20%<br>**Acceptance criteria is met.** | 0-18.5%<br>**Acceptance criteria is met.** | Successful cases showed no signs of pain |
| PMS data / safety databases | N/A | N/A | N/A | 1 complaint (0.0027%, also reported as adverse event) reporting on tooth extraction | 1 complaint (0.0027%, also reported as adverse events) where status of the patient tooth is unknown and therefore post-operative pain cannot be ruled out |
| **Radix-Anker® Long (Ti)** | | | | | |
| Scientific literature | N/F<br>Equivalent device:<br>80-93.9%<br>**Acceptance criteria is met.** | N/F<br>Equivalent device:<br>81.5-100%<br>**Acceptance criteria is met.** | N/F<br>Equivalent device:<br>6.1-20%<br>**Acceptance criteria is met.** | N/F<br>Equivalent device:<br>0-18.5%<br>**Acceptance criteria is met.** | N/F |
| PMS data / safety databases | N/A | N/A | N/A | No cases of tooth extraction reported | No complaints or adverse events reporting on pain |
| **Radix-Anker® Standard (SST)** | | | | | |

| | | | | | |
|---|---|---|---|---|---|
| Scientific literature [5-7] | 87-93.9%<br>**Acceptance criteria is met.** | 81.5-100%<br>**Acceptance criteria is met.** | 2.1-13%<br>**Acceptance criteria is met.** | 0-18.5%<br>**Acceptance criteria is met.** | No mention of post-operative pain reported |
| PMS data / safety databases | N/A | N/A | N/A | No cases of tooth extraction reported | No complaints or adverse events reporting on pain |
| **Unimetric™ X.X mm** | | | | | |
| Scientific literature [2,3,6] | 81.8-100%<br>**Acceptance criteria is met.** | 70.5-100%<br>**Acceptance criteria is considered met.*** | 0-18.2%<br>**Acceptance criteria is met.** | 0-29.5%<br>**Acceptance criteria is considered met.**** | No mention of post-operative pain reported |
| PMS data / safety databases | N/A | N/A | N/A | No cases of tooth extraction reported | No complaints or adverse events reporting on pain |

*The clinical survival rate associated with the use of Unimetric™ X.X mm posts as reported in study [2] is slightly lower when compared to the benchmark "clinical survival rate" derived from the SotA acceptance criteria. This discrepancy may be explained by the follow-up time in this study (18-22 years) being longer than the follow-up time of the benchmark value from the state-of-the-art

**The catastrophic failure (tooth loss) rate associated with the use of Unimetric™ X.X mm posts as reported in study [2] is slightly higher when compared to the benchmark "catastrophic failure (tooth loss) rate" derived from the SotA acceptance criteria. This discrepancy may be explained by the follow-up time in this study (18-22 years) being longer than the follow-up time of the benchmark value from the state-of-the-art

Conclusion:

A systematic literature review on clinical performance and safety of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)** (equivalent device to **Radix-Anker® Long (Ti)**), and **Radix-Anker® Long (Ti)** posts was performed. Results on clinical outcomes are summarized in Table 1. No published clinical data on **Radix-Anker® Long (Ti)** was identified by the systematic literature search.

The PMS data showed the broad use and acceptance of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts (which came to market between 1975 and 1989). Between April 2019 to September 2023 the number of complaints reported by the customers was low.

**RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts are sold worldwide. The main safety and recall databases in Europe (MHRA, BfArM, ANSM, SwissMedic) as well as the database MAUDE in United States (Dentsply Sirona, as a US company, has to report events to FDA) were therefore searched.

Safety databases searches identified 2 relevant events for **Radix-Anker® Standard (Ti)** posts, both of them related to post breakage. One event led to tooth extraction while the patient tooth's status is unknown in the second event. No relevant adverse events were identified for **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, and **Radix-Anker® Long (Ti)** posts neither in European adverse event databases nor in the US MAUDE database over the period of search (01.11.2019 – 20.12.2024).

No new or emerging risks, or significant increases in severity or frequency of known risks have been identified during this time period. No information has been received that would impact the benefit-risk conclusion identified during the risk management process. All existing risks were already addressed, mitigated and rated as acceptable in regard to the clinical benefits for the patients and are continuously monitored through Post Market Surveillance activities according to Dentsply Sirona Standard Operating Procedures in place. The benefits of Metallic Posts have been confirmed to outweigh the potential risks and the risk control measures are deemed effective.

Based on the outcome of the clinical evaluation, it can be stated that the performance and safety profile of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)** (equivalent device to **Radix-Anker® Long (Ti)**), and **Radix-Anker® Long (Ti)** posts are considered to be in line with the performance and safety profile of state-of-the-art similar devices, belonging to the same technology. All clinical data appraised confirm the performance and safety of the family of products when used as intended and demonstrate that the benefits of using these devices outweigh the risks, as required by the MDR 2017/745.

References:

[2] Caserio Valea, M. and V. Alonso de la Pena, Titanium posts and bonded amalgam core longevity: A 22-year clinical survival retrospective study. J Am Dent Assoc, 2017. 148(2): p. 75-80.

[3] Chekhani, U.N., A.A. Mikeli, and F.K.P. Huettig, All-ceramic prosthetic rehabilitation of a worn dentition: Use of a distal cantilever. Two-year follow-up. Dental research journal, 2013. 10(1): p. 126-131.

[4] Ellner, S., T. Bergendal, and B. Bergman, Four post-and-core combinations as abutments for fixed single crowns: a prospective up to 10-year study. Int J Prosthodont, 2003. 16(3): p. 249-54.

[5] Fokkinga, W.A., et al., Composite resin core-crown reconstructions: an up to 17-year follow-up of a controlled clinical trial. Int J Prosthodont, 2008. 21(2): p. 109-15.

[6] Mentink, A.G., et al., Clinical performance of different post and core systems--results of a pilot study. J Oral Rehabil, 1993. 20(6): p. 577-84.

[7] Ottl, P. and H.C. Lauer, Success rates for two different types of post-and-cores. J Oral Rehabil, 1998. 25(10): p. 752-8.

## 5.4 An overall summary of the clinical performance and safety

Based on available data, it has been judged that **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts are safe when used according to their intended purpose and perform as indicated by the manufacturer when used according to their corresponding IFUs. The data provided by the manufacturer is in line with the state of the art, and the information provided to users and patients is appropriate. The Post Market Surveillance activities confirm that these devices are used by dental professionals during endodontic treatments and the incident/complaint rate is very low when compared to the number of units sold. No new clinical risks were identified. Biological risks and risks related to substances were assessed in the biological safety assessment. The conclusion of the analysis showed that there is no substance detected from raw materials or chemical residues in sufficient concentration to be considered harmful to the end user. All other foreseeable hazards and risks related to the manufacturing process, degree of innovation and principle of operation of the device are taken into account in the Risk Management documents. Those risks have been addressed, mitigated appropriately and were considered acceptable. Moreover, Metallic Posts have been on the European market for more than a decade.

Conclusion:

The available clinical data provided sufficient evidence to support the safety and performance of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts. The benefit for the patient to restore tooth integrity and function as part of endodontic treatment (in case of insufficient residual tooth structure) when using the Metallic Posts clearly outweighs the potential risks associated with the endodontic treatment. All residual risks have been mitigated and are considered as acceptable. All possible mitigation actions have been implemented.

Further, considering that the alternative treatment would be a tooth extraction, the benefit of keeping the tooth and its function is significantly higher. Conclusively, the overall benefit/risk profile is judged to be acceptable.

## 5.5 Ongoing or planned post-market clinical follow-up

It has been judged that there is sufficient clinical evidence presented in the Clinical Evaluation Report to demonstrate that the safety and performance of the Metallic Posts are in line with the relevant General Safety and Performance Requirements. The current evaluation did not reveal any issue that cannot be identified by the post-market surveillance actions and procedures already in place. As part of ongoing PMS activities, the manufacturer will continue to collect evidence to monitor the safety and performance of **RS RADIX Root Post**, **Unimetric™ X.X mm**, **Radix-Anker® Standard (SST)**, **Radix-Anker® Standard (Ti)**, and **Radix-Anker® Long (Ti)** posts.

There are no uncertainties or unanswered questions towards performance, safety, and side effects. No risks remain unaddressed, and the overall residual risk is acceptable. Therefore, no additional Post-Market Clinical Follow-up (PMCF) studies are required and are currently planned.

PMCF Plan dedicated to the devices was established to define the methods and procedures for proactively collecting and evaluating clinical data of Metallic Posts. Based on the clinical data available, a literature vigilance search is planned, including a literature search and safety databases search for the Metallic Posts and similar devices.

## 6. POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

Restoring endodontically treated teeth with significant loss of coronal tooth structure often necessitates the use of a core retained with a post [271]. Posts enhance core retention but do not strengthen the tooth. They are classified as metallic or nonmetallic, with retention primarily attributed to the apical extension anchoring to the core material [36].

Intraradicular posts have been used for over 250 years to restore endodontically treated teeth, which are more prone to fractures due to changes in dentin properties, loss of proprioception, and structural loss from caries and endodontic access [36]. Since their introduction in the 1990s, fiber-reinforced composite (FRC) posts have become a popular alternative to traditional cast or prefabricated metal posts [37].

Prefabricated posts are selected to fit the available radicular space after conservative canal preparation, minimizing the removal of sound tooth structure. Custom-made posts, although requiring the removal of more dental tissue, have good survival outcomes [285, 286].

Fiber posts

Nonmetallic posts, especially glass fiber posts, have a closer elastic modulus to dentin, reducing stress concentrations and improving esthetics. They are advantageous because they can be readily retrieved for restoration replacement or endodontic retreatment [36]. However, they are prone to debonding at the dentin-cement interface. The primary function of fiber posts is to ensure the retention of restorations after significant dental structure loss, aiding in the secure filling and buildup of a prosthetic core [213]. However, fiber post debonding is influenced by several factors, including residual tooth tissue, occlusal scheme, periodontal status, presence of ferrule, and adhesion quality. The presence of at least 2 mm of coronal dentin, known as the "ferrule effect", significantly enhances restoration fracture resistance and supports better biomechanical performance and longevity of post-and-core restorations [37, 66].

Zirconium posts

Zirconium-oxide posts have been introduced to offer both esthetics and long-term durability. These posts are often combined with direct composite resin cores or indirect ceramic core materials, and can be either prefabricated or designed and milled using CAD/CAM technology for better adaptation to canal walls [182].

Metal cast posts

Cast post-and-cores made from precious or non-precious metal alloys, offer strong retention but can lead to stress concentration and root fractures due to their high elastic modulus. A cast post-and-core system is a (self-cutting) cast dowel or post in combination with a cast core build-up and was the previous gold standard in endodontic treatment before the emergence of fiber posts and composite core build-up materials. Currently, they only represent an option in case of teeth with very limited residual structure (one wall or less with/without ideal ferrule) [258]. This system relies more on macromechanical retention and the fit of the post within the canal. For its cementation, both adhesive and non-adhesive cements can be used.

| Type of Post | Advantages over Prefabricated Metallic posts | Disadvantages compared to Prefabricated Metallic posts |
|---|---|---|
| Fiber post | • Aesthetic advantage (improved esthetics for reproducing the natural aspect of the restored tooth) | • Similar survival rate<br>• Less radiopacity |

| | | |
|---|---|---|
| | • Elastic modulus similar to dentine | |
| Metal cast post | • Better geometric adaptation to excessively flared or elliptical canals | • Required 2 visits and laboratory fee<br>• Difficult and invasive to remove<br>• Requires minimum tooth structure removal |
| Zirconium post | • High flexural strength,<br>• High fracture toughness | • Retrieval of zirconium and ceramic posts is very difficult if endodontic retreatment is necessary or if the post fractures |
| No post | • Preserves remaining tooth structure; avoids risks associated with post placement<br>• Simpler and less invasive procedure | • Limited to cases with sufficient coronal tooth structure remaining<br>• Less retention for core and crown in severely damaged teeth<br>• Higher failure rates if minimal tooth structure remains |

References

[36] Giok, K.C., S.K. Veettil, and R.K. Menon, Comparative effectiveness of fiber and metal posts in the restoration of endodontically treated teeth: A systematic review with network meta-analysis. J Prosthet Dent, 2023.

[37] Yang, A., A. Lamichhane, and C. Xu, Remaining coronal dentin and risk of fiber-reinforced composite post-core restoration failure: a meta-analysis. Int J Prosthodont, 2015. 28(3): p. 258-64.

[66] Batista, V.E.S., et al., Influence of the ferrule effect on the failure of fiber-reinforced composite post-and-core restorations: A systematic review and meta-analysis. J Prosthet Dent, 2020. 123(2): p. 239-245.

[182] Al-Thobity, A.M., Clinical Outcomes of Zirconium-Oxide Posts: Up-to-Date Systematic Review. Compend Contin Educ Dent, 2016. 37(6): p. e17-20.

[213] Sorrentino, R., et al., Complications of endodontically treated teeth restored with fiber posts and single crowns or fixed dental prostheses-a systematic review. Clin Oral Investig, 2016. 20(7): p. 1449-57.

[258] Perdigão, J., Restoration of Root Canal-Treated Teeth. 2016: Springer International Publishing AG Switzerland.

[271] Tsintsadze, N., et al., Comparing survival rates of endodontically treated teeth restored either with glass-fiber-reinforced or metal posts: A systematic review and meta-analyses. J Prosthet Dent, 2024. 131(4): p. 567-578.

[285] Al-Dabbagh, R.A., et al., Effect of a circumferential ferrule on the survival and success of endodontically treated teeth restored with fiber posts: A systematic review and meta-analysis. J Prosthet Dent, 2024.

[286] Martins, M.D., et al., Is a fiber post better than a metal post for the restoration of endodontically treated teeth? A systematic review and meta-analysis. J Dent, 2021. 112: p. 103750.

## 7. SUGGESTED PROFILE AND TRAINING FOR USERS

The suggested profiles are:

- Endodontist: Qualified dental professionals specialized in Endodontic treatment. The Endodontist is the expert for complex cases. He/she practices endo treatment each day.
- General Practitioner Dentist: Qualified dental professionals not specialized in Endodontic treatment. The general practitioner practices endo treatment as well as any other kind of dental treatment.
- Dental Assistant: his/her role is processing and preparing the system for use. The dental assistants can also help the dentists during the treatment by preparing, holding or/and passing instruments.

No training is required by the manufacturer, the preparation and use of Metallic Posts is part of the common training of the dentist. Recommended processing instructions are described in the Instructions for Use.

## 8. REFERENCE TO ANY HARMONIZED STANDARDS AND CS APPLIED

| | | |
|---|---|---|
| **EN ISO 10993-9** | **2021** | Biological evaluation of medical devices<br>- Part 9: Framework for identification and quantification of potential degradation products (ISO 10993-9:2019) |
| **EN ISO 10993-10** | **2023** | Biological evaluation of medical devices<br>- Part 10: Tests for irritation and skin sensitization (ISO 10993-10:2021) |
| **EN ISO 10993-12** | **2021** | Biological evaluation of medical devices<br>- Part 12: Sample preparation and reference materials (ISO 10993-12:2021) |
| **EN ISO 10993-15** | **2023** | Biological evaluation of medical devices<br>- Part 15: Identification and quantification of degradation products from metals and alloys |
| **EN ISO 10993-17** | **2023** | Biological evaluation of medical devices<br>- Part 17: Establishment of allowable limits for leachable substances (ISO 10993-17:2023) |
| **EN ISO 10993-18/ A1** | **2023** | Biological evaluation of medical devices<br>- Part 18: Chemical characterization of medical device materials within a risk management process<br>- Amendment 1: Determination of the uncertainty factor (ISO 10993-18:2020) |
| **EN ISO 10993-23** | **2021** | Biological evaluation of medical devices<br>- Part 23: Tests for irritation (ISO 10993-23:2021) |
| **EN ISO 13485/ AC/ A11** | **2021** | Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016) |
| **EN ISO 14971/ A11** | **2021** | Medical devices - Application of risk management to medical devices (ISO 14971:2019) |
| **EN ISO 15223-1** | **2021** | Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied (ISO 15223-1:2021) |
| **EN ISO 17664-1** | **2021** | Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices (ISO 17664-1:2021) |

## 9. REVISION HISTORY

| SSCP REVISION NUMBER | DATE ISSUED | CHANGE DESCRIPTION | REVISION VALIDATED BY THE NOTIFIED BODY |
|---|---|---|---|
| *0* | *30.06.2025* | *First release* | ☐ Yes<br>Validation Language:<br>☒ No<br>Explanation: This document is not yet approved by the Notified Body. |

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