Dentsply Sirona

PathFile®

EN	Instructions For Use	2
FR	Instructions d'utilisation	7
IT	Istruzioni per l'uso	13
DE	Gebrauchsanweisung	19
ES	Instrucciones de uso	25
PT	Instruções de utilização	31
NL	Gebruiksinstructies	37
SV	Användarinstruktioner	43
DA	Brugsanvisning	48
NO	Brukerinstruksjoner	53
LT	Naudojimo instrukcija	58
LV	Lietošanas norādījumi	63
ET	Kasutusjuhised	68
CS	Návod k použití	73
SK	Návod na použitie	78
HU	Használati utasítás	84
RO	Instrucţiuni de utilizare	90
SR	Uputstva za upotrebu	96
BS	Uputstva za upotrebu	101
CNR	Uputstva za upotrebu	106
SQ	Udhëzimet për përdorim	111
EL	Οδηγίες χρήσης	117
TR	Kullanım Kılavuzu	123
PT-BR	Instruções de uso	128
KO	사용 지침	134

Dentsply Sirona

EN

PathFile®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED / REUSABLE – LIABLE FOR REPROCESSING
STERILE – STERILIZED BY RADIATION / NON-STERILE – STERILIZE BEFORE FIRST USE

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS
REUSABLE NON-STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 PathFile® are engine-driven endodontic tapered instruments used during the endodontic treatment for the preparation of the root canal (glide path) during endodontic therapy. They are designed with an operative part of nickel titanium alloy with cutting edges and a standardized shank. Multiple use PathFile® devices are available in non-sterile version.
- Clinical benefits:
 The clinical benefit for the patient being treated with PathFile® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (glide path).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

The PathFile® instruments are delivered in both versions: sterile blisters or non-sterile plastic boxes.

The PathFile® range comprises the following endodontic instruments:

Code	Product	ISO XXX YYY XXX-ISO size YYY-taper size	Available lengths	Packaging and content
A 0015	PathFile® STERILE	P1: ISO 013 002 P2: ISO 016 002 P3: ISO 019 002	21/25/31 mm	Blister with 6 instruments (sterile)
	PathFile® STERILE Assortment	P1 P2 P3	21/25/31 mm	Blister with 2 instruments of each ISO (sterile)
	PathFile®	P1: ISO 013 002 P2: ISO 0160 02 P3: ISO 019 002	21/25/31 mm	Plastic box with 6 instruments (non-sterile)
	PathFile® Assortment	P1 P2 P3	21/25/31 mm	Plastic box with 2 instruments of each ISO (non-sterile)
	PathFile® Assortment	P1 P2 P3	21/25 mm	Plastic box with 1 instrument of each ISO (non-sterile)

5 COMPOSITION

PathFile® instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with a colored ring.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 300 rpm.
- Set torque at: 0.6 - 1 Ncm.
- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.

- Cross-contamination.
- Heat generation due to insufficient lubrication and irrigation.
- Swallowing of working part of the instrument.
- Toxic or allergic reactions caused by processing residues.

- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, PathFile® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).

- Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
- Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
5) Scout the root canal with an ISO 010 hand file, determine the working length by using an electronic apex locator in combination with radiographs, irrigate and confirm patency.
6) Establish a manual Glide Path with hand file, like stainless steel K-File, ISO 008 and ISO 010.
7) Use first PathFile® ISO 013 002 (P1) and then PathFile® ISO 016 002 (P2) to working length.
8) Finally use PathFile® ISO 019 002 (P3) to working length (not needed before shaping procedure with ProTaper Next®).
9) Irrigate the expanded glide path.
10) Reconfirm your working length with a K-File ISO 015 combined with an apex locator, before starting the shaping procedure with NiTi ProTaper® files or another system.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- For instruments that are labelled “sterile”:
 - PathFile® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
 - Disinfection, cleaning and sterilization are not applicable.
- For instruments that are labelled “non-sterile”:
 - Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.

- Do not re-use silicone stops. Remove and discard silicone stops after each use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH